
# LABORGLASKATALOG LABORATORY GLASSWARE CATALOGUE

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*“I'D LIKE TO SAY THANKS FOR THE CONSTANT  
SERVICE, FOR THE RELIABILITY, FOR **THE MANY  
DIFFERENT OPTIONS** ... IN SHORT, CAN  
ANYBODY TELL ME HOW I CAN SAY THANKS  
TO SOME GLASSWARE?”*

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# WELCOME TO THE DURAN GROUP.

We are delighted with your interest in our chemical and heat-resistant DURAN® laboratory glassware. Based on the know-how of more than 650 employees at sites in Mainz and Wertheim in Germany and Pula in Croatia and with more than 100 years experience in manufacturing borosilicate glass 3.3, you will find we are the right contact for standard laboratory glassware as well as exacting, customized, special articles.

With our comprehensive range and our wide range of production and finishing options, we are sure we can find, in cooperation with you, the optimum laboratory glass for your special application.

This product catalogue provides an overview of our products and services within the field of DURAN® laboratory glassware. Our product managers and sales staff will be happy to discuss your individual needs as well as the properties of our speciality material, DURAN®.

You can find more information and details for your contact partner on our website [www.duran-group.com](http://www.duran-group.com).

Please do not hesitate to contact us.

Michael Merz  
General Manager Marketing & Sales

Tobias A. Thiele  
Vice President Marketing

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**DURAN GROUP. Magic of precision.**

## THERE ARE GOOD REASONS FOR GETTING TO KNOW US.

*The DURAN GROUP is one of the world's leading manufacturers of borosilicate glass, a special glass invented by Otto Schott in 1887 and registered in 1938 under the trade name DURAN®.*

*With sites in Mainz (DURAN Produktions GmbH & Co. KG), Wertheim (DURAN Group GmbH) and Pula (DURAN d.d.), the DURAN GROUP has a complete value-added chain from the glass melting process right through to glass forming and processing.*

*Over the last 70 years, DURAN® glass has proven itself both in laboratories as well as in industrial and domestic applications. At the same time our special industrial glass is used in numerous industries, amongst others, mechanical engineering, the electrical industry and medical technology. Well-known companies spread across all continents trust in the quality of our products and the creativity of our employees. Based on close dialogue with our customers, we achieve solutions, which meet individual requirements and permit highly precise applications in widely varying sectors.*

*After its spin-off from SCHOTT AG, the DURAN GROUP currently combines the strengths of its long history with the constantly present visionary spirit of the inventor Otto Schott.*

*The unique precision involved in the manufacturing and processing of DURAN® glass is the foundation of our corporate philosophy and guideline for our day-to-day operations. We offer our customers the flexibility of a medium-sized company with comprehensive know-how built up from many years of experience.*

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
## QUALITY WITHOUT COMPROMISE.

The requirement of our customers for our company to develop and manufacture reliable and safe products in accordance with the highest possible quality standards is at the very centre of our quality policy.

In close cooperation with all our staff and with the active involvement of our customers and suppliers, the DURAN GROUP has established a quality management system that conforms to DIN EN ISO 9001:2008 and which is integrated into daily practice.

This quality management system is the foundation for all the processes that our products have to go through: from the customer's initial enquiry, through to order processing and delivery and up to customer feedback.

We value the success of DURAN® products as a sign of confidence from our customers in our quality system, in our logistics and in our service.

Moreover, we are constantly working on improving our performance by assessment of the status quo and from it deriving new, ambitious goals, which we want to attain through process-orientated thinking, planning and action.

In everything we do, customer and staff satisfaction, as well as adherence to the principles of environmental sustainability are the important criteria that we constantly measure and aim to improve upon.

### Safety and environmental protection as corporate goals.

Human safety and protection of the environment are important goals for our company. To achieve this, we are constantly working on the development of more environmentally sustainable products and manufacturing processes while ensuring we use resources in a sparing and responsible manner and set high standards for the safe operation of our plant.

The DURAN GROUP is committed to responsible and fair behaviour towards society and the environment. Legislative guidelines and official ordinance are considered a minimum requirement. In our daily operations we go well beyond this minimum.

We commit ourselves further to reduce our demand for energy and at the same time to continuously improve our energy efficiency in the long run. To put our energy objectives into effect we introduce an energy management system according to DIN EN 16001:2009. Here we make sure that all requirements of this standard will be put into effect correctly and that the processes within this energy management system are subject to continuous improvement.

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## DURAN® LABORATORY GLASSWARE WITH BATCH IDENTIFICATION AND QUALITY CERTIFICATE.

*Increasing requirements and standards in terms of quality assurance as well as the traceability of primary packaging and auxiliary materials (EN-ISO 9000/8402, GMP, EU 178/2002) are becoming more and more significant.*

*To meet these requirements, all DURAN® laboratory bottles, beakers, Erlenmeyer flasks and premium caps are provided with a **Retrace Code**.*

*This eight-digit batch code (two digits for premium caps) allows the products to be retraced to the point of production and the matching batch. The retrace code is a contribution to a continuous documentation within the user's quality management system and it is therefore particularly important in the areas of medicine, the pharmaceutical and food industries.*

*By using this "Retrace Code" and the corresponding DURAN® item number (directly or via a selection menu), users can readily retrieve a batch and quality certificate via the internet. Besides the production date and the signature of the Quality Manager, the certificate also provides information on conformity with standards and USP/EP/JP.*

Please see [www.duran-group.com](http://www.duran-group.com) for more information about the Retrace Code.

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“MY FAVOURITE ASSISTANT ... **CHEMICALLY RESISTANT, PHYSICALLY RESISTANT ... SIMPLY: UNBEATABLE.**”

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## I | LABORATORY GLASS BOTTLES AND ACCESSORIES

*Due to their high chemical and thermal-shock resistance, DURAN® laboratory bottles are indispensable for routine work in the laboratory. As neutral glass with high hydrolytic resistance, DURAN® belongs to glass type I in accordance with EP, JP and USP.*

*The extensive range of original accessories includes caps for the widest possible range of applications. DURAN® laboratory bottles are completed by suitable pouring rings from different plastics, which enable drip-free working.*

*Alongside clear glass laboratory bottles, amber (absorbs light up to approx. 500 nm) and plastic-coated (absorbs light up to approx. 380 nm) bottles are also available to protect light-sensitive substances. The PU plastic coating also provides protection against scratching and, in the event of breakage, retains contents and glass splinters.*

### Usage tips:

- DURAN® laboratory bottles are, with the exception of the pressure-resistant DURAN® pressure plus bottles, not designed for use under pressure or vacuum conditions.
- When sterilising / autoclaving, the screw cap should

*only be loosely applied since with a closed bottle pressure equalisation does not take place. An ideal complement is the membrane cap (see pages 20, 34).*

- *When cleaning in a dishwasher, load such that the glass bottles, especially the threads, are not in contact with each other.*
- *The printed scale on many items of DURAN® laboratory glassware is indicative with an accuracy of  $\pm 10\%$ . Therefore the items are not suitable for use as volumetric glassware.*
- *DURAN® protect bottles can be autoclaved. The maximum permissible temperature is  $+135\text{ °C}$  (max. 30 Min.).*
- *When freezing media, store the bottles in an inclined position and do not fill more than three quarters full. Bear in mind that plastic caps and pouring rings have restricted service temperature ranges.*

*DURAN® laboratory glass and reagent bottles are provided with a retrace code. Using the eight-character code and the corresponding article number, a batch and quality certificate for the DURAN® laboratory glass bottles can be obtained at [www.duran-group.com](http://www.duran-group.com).*

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
*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Special thread means opening takes less than a turn. The 80 mm wide outer diameter of the bottle mouth permits easy filling and pouring out of powders and viscous substances.*

*Typical applications: storage, transport, safekeeping and sampling of substances, easy to use with granulated material, powders and viscous media, sampling of hot media.*

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**DURAN® GLS 80® LABORATORY BOTTLE, WIDE MOUTH**  
with GLS 80® thread

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**DURAN® GLS 80® LABORATORY**  
**BOTTLE, WIDE MOUTH, AMBER**  
**with GLS 80® thread**

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Alongside easy handling, UV protection up to 500 nm. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.*

*Typical applications: storage, transport and safekeeping of light-sensitive substances, easy to use with granulated material, powders and viscous media.*

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
*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue quick release closure (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of closure and pouring ring: +140 °C. Service temperature limit of the PU<sup>1</sup> plastic coating: -30 °C to +135 °C. The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 380 nm wavelength. High transparency. Suitable for microwaving.*

*Typical applications: storage, transport and safe handling of hazardous substances. Storage of high value viscous liquids, pastes and powder.*

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**DURAN® GLS 80® PROTECT**  
**LABORATORY BOTTLE, WIDE MOUTH**  
with GLS 80® thread, plastic coated

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

**DURAN® PRODUCTION AND STORAGE BOTTLES (CARBOYS)**

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*These larger sized bottles / carboys are ideal for bulk storage and handling of both liquid and solid intermediates and final formulations. Manufactured from Type I borosilicate 3.3 glass for durable performance and resistance to thermal stress. The glass conforms to American (USP), European (EP) and Japanese (JP) pharmacopoeia standards making the carboys ideal for pharmaceutical production applications. Manufactured with thickened, uniform side walls for higher mechanical strength. Available in two mouth sizes: The versatile DIN standard GL 45 or GLS80® wide mouth for easier cleaning, filling and emptying. Retrace Code for batch traceability and conformance certification. Manufactured from inorganic materials (Certified BST / TSE free). Suitable for high temperature sterilization, depyrogenisation or autoclaving. Feature large, permanent, easy-to-read, white enamel graduations marks. Also available with an external Polyurethane coating for enhanced scratch resistance, and to contain leakage in the event of damage. Available with customized logos, identification labeling or graduations. Supplied without screw cap or pouring ring, but can be used in conjunction caps and connector systems see pages 20 ff and 33 ff for options.*

*Typical applications: Flat robust base is ideal for mixing processes with large magnetic stir bars.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread GL/GLS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
11 602 00	10 000	45	227,50	410	1
11 601 00	20 000	45	288,50	505	1
11 602 20	10 000	80	227,50	385	1
11 601 10	20 000	80	288,50	480	1

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*Permits opening and closing of the DURAN® GLS 80® bottle with only a three-quarter turn. A matching PP<sup>1</sup> pouring ring is also available, permitting clean, drip-free use.*

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
*QUICK RELEASE CLOSURE FOR DURAN® GLS 80® LABORATORY BOTTLE from PP<sup>1</sup>, blue, with lip seal*

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
**QUICK RELEASE CLOSURE FOR DURAN® GLS 80® LABORATORY BOTTLE**  
**with sealing disc, (PSU compound<sup>1</sup>)**

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*The material used is a special compound based on a polyarylsulphone I. Consequently the chemical, thermal and mechanical properties of the material are noticeably improved and matched to laboratory requirements. Thanks to the thread, the DURAN® GLS 80® bottle can be opened or closed with only a three-quarter turn. The seal, which is coated on both sides with PTFE<sup>1</sup>, ensures the bottle can be tightly sealed. A matching PTFE<sup>1</sup> pouring<sup>2</sup> ring is also available, permitting clean, drip-free use.*

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
**MEMBRANE SCREW CAP FOR DURAN® GLS 80® LABORATORY BOTTLE**  
**from PP<sup>1</sup>, blue, with welded-in PTFE<sup>1</sup> membrane**  
**for pressure equalisation**

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*For GLS 80® thread. Ideal for autoclaving processes because the membrane permits pressure equalisation and screw tight sealing. Hence the risk of contamination is greatly reduced. Ingress of liquids or solids is prevented and the bottle contents remain sterile.*

*Typical applications: storage or transport of gas generating media, autoclaving of media.*

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*Materials used: PP<sup>1</sup> and PTFE<sup>1</sup>. Flexible modular system. Five different tubing diameters (3.2 mm; 6.0 mm; 8.0 mm; 10.0 mm and 12.0 mm) can be connected. Sterile pressure equalisation is possible by using the syringe filter. Unused ports can be closed with a blind cap.*


*Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).*

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**GLS 80<sup>®</sup> CONNECTION SYSTEM**  
**screw cap GLS 80<sup>®</sup> with four ports**  
**(GL 18 thread)**

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
**GLS 80® STIRRED REACTOR**

materials used: PP'/PTFE'/PEEK'/stainless steel

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*The GLS 80® stirred reactor is suitable for a wide range of laboratory mixing processes. The connections (2 x GL 14, 2 x GL 18) provided permit addition or removal of media from the bottle during the mixing process. The whole unit can be autoclaved and is therefore suitable for biological applications. By using components from the GLS 80® connection system, an additional media bottle (OD hose: 1.6 - 12.0 mm)<sup>2</sup> can be connected or a sterile pressure equalizer attached. Drive for the stirrer is provided by a standard commercial magnetic stirrer. The variable stirrer shaft can be used in DURAN® GLS 80® laboratory glass bottles (1 000 and 2 000 ml) and provides notably improved mixing in comparison with standard magnetic stirring bars. The agitating element is exchangeable and use up to 500 rpm is possible.*

*Typical applications: mixing of liquids, dissolving of solids, simple fermentation processes.*

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**DURAN® products with individual labelling**

Innovative laser marking, e.g.: Numbering, bar codes, logos or names

No effect on the product performance. More information, see page 229

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<sup>2</sup> Components for a secure hose connection, see page 21 and 36.

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**DURAN® LABORATORY BOTTLE**  
**with DIN thread**

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. With proven DURAN® properties. Complete with blue screw cap (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature level of screw cap and pouring ring: +140 °C. For screw cap assortment, see pages 33–42.*

*Typical applications: storage, sample preparation, transport. Autoclaving media.*

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<sup>2</sup> With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).

<sup>3</sup> Bottles of 100 – 2 000 ml capacity have a shoulder reinforcement ring corresponding to the bottles nominal capacity.

<sup>4</sup> Acceptance within ISO 4796 standard has been requested.

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
*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. UV protection up to approx. 500 nm wavelength. Unchanged DURAN® properties within the bottle, as colouration is only on the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology. For screw cap assortment, see pages 33–42.*

*Typical applications: storage and transport of light-sensitive substances.*

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**DURAN® LABORATORY BOTTLE, AMBER with DIN thread**

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<sup>1</sup> With specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required).

<sup>2</sup> Bottles of 100 – 2 000 ml capacity have a shoulder reinforcement ring corresponding to the bottles nominal capacity.

<sup>3</sup> Acceptance within ISO 4796 standard has been requested.

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**DURAN® PROTECT**  
**LABORATORY BOTTLE**  
 with DIN thread, plastic coated

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Service temperature limit of the PU<sup>1</sup> plastic coating: -30 °C to +135 °C (instructions for autoclaving, see page 217). The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. UV protection up to approx. 380 nm wavelength. High transparency. Suitable for microwaving. Caps and pouring rings see pages 33–42.*

*Typical applications: storage, transport and safe handling of hazardous or valuable substances.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Ohne Verschluss und Ausgießring, klar Without cap and pouring ring, clear					
10 926 76	25 <sup>2</sup>	25	36	70	10
10 926 77	50	32	46	87	10
21 805 24	100 <sup>3</sup>	45	56	100	10
21 805 29	150 <sup>4</sup>	45	62	105	10
21 805 36	250 <sup>3</sup>	45	70	138	10
21 805 44	500 <sup>3</sup>	45	86	176	10
21 805 51	750 <sup>4</sup>	45	95	203	10
21 805 54	1 000 <sup>3</sup>	45	101	225	10
21 805 63	2 000 <sup>3</sup>	45	136	260	10
21 805 69	3 500 <sup>4</sup>	45	160	295	1
21 805 73	5 000	45	182	330	1

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
*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable ceramic. Pressure resistance conforms to DIN ISO 1595, confirmed by GS marking (TÜV ID: 0000020716). Vacuum and/or pressure resistant from -1 to +1.5 bar due to modified geometry (based on ISO 4796-1). When pressure loaded the following apply: thermal shock resistance 30 K and maximum usage temperature +140 °C. Blue scale for visual differentiation from the standard laboratory bottle. Also available in amber. See pages 33-42 for spare caps and pouring rings.*

*Typical applications: safe working under pressure or vacuum, sampling under pressure, storage of gas generating media.*

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**DURAN® PRESSURE PLUS LABORATORY BOTTLE<sup>1</sup>**  
**with DIN thread, GL 45**

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Ohne Verschluss und Ausgießring, klar Without cap and pouring ring, clear					
10 922 34	250	45	70	138	10
10 922 35	500	45	86	176	10
21 810 54	1 000	45	101	225	10
Braun   Amber					
10 943 67	250	45	70	138	10
10 943 68	500	45	86	176	10
21 816 54	1 000	45	101	225	10

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<sup>1</sup> Bottle with plastic coating is available on demand.

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
**DURAN® HPLC BOTTLE**  
**with DIN thread, GL 45**

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable ceramic. Complete system comprising: DURAN® pressure plus laboratory bottle with 4-port screw cap (PP<sup>1</sup>); four screw connections (black, M8 thread) and silicone seals. Connection of different hose diameters (1.6 mm and 3.2 mm) as well as sterile pressure equalisation sets (syringe filter 0.2 µm) is possible. Unused ports can be sealed with silicone blanking seals.*

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<sup>3</sup> Replacement bottle, see page 27.

<sup>4</sup> Including silicone seals.

<sup>5</sup> GL 45 screw cap, 4 x M8 screw cap, 12 x silicone seal.

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
*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with blue screw cap (PP<sup>1</sup>, integral lip seal) and pouring ring (PP<sup>1</sup>) for drip-free pouring and clean, safe working. Service temperature limit of screw cap and pouring ring: +140 °C. Ergonomic handling due to angular shape, highly stable, good stackability. Alongside proven DURAN® properties, a space saving of 44% in comparison with standard laboratory bottles (example applies to 100 ml bottles). Screw caps are also available in the following colours: green, yellow and grey (see pages 33–42).*

*Typical applications: space-saving storage, space-saving transport.*

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**DURAN® LABORATORY BOTTLE, SQUARE with DIN thread**

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*With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Complete with pouring ring and cap from TpCh260<sup>1</sup> (similar to PFA). The premium cap with its PTFE<sup>1</sup> coated silicone seal is colourless and temperature resistant from -196 °C to +260 °C. Together with proven DURAN® properties, TÜV tested thermal shock resistance of 160 K, confirmed by GS-marking (TÜV ID: 0000020715). USP/FDA conformity for the entire system comprised of bottle, screw cap and pouring ring. DMF available upon request (DMF no. 19757). Accurate scale: ± 5%. Additional graduations as well as additional opposing scale simplify reading off.*

*Typical applications: Due to its properties, ideal for applications in the pharmaceutical industry, handling of aggressive media, sterilisation processes (hot air and dry sterilisation) and depyrogenation.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
11 270 75	100	45	56	105	10
11 270 76	250	45	70	143	10
11 270 77	500	45	86	181	10
11 270 78	1 000	45	101	230	10

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
*The unbreakable DURAN Group steel bottle is ideal for storage applications where glass is not applicable, due to the risk of breakage or chemical incompatibility. Manufactured from corrosion resistant AISI Type 316L (1.4404) stainless steel, and hygienically constructed with all welded seams. The bottle has a smooth inner surface finish that corresponds to IIIc (DIN 17441) with rounded inner edges for ease of cleaning. Polished and brushed durable exterior finish. The bottle has a GL 45 thread with a built-in pouring lip. Supplied without a cap, stainless steel cap is available. It is fully compatible with all the GL 45 caps (see page 33 ff).*

*Typical applications: Hazardous materials laboratory container for storage of liquids, intermediates, or solid product. Storage of precious materials, such as high purity fine chemicals, pharmaceutical or cosmetic products. Storage of light sensitive materials.*

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**DURAN GROUP STAINLESS STEEL LABORATORY BOTTLE**  
**with DIN thread, GL 45**

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
**DURAN GROUP STAINLESS STEEL SHIPPING BOTTLE**  
**with DIN thread GL 45, UN certified**

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*The unbreakable DURAN Group steel bottle is ideal for the storage and safe shipping of dangerous liquid goods such as solvents, and reagents, without the need for additional protective outer packaging<sup>1</sup>. Manufactured from corrosion resistant AISI Type 316L (1.4404) stainless steel, and hygienically constructed with all welded seams. The bottle has a smooth inner surface finish that corresponds to IIIc (DIN 17441) with rounded inner edges for ease of cleaning. Polished and brushed durable exterior finish. The bottle has a GL 45 thread with a built-in pouring lip. Supplied complete with GL 45 Stainless Steel screw cap and PTFE surfaced Platinum-catalyzed silicon liner, and UN certification number.*

*Typical applications: Container for the shipping of hazardous or non-hazardous liquids. Storage or shipping of precious liquids, such as high purity fine chemicals, pharmaceutical or cosmetic products. Storage of light sensitive liquids.*

Best.-Nr. Cat. No.		Inhalt Capacity ml	Verpackungs- einheit/ Stück Pack/Quantity
29 901 60 5	DURAN Group GL 45 (Transport-) Edelstahlfflasche, UN zertifiziert <sup>1</sup> , komplett DURAN Stainless Steel GL 45 Shipping Bottle, UN certified <sup>1</sup> , complete	1 500	1

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<sup>1</sup> Certified to UN standards for the carriage of liquids classified as dangerous goods in UN packing Groups II (medium danger) and III (low danger). It is suitable for liquids with a relative density of 2.0 or less. International regulations are subject to change, it is the user's responsibility for complying with all applicable laws and regulations.

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
Cap is moulded from a dye-free polymer to avoid risk of leaching. Wide usable temperature range from  $-196\text{ °C}$  to  $+260\text{ °C}$ . Very high chemical resistance. Complete with colourless PTFE<sup>1</sup> coated silicone sealing disk for high leak tightness. Registered as Type III Drug Master File in conjunction with Premium bottle (DMF no. 19757). Copy available on request. The matching TpCh260<sup>1</sup> pouring ring permits drip-free pouring is available separately. Replacement PTFE<sup>1</sup> coated silicone sealing disks<sup>2</sup> are available.

Typical applications: Due to its exceptional properties, it is ideal for applications in the pharmaceutical industry, storage of aggressive materials, and demanding sterilisation processes such as hot air sterilisation or depyrogenation.

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PREMIUM CAP  
from TpCh260<sup>1</sup> (similar to PFA) with  
PTFE<sup>1</sup> coated silicone seal<sup>2</sup>

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
The tamper-evident screw cap is available with lip seal or PTFE<sup>1</sup> coated silicone seal<sup>2</sup>. The integral plastic ring tears when the cap is first opened and is retained on the bottle neck. Therefore it can be clearly seen whether the bottle retains its original seal. After initial removal, can be used as “normal” screw cap.

Typical applications: secure storage and transport/shipping of valuable media.

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TAMPER-EVIDENT SCREW CAP  
from PP<sup>1</sup>, for laboratory bottles with DIN thread

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
**MEMBRANE SCREW CAP**  
 from PP<sup>1</sup>, blue with welded-in  
 PTFE<sup>1</sup> membrane for pressure equalisation

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For all GL 45, GL 32 and GL 25 threads.  
 Ideal for autoclaving processes because  
 the membrane permits pressure equalisation  
 and tight sealing, greatly reducing the risk of  
 contamination. Ingress of liquids or solids is  
 prevented and the bottle contents remain  
 sterile.

Typical applications: storage or transport  
 of gas generating media, autoclaving of  
 media.

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
Available in the colours blue, green, yellow  
 and grey with matching colour pouring  
 rings. Distinguishing media types is simplified  
 and interchanging of screw caps and carry  
 over of substances is practically eliminated.

Typical applications: visual aid to the safe  
 identification of different media.

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*High leak tightness through use of PTFE<sup>1</sup> coated silicone seal<sup>3</sup>. More chemically resistant than PP<sup>1</sup> screw cap. A matching ETFE<sup>1</sup> pouring ring is also available, permitting clean, drop-free use.*

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
<sup>2</sup> For screw connection caps with aperture, silicone seals with aperture and septa, see pages 179–181.

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
**CONNECTION SYSTEM**

screw cap GL 45 with two or three ports (GL 14 thread)

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Materials used: PP<sup>1</sup> and PTFE<sup>1</sup>. Flexible modular system. Four different tubing diameters (1.6 mm; 3.0 mm; 3.2 mm and 6.0 mm) can be connected. Sterile pressure equalisation is possible through use of syringe filter. Unused ports can be provided with a blind cap.

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Typical applications: safe transfer of liquid media within a closed and sterile system (evaporation is reduced).

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
*GL45 cap connection system for the easy connection of flexible tubing to the DURAN® GL45 bottles. The two polypropylene connectors have angled top and straight underside connectors. The cap is ideal for use with soft elastic tubing that has an inner diameter of 6 – 9 mm (such as Viton®, Tygon®, or other silicone tubing). An optional venting connector is available with, or without a syringe filter.*

*Highly versatile as the screw cap is based on the standard GL 45 thread. Grey polypropylene cap center rotates freely, allowing the bottle to be unscrewed without the need to disconnect the tubing. The angled top connectors prevent kinks in the tubing. Temperature resistant up to +140 °C. Fully autoclavable and washer-proof.*

*Typical applications: Possible biotech applications include the transfer of sterile media from one container to another using a peristaltic pump.*

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

*Self contained stirring system ideal for mixing processes. Stirrer shaft can be adjusted to the optimal position in either a 500 or 1,000 ml DURAN® GL 45 bottles. Drive with standard magnetic stirres up to 500 rpm. Significantly improved mixing performance compared to standard magnetic stirring bars. Multi-connector cap is compatible with the proven DURAN® connection systems; tubing with 1.6 mm to 6.0 mm bores can be used to add or remove liquid or gas. Fully autoclavable. Parts in contact with media conform to FDA requirements. Available separately or with bottle.*

*Typical applications: mixing chemicals, media or cultures in the DURAN® GL 45 laboratory bottle.*

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*DURAN Group Bottle-top dispensers make the serial dispensing of liquids directly from DURAN® laboratory glass bottles safer and quicker.*


*Reliable and reproducible liquid dispensing due to a very precise borosilicate 3.3 glass cylinder. Volume settings can be quickly and accurately adjusted.*

*Simple dispensing operation, even with one hand. Dispenser rotates freely 360° on the bottle to a convenient position. Materials of construction offer good chemical resistance. Adaptable to different bottles threads (adaptors for GL 32, GL 38, S 40 or GL 28, GL 45, S 40) and bottle heights with a telescopic intake tube. PFA sealing lip prevents crystallization build-up. The modular design is easy to dismantle and clean.*

*Supplied with certificate of performance, and multi-lingual instruction manual.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread GL	Adapter Adaptor	Verpackungseinheit/ Stück Pack/Quantity
29 115 03	0,25 – 2,5	32	GL 28, GL 45, S40	1
29 115 06	0,5 – 5,0	32	GL 28, GL 45, S40	1
29 115 09	1,0 – 10,0	32	GL 28, GL 45, S40	1
29 115 14	2,5 – 25,0	45	GL 32, GL 38, S40	1
29 115 17	5,0 – 50,0	45	GL 32, GL 38, S40	1
29 115 24	10,0 – 100,0	45	GL 32, GL 38, S40	1

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**DG SAFETY CAP**  
**screw cap GL 45 with one, two or three ports**

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*High-quality and flexible connection system, specifically designed for handling aggressive, gas-generating media. Materials used: PP<sup>1</sup>, PBT<sup>1</sup> and PTFE<sup>1</sup>. All parts coming into contact with media are from PTFE<sup>1</sup>. Three different hose diameters (1.6 mm; 2.3 mm; 3.2 mm) can be connected. Pressure equalisation is possible through use of membrane filter. Outgassing of media is reduced by use of a one-way valve that decreases media loss and protects the user from vapours. A free moving PTFE<sup>1</sup> core prevents rotation of the hoses when screwing onto the DURAN<sup>®</sup> bottle. Unused ports can be provided with a blind cap. Optionally available with stopcock to prevent hoses running dry or air intakes during bottle changeover. DG safety caps are also available for bottles with standard ground joint.*

*Typical applications: safe transfer of aggressive media within a closed system (evaporation is reduced).*

Best.-Nr. Cat. No.	Bezeichnung Description	DIN-Gewinde DIN Thread GL	Verpackungseinheit/ Stück Pack/Quantity
11 684 01 <sup>2</sup>	DG Safety Cap I (1 Anschluss) DG Safety Cap I (1 port)	45	1
11 684 02 <sup>2</sup>	DG Safety Cap II (2 Anschlüsse) DG Safety Cap II (2 ports)	45	1
11 684 13 <sup>2</sup>	DG Safety Cap III (3 Anschlüsse) DG Safety Cap III (3 ports)	45	1
11 684 14 <sup>2</sup>	DG Safety Cap I (1 Anschluss), mit Absperrhahn DG Safety Cap I (1 port), with stopcock	45	1
11 684 15 <sup>2</sup>	DG Safety Cap II (2 Anschlüsse), mit Absperrhähnen DG Safety Cap II (2 ports), with stopcocks	45	1
11 684 16 <sup>2</sup>	DG Safety Cap III (3 Anschlüsse), mit Absperrhähnen DG Safety Cap III (3 ports), with stopcocks	45	1
11 684 17 <sup>3</sup>	DG Safety Cap für NS 29/32 DG Safety Cap for NS 29/32		1

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<sup>2</sup> Complete with GL 45 screw cap, connection screws, conical nipple (ID 3.2 mm, PTFE), ventilation valve and membrane filter.

<sup>3</sup> Complete with screw cap for NS 29/32, connection screws, conical nipple (ID 3.2 mm, PTFE), ventilation valve and membrane filter.

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
*The activated carbon filter prevents the escaping of hazardous gases and thus increases safety in the laboratory. Materials used: PBT<sup>1</sup> and PTFE<sup>1</sup>. Optional special connection for installations with a leakage system. A free moving PTFE<sup>1</sup> core prevents rotation of the hoses when screwing onto the DURAN® bottle. Unused ports can be provided with a blind cap.*

*Typical applications: safe storage and collection of liquid laboratory waste.*

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**DG SAFETY WASTE CAP**  
 additional connection option for activated carbon filter and leak hose, GL 45

Best.-Nr. Cat. No.	Bezeichnung Description	DIN-Gewinde DIN Thread GL	Verpackungseinheit/ Stück Pack/Quantity
11 684 18 <sup>2</sup>	DG Safety Waste Cap DG Safety Waste Cap	45	1
11 684 20 <sup>3</sup>	DG Safety Waste Cap für Schlauchverbindung DG Safety Waste Cap for hose connector	45	1

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<sup>2</sup> Complete with GL 45 screw cap, connection screws and conical nipple (ID 3.2 mm, PTFE).

<sup>3</sup> Complete with GL 45 screw cap, connection screws, conical nipple (ID 3.2 mm, PTFE) and connection for leak hose.

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
*Matching cap, see pages 44–45.  
A maximum usage temperature of +100 °C is recommended. Thermal shock resistance 30 K. Hydrolytic class 3.*

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SCREW CAP BOTTLE, SQUARE,  
FROM SODA-LIME GLASS  
narrow neck with thread

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
*Without pouring ring and dust cover cap<sup>1</sup>.  
A maximum usage temperature of +100 °C is recommended. Thermal shock resistance 30 K. Hydrolytic class 3.*

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SCREW CAP BOTTLE, SQUARE  
FROM SODA-LIME GLASS  
wide neck with thread

Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Klartglas   <i>Clear glass</i>					
23 820 17	50	32 <sup>2</sup>	48	70	10
23 820 24	100	32 <sup>2</sup>	49	111	10
23 820 36	250	45 <sup>2</sup>	64	146	10
23 820 44	500	54 <sup>2</sup>	76	173	10
23 820 54	1 000	60 <sup>2</sup>	97	213	10
Braunglas   <i>Amber glass</i>					
23 826 24	100	32 <sup>2</sup>	49	111	10
23 826 36	250	45 <sup>2</sup>	64	146	10
23 826 44	500	54 <sup>2</sup>	76	173	10
23 826 54	1 000	54 <sup>2</sup>	97	213	10

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**SCREW CAP BOTTLE, ROUND, AMBER**  
**FROM SODA-LIME GLASS**  
 with thread

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Complete with pouring ring and dust cover  
 cap. A maximum usage temperature of  
 +100 °C is recommended. Thermal shock  
 resistance 30 K. Hydrolytic class 3.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Braunglas   Amber glass					
23 835 66 5	2 500	45 <sup>1</sup>	139	283	1

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**TAMPER-EVIDENT SCREW CAP**  
 from PP<sup>1</sup>, blue, for soda-lime square screw  
 cap bottle, narrow neck

Best.-Nr. Cat. No.	Gewinde Thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Originalitätsverschluss   Tamper evident screw cap				
29 301 19	32 <sup>2</sup>	45	32	10
29 301 28	45 <sup>2</sup>	60	35	10
Ausgießring   Pouring ring				
29 251 19	32 <sup>2</sup>	45	32	10
29 251 28	45 <sup>2</sup>	60	35	10

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**SECURITY SCREW CAP**  
 from PP<sup>1</sup> with valve<sup>2</sup>, red, for soda-lime  
 square screw cap bottles, narrow neck

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

**TAMPER-EVIDENT SCREW CAP**  
 from PP<sup>1</sup>, blue, for soda-lime  
 square screw cap bottles, wide neck

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*With precision ground neck. All glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts.*

*Typical application: storage of powders.*

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**DURAN® REAGENT BOTTLE, WIDE NECK**

neck with standard ground joint,  
with standard ground glass flat-head stopper

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*With precision ground neck. All glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts. Also available in amber. Unchanged DURAN® properties within the bottle, as colouration is only to the outer surface. Very uniform amber colour thanks to innovative technology, moreover long durability and high chemical resistance of the amber colouring.*

*Typical applications: storage of powders.*

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<sup>1</sup> Narrow neck bottles with standard ground joint, 1 000 ml, also suitable for use as burette bottles, see page 124.

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*With precision ground neck. All glass components, therefore also suitable for storage of aggressive media, which could attack plastic parts. Also available in amber. Unchanged DURAN® properties within the bottle, as colouration is only on the outer surface. Very uniform, durable and chemically resistant amber colour due to use of innovative technology.*

*Typical application: storage of liquids.*

Best.-Nr. Cat. No. klar/clear	Best.-Nr. Cat. No. braun/amber	Inhalt Capacity ml	d mm	h mm	Hals Neck NS	Verpackungseinheit/ Stück Pack/Quantity
Hals mit Normschliff und NS-Glasstopfen Neck with standard ground joint and glass flat-head stopper						
21 165 08 <sup>1</sup>		10	28	52	10/19	10
21 165 14	21 168 14	25	36	64	12/21	10
21 165 17	21 168 17	50	42	80	14/15	10
21 165 24	21 168 24	100	52	96	14/15	10
21 165 36	21 168 36	250	70	130	19/26	10
21 165 44	21 168 44	500	86	164	24/29	10
21 165 54	21 168 54	1 000	107	200	29/32	10
21 165 63	21 168 63	2 000	134	248	29/32	10
21 165 73	21 168 73	5 000	182	323	45/40	1
21 165 86	21 168 86	10 000	227	398	60/46	1
21 165 91	21 168 91	20 000	288	492	60/46	1

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


DURAN® REAGENT BOTTLE, NARROW NECK

neck with standard ground joint, with standard ground glass flat-head stopper

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REAGENT BOTTLE, WIDE NECK  
 FROM SODA-LIME GLASS  
 neck with standard ground joint

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Without stopper<sup>1</sup>. A maximum usage temperature of +100 °C is recommended. Thermal shock resistance 30 K. Hydrolytic class 3.

Best.-Nr. Cat. No. klar/clear	Best.-Nr. Cat. No. braun/amber	Inhalt Capacity ml	d mm	h mm	Hals Neck NS	Verpackungseinheit/ Stück Pack/Quantity
Hals mit Normschliff, ohne Stopfen Neck with standard ground joint, without stoppers						
23 184 17		50	44	79	24/20	10
23 184 24	23 187 24	100	52	97	29/22	10
23 184 36	23 187 36	250	71	129	34/24	10
23 184 44	23 187 44 <sup>2</sup>	500	86	164	45/40	10
23 184 54	23 187 54	1 000	107	200	60/46	10

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
**Note** on using DURAN® glass stoppers with soda-lime glass bottles: If the bottle and the stopper have a temperature difference greater than 30 °C, the glass stoppers can become stuck!

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REAGENT BOTTLE, WIDE NECK  
 FROM SODA-LIME GLASS  
 neck with standard ground joint,  
 with standard ground glass flat-head stopper<sup>1</sup>

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A maximum usage temperature of +100 °C is recommended. Thermal shock resistance 30 K. Hydrolytic class 3.

Best.-Nr. Cat. No. klar/clear	Best.-Nr. Cat. No. braun/amber	Inhalt Capacity ml	d mm	h mm	Hals Neck NS	Verpackungseinheit/ Stück Pack/Quantity
Hals mit Normschliff und NS-Glasstopfen Neck with standard ground joint and glass flat-head stopper						
23 185 17		50	44	79	24/20	10
23 185 24	23 188 24	100	52	97	29/22	10
23 185 36	23 188 36	250	71	129	34/24	10
23 185 44	23 188 44 <sup>2</sup>	500	86	164	45/40	10
23 185 54	23 188 54	1 000	107	200	60/46	10

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REAGENT BOTTLE, NARROW NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint


Best.-Nr. Cat. No.	Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Hals Neck NS	Verpackungseinheit/ Stück Pack/Quantity
klar/clear	braun/amber					
Hals mit Normschliff Neck with standard ground joint						
	23 167 17	50	42	80	14/15	10
23 164 24	23 167 24	100	52	96	14/15	10
23 164 36	23 167 36	250	72	130	19/26	10
23 164 44	23 167 44	500	89	165	24/29	10
23 164 54 <sup>1</sup>	23 167 54 <sup>1</sup>	1 000	110	200	29/32	10

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
REAGENT BOTTLE, NARROW NECK  
FROM SODA-LIME GLASS  
neck with standard ground joint,  
with standard ground glass flat-head stopper

Best.-Nr. Cat. No.	Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Hals Neck NS	Verpackungseinheit/ Stück Pack/Quantity
klar/clear	braun/amber					
Hals mit Normschliff und NS-Glasstopfen Neck with standard ground joint and glass flat-head stopper						
	23 168 17	50	42	80	14/15	10
23 165 24	23 168 24	100	52	96	14/15	10
23 165 36	23 168 36	250	72	130	19/26	10
23 165 44	23 168 44	500	89	165	24/29	10
23 165 54	23 168 54	1 000	110	200	29/32	10

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*GLASS STOPPER, AMBER*  
*with standard ground joint, octagonal*

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
*PLASTIC STOPPER*  
*from polyethylene<sup>1</sup>*

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*GLASS STOPPER*  
*with standard ground joint, octagonal*

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

**DURAN® GLASS STOPPER**

ground conical, for reagent bottles  
(oxygen bottles according to Winkler)

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
**Note** on using DURAN® glass stoppers with **soda-lime** glass bottles: If the bottle and the stopper have a temperature difference greater than 30 °C, the glass stoppers can become stuck!

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**DURAN® ASPIRATOR (LEVELLING)**  
**BOTTLE**  
*with plain neck and bottom sidearm*

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
*Useful as a delivery or storage container for solutions. Tubular sidearm (Diameter: Outer  $d_1$ : 11 mm, Inner  $d_2$ : 5 mm). Outlet facilitates attachment of flexible tubing.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 708 36	250	73	130	10
24 708 44	500	89	164	10
24 708 54	1 000	111	200	1

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**DURAN® ACID BOTTLE<sup>1</sup>**  
*with standard ground "pennyhead" stopper,*  
*conical shoulders, interchangeable glass cap*

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*In addition to the ground stopper, a glass cap with ground joint is supplied. This provides an improved seal and protection against acid vapours.*

Best.-Nr. Cat. No.	Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	NS	Verpackungseinheit/ Stück Pack/Quantity
21 275 24 klar   clear	21 275 24 6 braun   amber	100	55	145	19/17	10
21 275 36	21 275 36 6	250	75	180	19/17	10
21 275 44	21 275 44 6	500	82	220	24/20	10
21 275 54	21 275 54 6	1 000	109	260	29/32	10

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**DROPPING BOTTLE**  
**FROM SODA-LIME GLASS**  
with interchangeable clear glass pipette<sup>1</sup> standard  
ground joint, complete with rubber teats<sup>2</sup>

Best.-Nr. Cat. No.	Best.-Nr. Cat. No.	Inhalt Capacity ml	h mm	NS	Verpackungseinheit/ Stück Pack/Quantity
klar   clear	braun   amber				
23 270 17	23 270 17 6	50	79	14/15	10
23 270 24	23 270 24 6 <sup>3</sup>	100	105	14/15	10

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

*Dosing of liquids is possible via an outlet.*

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**DURAN® ASPIRATOR (LEVELLING)**  
**BOTTLE**  
tubulated with standard ground joint,  
without stoppers, neck unground


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Bodentubus Tubulature NS	Verpackungseinheit/ Stück Pack/Quantity
24 701 44	500	86	164	19/26	10
24 701 54	1 000	107	200	19/26	10
24 701 63	2 000	134	249	19/26	10
24 701 73	5 000	182	320	29/32	1
24 701 86	10 000	228	398	29/32	1
24 701 91	20 000	289	492	29/32	1

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**DURAN® ASPIRATOR (LEVELLING)**  
**BOTTLE**

*tubulated with standard ground joint,  
 complete with standard ground stopcock  
 and standard ground stopper*

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*Dosing of liquids is possible via a stopcock.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Bodentubus Tubulature NS	Hals Neck NS	Orig. Pack./Stück Stand. pack/quantity
24 702 44	500	86	164	19/26	24/29	10
24 702 54	1 000	107	200	19/26	29/32	10
24 702 63	2 000	134	249	19/26	29/32	10
24 702 73	5 000	182	320	29/32	45/40	1
24 702 86	10 000	228	398	29/32	60/46	1
24 702 91	20 000	289	492	29/32	60/46	1

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**STOPCOCK WITH STANDARD**  
**GROUND JOINT FOR ASPIRATOR**  
**BOTTLE**

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
**DURAN® ASPIRATOR (LEVELLING)**  
**BOTTLE**

*neck with DIN thread GL 45,  
 tubulator with GL 32*

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*Complete with screw connection cap<sup>1</sup>,  
 silicone seal and stopcock with PTFE<sup>2</sup>  
 spindle. Dosing of liquids is possible via a  
 stopcock.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 703 54	1 000	45	101	225	1
24 703 63	2 000	45	136	260	1
24 703 73	5 000	45	182	330	1
24 703 86	10 000	45	230	410	1

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<sup>2</sup> For chemical and thermal resistance, see page 222.


Best.-Nr. Cat. No.	Bohrung Hole mm	für Inhalt for Capacity ml	Orig. Pack./Stück Stand. pack/quantity
24 147 03	6	1 000 + 2 000	1
24 147 04	8	5 000 + 10 000	1

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**STOPCOCK FOR ASPIRATOR BOTTLE**  
with PTFE<sup>1</sup> spindle, for GL 32 screw thread

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
**PTFE-SILICONE<sup>1</sup> SEAL**  
for GL 32 screw thread

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“LONDON, BOSTON, SHANGHAI ... IN REALITY,  
MASTERMINDS ARE VALUED ANYWHERE IN THE  
WORLD ... AS IS **GLASSWARE MADE IN GERMANY.**”

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## 2 | BOILING FLASKS AND GENERAL LABORATORY GLASSWARE

*DURAN® laboratory glassware, including heating vessels, has very good thermal-shock resistance ( $\Delta T=100\text{ K}$ ) and a high operating temperature ( $+500\text{ °C}$ ). Not only the glass type but also its uniform wall thickness distribution are critical in preventing uneven expansion and stressing of the glass which could result in failure. For this reason, wall thickness distribution is, as a vital quality characteristic, continuously checked during the production process.*

*The beakers are primarily used as heating vessels. The tall shape is particularly suited to heating in liquid baths where the beaker contents are protected against the surrounding medium.*

*Erlenmeyer flasks are well suitable for mixing, because of their conical shape.*

*Weighing bottles are used when accurately weighing out substances. Close fitting lids with moulded grips are used to prevent the substances from being lost, e.g. during transport within the laboratory.*

*Watch glass dishes can be used both for covering beakers and Erlenmeyer flasks as well as for weighing small quantities of substances.*

*Our product range also includes a wide range of test tubes. In addition to DURAN® glass, other glass types are available (FIOLAX®, soda-lime). The characteristics of each glass type may be found on pages 217-221, 234.*

### *Usage tips:*

- *Due to the uniform wall thickness distribution suitable for very high temperature changes.*
- *The printed scale on many items of DURAN® laboratory glassware is indicative with an accuracy of  $\pm 10\%$ . Therefore the items are not suitable for use as volumetric glassware.*
- *The products are not designed for use under differential pressure or vacuum conditions.*

*DURAN® beakers and Erlenmeyer flasks are provided with a retrace code. Using the eight-character code and the corresponding article number, a batch and quality certificate can be obtained at [www.duran-group.com](http://www.duran-group.com).*

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*With easy-to-read scale and large labeling field for easy marking in fired-on, highly durable, white ceramic. Spout for clean pouring. Uniform wall thickness distribution makes these beakers ideal for heating applications.*

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
*Better mechanical stability due to increased glass content. Reinforced rim increases shock resistance and reduces risk of breakage.*

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
DURAN® SUPER DUTY BEAKER  
low form with spout

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*Typical applications: Working under mechanical load.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 107 29 09	150	60	80	10
21 107 36 05	250	70	95	10
21 107 41 04	400	80	110	10
21 107 48 07	600	90	125	10
21 107 54 09	1 000	105	145	10
21 107 63 02	2 000	132	185	10
21 107 73 07	5 000	170	270	1

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
**Note**

*Uniform and slow heating is recommended. The standard DURAN® beakers should be used when working at very high temperatures or if rapid temperature changes are expected.*

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### DURAN® BECHER dickwandig (Filtrierbecher)

### DURAN® BEAKER heavy-wall (filtering beaker)

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*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Has, due to the increased wall thickness, better mechanical properties than the standard beaker. Thermal shock resistance, however, is reduced so only limited application for heating. Spout for clean pouring.*


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 131 24	100	52	85	10
21 131 29	150	54	93	10
21 131 36	250	70	94	10
21 131 44	500	89	124	10
21 131 54	1 000	105	160	10
21 131 63	2 000	135	195	10
21 131 68	3 000	157	205	4
21 131 73	5 000	182	256	1
21 131 86 <sup>1</sup>	10 000	225	340	1
21 131 88 <sup>1</sup>	15 000	260	390	1
21 131 91 <sup>1</sup>	20 000	285	430	1

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*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Spout for clean pouring. Uniform wall thickness distribution makes these beakers ideal for heating applications.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 116 17	50	38	70	10
21 116 24	100	48	80	10
21 116 29	150	54	95	10
21 116 36	250	60	120	10
21 116 41	400	70	130	10
21 116 48	600	80	150	10
21 116 53	800	90	175	10
21 116 54	1 000	95	180	10
21 116 63	2 000	120	240	10
21 116 68	3 000	135	280	2

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
*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Uniform wall thickness distribution makes these beakers ideal for heating applications.*

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**DURAN® BEAKER**  
**high form, without spout**

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 117 17	50	38	70	10
21 117 24	100	48	80	10
21 117 29	150	54	95	10
21 117 36	250	60	120	10
21 117 41	400	70	130	10
21 117 48	600	80	150	10
21 117 54	1 000	95	180	10

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

**DURAN® BERZELIUS BEAKER**  
**without spout**

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
*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Due to conical form, suited to the mixing of liquids. Uniform wall thickness distribution makes these flasks ideal for heating applications.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d	d <sub>i</sub>	h	Verpackungseinheit/ Stück Pack/Quantity
21 216 14	25 <sup>2</sup>	42	22	75	10
21 216 17	50	51	22	90	10
21 216 24	100	64	22	105	10
21 216 28	125	67	28	112	10
21 990 27 <sup>1</sup>	150	74	28	118	10
21 216 32 <sup>1</sup>	200	79	34	131	10
21 216 36	250	85	34	145	10
21 216 39 <sup>1</sup>	300	87	34	156	10
21 216 44	500	105	34	180	10
21 216 54	1 000	131	42	220	10
21 216 63	2 000	166	50	280	10
21 216 68	3 000	187	52	310	2
21 216 73	5 000	220	52	365	1

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
#### Note

*Uniform and slow heating is recommended. The standard DURAN® Erlenmeyer flasks should be used when working at very high temperatures or if rapid temperature changes are expected.*

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

*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. Due to conical form, suited to the mixing of liquids. Uniform wall thickness distribution makes these flasks ideal for heating applications. The wide neck enables easy filling and cleaning.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 226 14 <sup>1</sup>	25 <sup>2</sup>	43	31	70	10
21 226 17	50	51	34	85	10
21 226 24	100	64	34	105	10
21 226 32 <sup>1</sup>	200	79	50	131	10
21 226 36	250	85	50	140	10
21 226 39 <sup>1</sup>	300	87	50	156	10
21 226 44	500	105	50	175	10
21 226 54	1 000	131	50	220	10
21 226 63 <sup>1</sup>	2 000	153	72	276	10

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
*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. The flask can be closed with a PBT cap<sup>1</sup> or membrane cap<sup>2</sup> (permits gas exchange).*

*Typical applications: The flask is suitable for storage, media preparation and cultivation.*

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
*DURAN® IODINE FLASK*  
*Erlenmeyer shape, with standard ground joint*  
*and glass stopper*

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*With easy-to-read scale and large labelling field for easy marking in fired-on, highly durable, white ceramic. The flask can be closed with a glass stopper.*


*Typical applications: the iodine flask is suitable for determining the iodine number, i.e. the content of unsaturated fatty acids in oils and fats.*

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*DURAN® CONICAL FLASK*  
*Erlenmeyer shape, wide neck*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 227 68 <sup>1</sup>	3 000	190	106	285	1
21 227 73 <sup>1</sup>	5 000	220	108	322	1
21 227 86 <sup>1</sup>	10 000	285	147	420	1

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*Uniform wall thickness distribution makes these flasks ideal for heating applications. The geometry permits very uniform heating. Flasks with a neck diameter of 65 mm or more have a reinforced rim.*

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DURAN® ROUND BOTTOM FLASK  
narrow neck, with beaded rim


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 721 17	50	51	26	95	10
21 721 24	100	64	26	110	10
21 721 36	250	85	34	144	10
21 721 44	500	105	34	168	10
21 721 54	1 000	131	42	200	10
21 721 64 <sup>1</sup>	2 000	166	42	260	10
21 721 68 <sup>1</sup>	3 000	185	50	260	1
21 721 71	4 000	207	52	290	1
21 721 73 <sup>1</sup>	5 000	223	50	305	1
21 721 77 <sup>1</sup>	6 000	236	51	355	1
21 721 86	10 000	279	65	380	1
21 721 87 <sup>1,2</sup>	12 000	295	65	380	1
21 721 91 <sup>2</sup>	20 000	345	76	515	1

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*Uniform wall thickness distribution makes these flasks ideal for heating applications. The geometry permits very uniform heating. The wide neck permits easy filling and removal of flask contents. Flasks with a neck diameter of 76 mm or more have a reinforced rim.*

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
DURAN® ROUND BOTTOM FLASK  
wide neck, with beaded rim

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 741 17 <sup>1</sup>	50	51	34	105	10
21 741 24	100	64	35	110	10
21 741 36	250	85	51	143	10
21 741 44	500	105	50	168	10
21 741 54	1 000	131	50	200	10
21 741 55 <sup>1</sup>	1 000	131	65	200	10
21 741 63	2 000	165	76	240	10
21 741 64 <sup>1</sup>	2 000	166	50	240	10
21 741 68 <sup>1</sup>	3 000	185	65	260	1
21 741 71	4 000	206	76	290	1
21 741 73 <sup>1</sup>	5 000	223	65	310	1
21 741 76	6 000	236	89	330	1
21 741 77 <sup>1</sup>	6 000	236	65	330	1
21 741 86 <sup>1</sup>	10 000	279	89	420	1
21 741 91 <sup>1</sup>	20 000	345	89	520	1

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
**DURAN® FLAT BOTTOM FLASK**  
**narrow neck, with beaded rim**

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*Uniform wall thickness distribution makes these flasks ideal for heating applications. Flat base means flasks can be set down without a supporting ring. Flasks with a neck diameter of 65 mm have a reinforced rim.*


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 711 17	50	51	26	90	10
21 711 24	100	64	26	105	10
21 711 36	250	85	34	138	10
21 711 44	500	105	34	163	10
21 711 54	1 000	131	42	190	10
21 711 64 <sup>1</sup>	2 000	166	42	250	10
21 711 68 <sup>1</sup>	3 000	185	50	250	1
21 711 71	4 000	207	50	275	1
21 711 73 <sup>1</sup>	5 000	223	50	290	1
21 711 76	6 000	237	65	315	1
21 711 86	10 000	280	65	360	1

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**DURAN® FLAT BOTTOM FLASK**  
**wide neck, with beaded rim**

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*Uniform wall thickness distribution makes these flasks ideal for heating applications. Flat base means flasks can be set down without a supporting ring. The wide neck permits easy filling and removal of flask contents. Flasks with a neck diameter of 76 mm have a reinforced rim.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 731 17	50	51	34	90	10
21 731 24	100	64	34	105	10
21 731 36	250	85	50	138	10
21 731 44	500	105	50	163	10
21 731 54	1 000	131	50	190	10
21 731 63 <sup>1</sup>	2 000	166	76	230	10
21 731 64	2 000	166	50	230	10

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
*Uniform wall thickness distribution makes these flasks ideal for heating applications and distillations.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>i</sub> mm	h mm	Seitenrohr Side arm		Verpackungseinheit/ Stück Pack/Quantity
					d <sub>2</sub> mm	l mm	
21 653 24	100	66	20	215	6	100	10
21 653 28	125	68	22	215	7	100	10
21 653 29	150	73	20	223	6	100	10

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*DURAN® ENGLER FLASK*  
*with side outlet*

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
*DURAN® ENGLER FLASK*  
*with standard ground joint 19/26, side outlet*


**GLASKERAMIK LABOR-SCHUTZPLATTE**

**GLASS CERAMIC LABORATORY PROTECTION PLATE**

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*Due to low thermal expansion stresses, these glass ceramic plates are well suited to heating glassware with a Bunsen burner.*

Best.-Nr. Cat. No.	Abmessungen der Platte Plate dimensions a x b mm	Verpackungseinheit/ Stück Pack/Quantity
23 821 53	135 x 135	10
23 821 57	155 x 155	10
23 821 58	175 x 175	10

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**SQUARE QUADRUPOD  
for glass ceramic laboratory protection plate**

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*Plate holder for glass ceramic plates. Made from heat-resistant chrome-nickel steel with four legs for extra stability.*


Best.-Nr. Cat. No.	Für Platte For plates mm	Abmessungen der Platte Plate dimensions a x b mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
29 077 53	23 821 53	135 x 135	210	5
29 077 57	23 821 57	155 x 155	210	5
29 077 58	23 821 58	175 x 175	220	5

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**PLATE HOLDER  
for glass ceramic laboratory protection plate**

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*Plate holder for glass ceramic plates. Made from heat-resistant chrome-nickel steel.*


Best.-Nr. Cat. No.	Für Platte For plates mm	Abmessungen der Platte Plate dimensions a x b mm	Verpackungseinheit/ Stück Pack/Quantity
29 078 53	23 821 53	135 x 135	5
29 078 57	23 821 57	155 x 155	5
29 078 58	23 821 58	175 x 175	5

Ideal zur Reinigung von Glaskeramikplatten. *Ideal for cleaning glass ceramic plates.*

Best.-Nr. Cat. No.	Verpackungseinheit/Stück Pack/Quantity
29 079 01	10

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
CLEANING SCRAPER  
for glass ceramic laboratory protection plate

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SPIRIT LAMP  
FROM SODA-LIME GLASS  
without filling tubulature, with ground over-cap


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 301 32 <sup>1</sup>	15	50	25	10
21 301 34 <sup>1</sup>	45	60	30	10
21 301 38 <sup>2</sup>	60	70	35	10
21 301 41 <sup>2</sup>	90	80	45	10
21 301 44	170	95	55	10
21 301 49	320	115	65	10
21 301 54	600	140	80	10
21 301 59	1 500	190	100	10
21 301 63	2 500	230	130	10

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<sup>1</sup> Without print.  
<sup>2</sup> Without labelling field.

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DURAN® EVAPORATING DISHES  
with spout

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
**DURAN® CRYSTALLIZING DISH**  
*with and without spout*


Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
Mit Ausguss DIN 12 338 With spout DIN 12 338				
21 311 24 <sup>1</sup>	20	40	25	10
21 311 32 <sup>1</sup>	40	50	30	10
21 311 34 <sup>1</sup>	60	60	35	10
21 311 38	100	70	40	10
21 311 41	150	80	45	10
21 311 44	300	95	55	10
21 311 49	500	115	65	10
21 311 54	900	140	75	10
21 311 59	2 000	190	90	10
21 311 63	3 500	230	100	10
Ohne Ausguss DIN 12 337 Without spout DIN 12 337				
21 313 24 <sup>1</sup>	20	40	25	10
21 313 32 <sup>1</sup>	40	50	30	10
21 313 34 <sup>1</sup>	60	60	35	10
21 313 38	100	70	40	10
21 313 41	150	80	45	10
21 313 44	300	95	55	10
21 313 49	500	115	65	10
21 313 54	900	140	75	10
21 313 59	2 000	190	90	10
21 313 63	3 500	230	100	10

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**DURAN® WEIGHING BOTTLE**  
*with ground lid*

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*Close-fitting lid prevents any sample loss during transport after weighing. Available in low and high forms.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d x h mm	Verpackungseinheit/ Stück Pack/Quantity
Niedrige Form Flat shape			
24 210 13	5	28 x 25	10
24 210 23	15	38 x 30	10
24 210 32	30	54 x 30	10
24 210 41	80	85 x 30	10
Hohe Form Tall shape			
24 211 13	10	28 x 40	10
24 211 18	20	32 x 50	10
24 211 23	45	38 x 70	10
24 211 24	70	44 x 80	10

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
Available in DURAN® and also in *soda-lime glass* (for properties, see pages 217–222, 234).

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
DURAN® ORGANGLAS  
ohne Stopfen

DURAN® ORGAN STORAGE JAR  
without stopper

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
*DURAN® CYLINDER*  
 with knobbed lid, polished rim


Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Verpackungseinheit/ Stück Pack/Quantity
24 205 01	80 x 80	250	10
24 205 03	100 x 100	500	10
24 205 05	120 x 120	1 000	1
24 205 10	150 x 150	2 000	1
24 205 21	210 x 210	6 000	1
24 205 32	260 x 260	12 000	1

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	Schliffgröße Grinding height mm	h mm	h <sub>1</sub> mm	Verpackungseinheit/ Stück Pack/Quantity
21 580 24	100	52	34,5	135	39	10
21 580 39	300	69	45	163	48	10
21 580 48	600	81	50	214	50	10
21 580 51	750	90	60	240	57	10
21 580 58	1 200	100	60	253	57	10

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
*The precise grinding of the knobbed lid and base vessel enables a very tight seal.*

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**DURAN® SPECIMEN JAR  
with ground, knobbed lid**


Best.-Nr. Cat. No.	d x h mm	Volumen ca. Volume approx. ml	Verpackungseinheit/ Stück Pack/Quantity
24 209 02	65 x 63	125	10
24 209 09	65 x 103	250	10
24 209 11	110 x 103	800	10
24 209 16	85 x 123	600	10
24 209 17	127 x 123	1 400	1
24 209 24	85 x 153	750	1
24 209 26	110 x 153	1 250	1
24 209 28	157 x 153	2 500	1
24 209 38	110 x 203	1 750	1
24 209 39	157 x 203	4 000	1
24 209 49	110 x 253	2 000	1
24 209 50	157 x 253	4 500	1
24 209 57	120 x 300	2 250	1
24 209 59	252 x 303	13 000	1

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*The precise grinding of the glass plate and base vessel enables a very tight seal.*

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
**DURAN® MUSEUM JAR  
with ground glass plate**

Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/ Stück Pack/Quantity
21 363 05	60 x 50 x 100	10
21 363 11	100 x 50 x 120	1
21 363 13	130 x 50 x 130	1
21 363 19	150 x 50 x 150	1
21 363 28	120 x 60 x 180	1
21 363 47	210 x 100 x 210	1
21 363 58	250 x 140 x 250	1

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**DURAN® MULTI-PURPOSE CYLINDER**  
**with round base**

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
**DURAN® STANDING CYLINDER**  
**with round base**

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*Wall thickness and geometry designed  
 to suit vacuum applications.*

Best.-Nr. Cat. No.	h + h <sub>1</sub> mm	h x d mm	Verpackungseinheit/ Stück Pack/Quantity
24 460 59	300	250 x 185	1
24 460 66	305	255 x 260	1
24 460 69	350	300 x 315	1

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
Wall thickness and geometry designed to suit vacuum applications. Neck aperture, standard ground joint 34/35.

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**DURAN® BELL JAR with aperture in neck**

Best.-Nr. Cat. No.	h + h <sub>1</sub> mm	h x d mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 465 59	300	250 x 185	34/35	1
24 465 61	350	300 x 215	34/35	1
24 465 69	550	500 x 315	34/35	1

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
For scale divisions and accuracy limits, see table:

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**DURAN® SEDIMENTATION CONE Imhoff type, graduated**

Best.-Nr. Cat. No.	Inhalt Capacity ml	Größter d Max. d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 401 54	1 000	120	470 max.	10

Graduierung Scale ml	Teilung Division ml	Toleranz +/- Tolerance +/- ml
0-2	0,1	0,1
2-10	0,5	0,5
10-40	1	1
40-100	2	2
1 000	Ringmarke / Circular marking	10

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
#### DURAN® TEST TUBE with beaded rim or straight rim

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
The test tubes are thick-walled and therefore mechanically very resistant, yet still retain good thermal shock resistance.

Best.-Nr. Bördelrand Cat. No. Beaded rim	Best.-Nr. Rand gerade Cat. No. Straight rim	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
26 130 01	26 131 01	8 x 70	2	0,8-1,0	100
26 130 03	26 131 03	10 x 75	4	0,8-1,0	100
26 130 06	26 131 06	10 x 100	5	0,8-1,0	100
26 130 08	26 131 08	12 x 75	6	0,8-1,0	100
26 130 11	26 131 11	12 x 100	8	0,8-1,0	100
26 130 12	26 131 12	13 x 100	9	0,8-1,0	100
26 130 13	26 131 13	14 x 130	16	0,8-1,0	100
26 130 16	26 131 16	16 x 130	17	1,0-1,2	100
26 130 21	26 131 21	16 x 160	21	1,0-1,2	100
26 130 23	26 131 23	18 x 180	32	1,0-1,2	100
26 130 26	26 131 26	20 x 150	34	1,0-1,2	100
26 130 28	26 131 28	20 x 180	40	1,0-1,2	100
26 130 33	26 131 33	25 x 150	55	1,0-1,2	50
26 130 36	26 131 36	25 x 200	70	1,0-1,2	50
26 130 38	26 131 38	30 x 200	100	1,0-1,4	50

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#### FIOLAX® BOROSILICATE TEST TUBE with beaded rim

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Thin-walled test tubes suited to rapid temperature changes or localized heating.


Best.-Nr. Cat. No.	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
26 110 01	8 x 70	2	0,4-0,5	100
26 110 03	10 x 75	4	0,4-0,5	100
26 110 06	10 x 100	6	0,4-0,5	100
26 110 08	12 x 75	6,5	0,4-0,5	100
26 110 11	12 x 100	9	0,4-0,5	100
26 110 13	14 x 130	16	0,4-0,5	100
26 110 16	16 x 130	20	0,5-0,6	100
26 110 21	16 x 160	25	0,5-0,6	100
26 110 23	18 x 180	35	0,5-0,6	100
26 110 26	20 x 150	39	0,5-0,6	100
26 110 28	20 x 180	45	0,5-0,6	100
26 110 33	25 x 150	60	0,6-0,7	50
26 110 36	25 x 200	80	0,6-0,7	50
26 110 38	30 x 200	110	0,7-0,8	50

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*NMR tubes are available, according to requirement, in three accuracy classes. The correct tube can be selected depending on magnetic field and spin. These tubes are noteworthy for their close tolerances and accuracy, especially to their straightness, wall thickness and wall thickness distribution. Consequently, quick and accurate test results are achievable.*

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*“IS EVERYTHING MADE CORRECTLY? ... TO START WITH, WE HAD SO MANY QUESTIONS ... BUT AT LEAST ONE THING IS CLEAR: IT'S GOOD TO KNOW WE CAN RELY ON SOMEONE WHO HAS BEEN AROUND FOR 100 YEARS.”*

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### 3 | DESICCATORS

*DURAN® desiccators are used for drying moist substances or as storage vessels for moisture-sensitive products. To accelerate the drying process, the desiccators can be used under vacuum. Due to the high wall-thickness of the vessels and the exact machining of the vacuum-tight ground joints on the lid and base, storage under vacuum is possible - even over extremely long periods.*

*All individual parts and a wide range of accessories such as lids, stopcocks, bases, etc. are compatible and can be interchanged as required. Always ensure the individual parts have the same DN (nominal diameter).*

*For desiccators, the DN is based on the diameter of the sieve plate; this can be measured directly. For lids, measure the outside diameter of the flange and cross-reference with the tables on the product pages.*

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*An overview on page 84 indicates which individual parts are required to assemble the desired desiccator.*

#### Usage tips:

- Designed for use under full vacuum.
- Due to the high wall thickness and the reduced thermal shock resistance under pressure loading, the desiccators must not be heated on one side only or heated using a naked flame.
- Before evacuation, it is recommended that the glass surfaces of the desiccator be checked for damage such as scratches, cracks or nicks. Damaged desiccators must not be used for safety reasons.
- Never expose desiccators to abrupt pressure changes (do not suddenly ventilate evacuated vessels).

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*Vacuum-tight, made from DURAN® borosilicate glass 3.3. To accelerate drying, a vacuum can be applied via the stop-cock. Spare parts such as lids, bases, stopcocks and caps can be interchanged (observe DN).*


*Typical applications: drying of moist samples and storage of moisture-sensitive substances.*

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*DURAN® VACUUM DESICCATOR  
with NOVUS standard ground joint (NS 24/29)  
junction tube in the lid, stopcock and flat flange*

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*Made from DURAN® borosilicate glass 3.3. Components are vacuum tight (DIN 12491). Spare parts such as lids and bases can be interchanged (observe DN).*


*Typical applications: drying of moist products and storage of moisture-sensitive substances.*

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*DURAN® DESICCATOR  
with flat flange and knobbed lid,  
no connection*

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**Note:** The desiccators on this page are no longer available as sets, but can be assembled by selecting components from the accompanying table.

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*DURAN® VACUUM DESICCATOR*  
flange with groove and NOVUS standard ground joint (NS 24/29) junction tube in the lid, with O-ring and stopcock

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*DURAN® VACUUM DESICCATOR*  
base with NOVUS standard ground joint (NS 24/29) and stopcock, flat flange and knobbed lid

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*DURAN® VACUUM DESICCATOR*  
base with flat flange, lid with WERTEX outlet and safety stopper

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*DURAN® VACUUM DESICCATOR*  
with MOBILEX screw thread (GL 32) in lid and base, with stopcock

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
<sup>1</sup> The characters XX in the 6th and 7th position should be replaced, dependent on DN size: DN 100 = 46, DN 150 = 57, DN 200 = 61, DN 250 = 66, DN 300 = 69

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
24 770 46	100	112	58	119	153 +/- 2	0,7	1
24 770 57	150	154	81	172	215 +/- 2	2,4	1
24 770 61	200	202	115	224	270 +/- 2	5,8	1
24 770 66	250	235	120	274	320 +/- 2	10,5	1
24 770 69	300	283	150	332	380 +/- 2	18,5	1

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*DURAN® DESICCATOR BASE  
with flat flange, no outlet,  
suitable for all types of lids*

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
*DURAN® DESICCATOR BASE  
with ring-grooved flange,  
suitable for all types of lids*

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
DURAN® DESICCATOR BASE  
with flat flange, standard ground outlet (24/29),  
type NOVUS, suitable for all types of lids

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
24 771 46	100	112	58	119	153 +/- 2	0,7	1
24 771 57	150	154	81	172	215 +/- 2	2,4	1
24 771 61	200	202	118	224	270 +/- 2	5,8	1
24 771 66	250	235	122	274	320 +/- 2	10,5	1
24 771 69	300	283	154	332	380 +/- 2	18,5	1

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
DURAN® DESICCATOR BASE  
with flat flange, screw thread outlet,  
type MOBILEX (GL 32), suitable for all  
types of lids

Best.-Nr. Cat. No.	DN	h mm	h <sub>1</sub> mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	ca. Volumen approx. Vol. l	Verpackungseinheit/ Stück Pack/Quantity
24 772 57	150	154	81	172	215 +/- 2	2,4	1
24 772 61	200	202	118	224	270 +/- 2	5,8	1
24 772 66	250	235	122	274	320 +/- 2	10,5	1
24 772 69	300	283	154	332	380 +/- 2	18,5	1

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
*DURAN® DESICCATOR LID  
with knob, suitable for all types of bases*

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
*DURAN® DESICCATOR LID  
for standard ground joint stopcocks (NS 24/29),  
type NOVUS, suitable for all types of bases*

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
*DURAN® DESICCATOR LID*  
*with special tube (NS 24/29 type WERTEX),*  
*with ring grooved flange, suitable for all types*  
*of bases*


Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 430 57	150	85	172	215 +/- 2	24/29	1
24 430 66	250	109	274	320 +/- 2	24/29	1

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*DURAN® DESICCATOR LID*  
*with threaded outlet, type MOBILEX (GL 32),*  
*suitable for all bases*

Best.-Nr. Cat. No.	DN	h mm	ID Flansch ID Flange mm	AD Flansch OD Flange mm	GL	Verpackungseinheit/ Stück Pack/Quantity
24 440 57	150	85	172	215 +/- 2	32	1
24 440 61	200	94	224	270 +/- 2	32	1
24 440 66	250	109	274	320 +/- 2	32	1
24 440 69	300	137	332	380 +/- 2	32	1

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


O-RING<sup>1</sup>  
suitable for desiccators, from silicone (VMQ)

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Ordering advice: the O-ring is dimensioned according to the groove diameter.

Best.-Nr. Cat. No.	DN	d mm	s mm	Verpackungseinheit/ Stück Pack/Quantity
29 224 57	150	174	5,7	1
29 224 60	200	229	5,7	1
29 224 66	250	274	6,5	1
für Modelle ab 1996 for articles since 1996				
29 215 57	150	185	5,3	1
29 215 61	200	236	5,3	1
29 215 66	250	290	5,3	1

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DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE  
for desiccator base side outlets, type NOVUS  
(NS 24/29)


Best.-Nr. Cat. No.	l mm	d mm	NS	Verpackungseinheit/ Stück Pack/Quantity
24 798 03 <sup>2</sup>	85	8	24/29	1

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
<sup>2</sup> With self-sealing valve stopcock and PP tightening torque limiter.

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**DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE**  
for desiccator lid outlets, type NOVUS (NS 24/29)

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Matching screw cap and seal, see page 92.

Best.-Nr. Cat. No.	l mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
24 797 03 <sup>2</sup>	160	8	1

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
**DURAN® STOPCOCK WITH  
PTFE<sup>1</sup> SPINDLE**  
for threaded outlets, type MOBILEX (GL 32)

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
**SCREW CAP WITH APERTURE GL 32**

from PBT<sup>1</sup>, red, with pierced aperture for  
MOBILEX type threaded outlets

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Plastic (PBT<sup>1</sup>). Suitable seal: Cat. No.  
29 236 12, see pages 180–181.

Best.-Nr. Cat. No.	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
29 227 08	41	26	10

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
**SCREW CAP GL 32**

from PBT<sup>1</sup>, red, with PTFE<sup>1</sup> coated seal for  
MOBILEX type threaded outlets

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Plastic. Cap for blanking unused outlets  
(type MOBILEX).

Best.-Nr. Cat. No.	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
29 240 19	41	26	10

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
**SAFETY STOPCOCK**  
for safety outlets type Wertex, NS 24/29

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*“GROWTH... DEPENDS ON MANY FACTORS,  
NOT JUST FOR MY CULTURES... IT'S GOOD TO  
KNOW, THAT IN A CHANGING WORLD, DURAN®  
IS ONE **RELIABLE CONSTANT**.”*

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## 4 | GLASSWARE FOR MICROBIOLOGY

*Due to its high thermal-shock resistance, DURAN® microbiology glassware is ideal for autoclaving and sterilisation processes and shows, even after multiple use no signs of wear. Unlike plastic items, it is very resistant to mechanical wear even after repeated use and sterilisation cycles.*

*Due to the nearly inert behaviour, there are no interactions (e.g. ion exchange) between medium and glass and any spurious influence on experiments is thereby effectively excluded.*

*DURAN® products are completely transparent in visible light and unlike many plastic products are ideal for use under the microscope.*

*DUROPLAN® Petri dishes are outstanding due to their distortion-free transparency and high planarity. These*

*excellent geometrical properties enable uniform agar distribution and reproducible culture growth.*

*Alongside the Petri dishes, the DURAN® range includes a wide range of culture bottles, culture flasks, roller bottles and spot plates.*

*In addition, there are various types of staining dishes.*

### Usage tips:

- Only autoclave products, which are free from damage such as scratches, cracks or nicks.
- The outstanding thermal properties (max. operating temperature of +500 °C, thermal shock resistance  $\Delta T = 100$  K) enable high temperature processes, such as hot air sterilisation.

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*These Petri dishes are made from DURAN® borosilicate glass 3.3 using a special manufacturing process, which permits the uniform distribution of agar and guarantee distortion-free viewing.*


*Typical applications: biological and medical work, cultivation of microorganism, microscopy.*

Best.-Nr. Cat. No.	$d_1 \times h_1$ mm	$d_2$ mm	$h$ mm	Verpackungseinheit/ Stück Pack/Quantity
21 755 41	60 × 20	54	22	10
21 755 43	80 × 20	74	22	10
21 755 46	100 × 15	94	17	10
21 755 48	100 × 20	94	22	10
21 755 51	120 × 20	114	22	10
21 755 53	150 × 27	143	32	10

Best.-Nr. Cat. No.	$d \times h$ mm	Verpackungseinheit/Stück Pack/Quantity
23 755 39	40 × 12	10
23 755 40	60 × 15	10
23 755 42	80 × 15	10
11 840 71	90 × 15	10
23 755 45	100 × 10	10
23 755 46	100 × 15	10
23 755 48	100 × 20	10
23 755 51	120 × 20	10
23 755 52	150 × 25	10
23 755 56	180 × 30	10
23 755 59	200 × 30	10
23 755 61	200 × 45	10

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
*DUROPLAN® PETRI DISH*  
base and lid are flat inside and out, and free  
from bubbles and streaks

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
*STERIPLAN® PETRI DISH*  
MADE OF SODA-LIME GLASS

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*SPOT PLATE, TYPE FEIGL*  
**MADE OF SODA-LIME GLASS**

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
*Typical application: detection reactions.*

Best.-Nr. Cat. No.	l mm	b mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
23 671 52	130	100	14	10

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*Bottom of the base has three rests for stability.*

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*Typical applications: biological and medical work, preparation of agars, microscopy.*

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*Centrifuge tubes are thick-walled and very resistant to mechanical loading (For further details needed for calculation of a maximum speed in rpm, see page 234). The higher density fraction collects in the bottom. Consequently solids can be collected and separated.*

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Verpackungseinheit/ Stück Pack/Quantity
21 601 10	12 x 100	6	50
21 601 11	16 x 100	12	50
21 601 14	24 x 100	25	10
21 601 17	34 x 100	50	10
21 601 24 <sup>1</sup>	40 x 115	80	10
21 601 26	44 x 100	80	10
21 601 36 <sup>1</sup>	56 x 147	250	10

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*Centrifuge tubes are thick-walled and very resistant to mechanical loading (For further details needed for calculation of a maximum speed in rpm, see page 234). The higher density fraction collects in the pointed centre of the bottom. Consequently even small amounts of solids can be collected and separated.*

Best.-Nr. Cat. No.	d x h mm	Nenninhalt Nominal capacity ml	Verpackungseinheit/ Stück Pack/Quantity
24 263 09	16 x 100	12	50

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



*The straight rim permits the use of Kapsenberg caps; tubes are therefore well suited to the culture of micro-organisms.*

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
*Typical applications: growth and storage of sterile cultures.*

Best.-Nr. Cat. No.	Best.-Nr für passende Kapsenberg-Kappen Cat. No. for matching Kapsenberg cap	d x h mm	ca. Volumen approx. vol. ml	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
26 132 21	29 010 09	16 x 160	20	1,0-1,2	100
26 132 23	29 010 11	18 x 180	30	1,0-1,2	100

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*DURAN® CULTURE TUBE*  
*with DIN thread and screw cap from PBT<sup>1</sup>*

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*The DIN thread permits the use of PBT<sup>1</sup> screw caps; tubes are therefore well suited to the culture of micro-organisms. The contents only come into contact with the glass and PTFE<sup>1</sup> coating of the seal. For details of the PBT<sup>1</sup> screw cap, see pages 179–181.*

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*Typical applications: growing and storage of sterile cultures.*

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<sup>1</sup> For chemical and thermal properties, see page 222.

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
*DISPOSABLE CULTURE TUBE*  
*FROM SODA-LIME GLASS*  
*with DIN thread and PP<sup>1</sup> screw cap*

Best.-Nr. Cat. No.	d x h mm	ca. Volumen approx. vol. ml	DIN-Gewinde DIN-thread GL	Wanddicke Wall thickness mm	Verpackungseinheit/ Stück Pack/Quantity
Mit TPE Dichtung With TPE seal					
23 175 11 5	12 x 100	6	14	1	100
23 175 14 5	16 x 100	12	18	1	100
23 175 21 5	16 x 160	22	18	1	100
23 175 23 5	18 x 180	32	18	1	100

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
*SCREW CAP FOR CULTURE TUBES*  
*FROM SODA-LIME GLASS*  
*from PP<sup>1</sup> with seal*

Best.-Nr. Cat. No.	DIN-Gewinde DIN-thread GL	Verpackungseinheit/Stück Pack/Quantity
29 990 12	14	100
29 990 13	18	100

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
Rand gerade

**DISPOSABLE CULTURE TUBE**  
**FROM SODA-LIME GLASS**

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<sup>1</sup> Other dimensions available upon request.

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*Baffled flasks disrupt the laminar flow and cause a turbulent flow. The baffles increase the surface area of the liquid, the gas-exchange area and the oxygen intake. Complete with membrane screw cap<sup>1</sup>. For further information please see page 235.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN-thread	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 283 36 0	250	45	85	145	4
21 283 36 5 <sup>1</sup>	250	45	85	145	4
21 283 44 0	500	45	105	149	4
21 283 44 5 <sup>1</sup>	500	45	105	149	4

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
*Large, flat, bottom surface allows uniform culture thickness.*

*Typical application: preparation of cultures in nutrient media.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 501 43	450	117	29	100	10

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

**DURAN® CULTURE FLASK,  
FERNBACH TYPE  
bulbous shape**

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*DURAN® CULTURE FLASK,*  
*FERNBACH TYPE*  
*conical shape*

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
*DURAN® CULTURE FLASK,*  
*FERNBACH TYPE*  
*conical shape, straight neck for metal caps*

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*Large, flat, bottom surface allow uniform culture thickness. Metal cap see page 110.*

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
*DURAN® CULTURE MEDIA BOTTLE*  
*straight rim, for glass caps*

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

*DURAN® SQUARE BOTTLE*  
*after Breed-Demeter*

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
**DURAN® CULTURE MEDIA BOTTLE with beaded rim**

Best.-Nr. Cat. No.	Inhalt Capacity ml	d	d <sub>1</sub>	h	Verpackungseinheit/ Stück Pack/Quantity
21 451 24	100	50	29	115	10
21 451 39	300	70	42	168	10
21 451 44	500	83	42	207	10
21 451 54	1 000	105	46	237	10
21 451 66	2 500	150	50	315	1
21 451 73	5 000	185	54	390	1

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*Conical geometry makes the flasks particularly suitable for shaking experiments (e.g. media optimisation). Erlenmeyer flasks with screw thread are also available. These flasks can be closed with a PBT cap<sup>1</sup> or membrane cap<sup>2</sup> (permits gas exchange).*

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
**DURAN® CULTURE FLASK, ERLENMEYER SHAPE straight neck for metal caps**

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
*DURAN® CULTURE MEDIA BOTTLE*  
*straight neck for metal caps*

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
*Suitable for: culture flask no. 21 771 XX, 21 774 62 and culture media bottle Cat. No. 21 773 XX.*

Best.-Nr. Cat. No.	Material Material	Für Hals For neck Ø mm	Verpackungseinheit/ Stück Pack/Quantity
29 012 24	Edelstahl Stainless steel	38	10
29 013 24	Aluminium, blau eloxiert Aluminium, anodised blue	38	10

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*Innovative laser marking, e.g.: Numbering, bar codes, logos or names*

*No effect on the product performance. More information, see page 229*

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*Typical applications: sampling and  
 cultivation.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>i</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 461 24	100	45	17	135	10
21 461 36	250	57	17	182	10
21 461 44	500	74	17	218	10
21 461 54	1 000	95	17	265	10

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*In addition, the following individual parts  
 are available: porcelain clamp closure  
 (Cat. No. 29 701 08), replacement rubber  
 seal (Cat. No. 29 990 31) or replace-  
 ment silicone seal, autoclavable (Cat. No.  
 29 990 10).*

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
*Typical applications: sampling and cultivation.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>i</sub> mm	Verpackungseinheit/ Stück Pack/Quantity
21 465 24	100	45	17	10
21 465 36	250	57	17	10
21 465 44	500	74	17	10
21 465 54	1 000	95	17	10

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**DURAN® ROLLED FLANGE BOTTLE**  
**without closure**

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
**DURAN® ROLLED FLANGE BOTTLE**  
**with clamp closure**

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**DURAN® ROLLER BOTTLE  
 FOR CELL CULTURES**  
 with DIN thread, GL 45

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*With screw cap and pouring ring (blue, PP<sup>1</sup>). Spare caps and pouring rings see pages 33–35.*

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**STAINING JAR, COPLIN TYPE  
FROM SODA-LIME GLASS<sup>1</sup>**


Best.-Nr. Cat. No.	AD OD mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
23 319 00	66	108	10

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**Note:** Do not clean staining dishes and staining jars at temperatures above 60 °C (glass corrosion is possible).

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**STAINING DISH,  
HELLENDAHL TYPE  
FROM SODA-LIME GLASS<sup>1</sup>**


Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
23 314 00	90 x 40 x 90	10

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For 16 microscope slides 76 x 26 mm,  
with widening towards the top.

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
*STAINING DISH, SCHIEFFERDECKER TYPE*  
**FROM SODA-LIME GLASS<sup>1</sup>**

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

*GLASS BOX*  
**FROM SODA-LIME GLASS<sup>1</sup>**  
 for staining tray 21 317 00

Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
23 318 00	108 x 90 x 70	10

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
For 10 microscope slides 76 x 26 mm or each width up to 52 mm.


Best.-Nr. Cat. No.	l x b x h mm	Verpackungseinheit/Stück Pack/Quantity
21 317 00	88 x 40 x 70	10

Best.-Nr. Cat. No.	Verpackungseinheit/Stück Pack/Quantity
29 075 00	10

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ACCESSORIES:  
STAINLESS STEEL HANDLE  
for staining tray

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“WHERE WORK'S CONCERNED, I TAKE THINGS  
PRETTY SERIOUSLY ... **EXACT RESULTS**,  
EVERYTHING UNDER CONTROL ... BUT THERE'S  
AN EVEN GREATER PERFECTIONIST IN MY TEAM –  
MY MEASURING INSTRUMENT.”

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## 5 | VOLUMETRIC GLASSWARE

*DURAN® volumetric products have closely calibrated scales that permit very accurate determination and measurement of volumes. DURAN® products are available in two accuracy classes: class A/AS and class B (see the corresponding product descriptions on page 243). The two classes differ in the accuracy of measurement with class A being the highest accuracy class and class B being approximately half that of class A. Class AS has the same tolerances as class A, but is designed to permit more rapid outflow.*

*Volumetric flasks and cylinders are calibrated to measure the amount of fluid contained ("In"). Thus, for example the desired concentration can be precisely set.*

*Pipettes and burettes are calibrated to measure the amount of fluid delivered ("Ex"). This calibration takes into account surface adhesion to the glass / capillary effects; specified waiting times must be observed.*

*Volumetric flasks, volumetric and mixing cylinders as well as burettes are manufactured from DURAN® borosilicate glass 3.3 and have excellent chemical and thermal resistance. Measurement and bulb pipettes are made from soda-lime glass (see page 234 for more information on soda-lime glass).*

### Usage tips:

- *To ensure a long service life for your volumetric glassware and to exclude possible volume changes, these products should not be heated above +180 °C in drying cabinets or sterilisers.*
- *Never heat volumetric glassware on a hot plate.*
- *Always heat up and cool down volumetric glassware gradually, to avoid thermal stresses and thus any possible breakage of the glass.*

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*Calibration is based on the poured in volume (“In”) at a +20 °C reference temperature. The volume content tolerances conform to accuracy class A, the accuracy limits of the German weights and measures regulations and DIN and ISO specifications.*

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*Typical applications: precise measurement of specified liquid amounts, preparation and storage of standard solutions.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	h mm	d mm	Stopfengröße Stopper size	Verpackungseinheit/ Stück Pack./Quantity
21 678 07	5	0,025	70	24	7/16	2
21 678 08	10	0,025	90	28	7/16	2
21 678 12 <sup>3</sup>	20	0,04	110	37	10/19	2
21 678 14 <sup>3</sup>	25	0,04	110	39	10/19	2
21 678 17 <sup>3</sup>	50	0,06	140	49	12/21	2
21 678 24 <sup>3</sup>	100	0,1	170	61	12/21	2
21 678 25 <sup>3</sup>	100	0,1	170	61	14/23	2
21 678 32 <sup>3</sup>	200	0,15	210	76	14/23	2
21 678 36 <sup>3</sup>	250	0,15	220	81	14/23	2
21 678 44 <sup>3</sup>	500	0,25	260	101	19/26	2
21 678 54 <sup>3</sup>	1 000	0,4	300	127	24/29	2
21 678 63 <sup>3</sup>	2 000	0,6	370	161	29/32	2
21 678 73 <sup>2,3</sup>	5 000	1,2	470	217	34/35	1

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<sup>1</sup> For chemical and thermal resistance, see page 222.

<sup>2</sup> Non-ISO size.

<sup>3</sup> For suitable glass stoppers, see pages 50–51.

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*The large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs which increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on the poured in volume (“In”) at a +20 °C reference temperature. Mixing cylinder accuracy limits conform to DIN and ISO standards.*

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*Typical applications: diluting solutions, mixing several components with specified proportions.*


Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	h mm	d mm	NS	Verpackungseinheit/ Stück Pack./Quantity
21 618 08	10	0,2	0,2	156	16	10/19	2
21 618 14	25	0,5	0,5	190	22	14/23	2
21 618 17	50	1	1	222	25	19/26	2
21 618 24	100	1	1	285	29	24/29	2
21 618 36	250	2	2	363	39	29/32	2
21 618 44	500	5	5	395	54	34/35	2
21 618 54	1 000	10	10	500	66	45/40	1
21 618 63	2 000	20	20	540	85	45/40	1

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
*DURAN® VOLUMETRIC FLASK  
with circular graduation mark and plastic stopper from PE<sup>1</sup>, finely ground stopperbed, suitable for standard ground stoppers*

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
*DURAN® MIXING CYLINDER  
WITH HEXAGONAL BASE  
with graduation, standard ground joint  
and plastic PE<sup>1</sup> stopper*

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**DURAN® MEASURING CYLINDER WITH HEXAGONAL BASE with graduation**

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*The large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs which increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on the poured in volume („In“) at a + 20 °C reference temperature. Measuring cylinder accuracy limits conform to DIN and ISO standards (class B).*

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
*Typical applications: holding and simultaneous measurement of varying liquid amounts.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
21 396 07	5	0,1	0,1	112	13	2
21 396 08	10	0,2	0,2	137	14	2
21 396 14	25	0,5	0,5	167	21	2
21 396 17	50	1	1	196	25	2
21 396 24	100	1	1	256	29	2
21 396 36	250	2	2	331	39	2
21 396 44	500	5	5	360	53	2
21 396 54	1 000	10	10	460	65	1
21 396 63	2 000	20	20	500	85	1

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**DURAN® MEASURING CYLINDER, GRADUATED LOW FORM with hexagonal base**

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
21 395 08	10	0,2	1	90	21	10
21 395 14	25	0,5	1	115	25	10
21 395 17	50	1	2	145	29	10
21 395 24	100	1	2	165	39	10
21 395 36	250	2	5	195	54	10
21 395 44	500	5	10	250	65	10
21 395 54	1 000	10	20	285	85	1
21 395 63	2 000	20	50	340	105	1

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*With Schellbach stripe and main graduations as circular divisions. Calibration is based on the poured out volume („Ex“) at a +20 °C reference temperature. Volume content tolerances conform to DIN.*

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
*DURAN® BURETTE*  
*with straight standard ground stopcock and with glass key, class AS, 30 seconds waiting time*

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
*Work in the laboratory is simplified by the fact that unlike glass keys, the PTFE keys do not have to be lubricated.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	h mm	Ablaufzeit Run-out time s	Verpackungseinheit/ Stück Pack/Quantity
24 330 27 02	10	0,02	0,02	750	35-45	2
24 330 33 04	25	0,03	0,05	750	35-45	2
24 330 36 04	50	0,05	0,1	750	35-45	2

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
*DURAN® BURETTE*  
*with straight standard ground stopcock and with PTFE key, class AS, 30 seconds waiting time*

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**DURAN® BURETTE**  
with straight standard ground stopcock, class B

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*Calibration is based on the poured out volume („Ex“) at a + 20 °C reference temperature. Volume content tolerances conform to DIN and ISO. The Class B accuracy limit is roughly one and a half times wider than for Class AS. The tolerances are thus more strict than specified by DIN.*

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
*DURAN® MICRO-BURETTE*  
with straight standard ground stopcock,  
class AS, 30 seconds waiting time

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*With Schellbach stripe and main graduations as circular divisions, reservoir bottle (2,000 ml) and blowball.*

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*DURAN® AUTOMATIC BURETTE*  
Pellet-type, side-positioned standard ground  
stopcock, class AS, 30 seconds waiting time

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Replacement bottle for automatic burettes.


Best.-Nr. Cat. No. klar   clear	Best.-Nr. Cat. No. braun   amber	Inhalt Capacity ml	h mm	d mm	Verpackungseinheit/ Stück Pack/Quantity
Hals ungeschliffen Neck unground					
21 150 63 <sup>1</sup>	21 150 63 6	2 000	200	160	1
Mit Normalschliff NS 29/32 With standard ground NS 29/32					
21 159 63	21 159 63 6	2 000	200	160	1

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**MEASURING PIPETTE**  
**FOR PARTIAL OUTFLOW**  
**FROM SODA-LIME GLASS**  
**class B, graduated**

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Numbering from the top down. Calibration is based on the poured out volume („Ex“) at a +20 °C reference temperature. Due to the scale, variable volumes can be held and then dispensed in the same or differing increments.

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Typical applications: accurate measurement and decanting of liquids.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Teilung Graduation divisions ml	Farbkenn- zeichnung Colour code DIN 12 621	Gesamtlänge Overall length mm	Verpackungs- einheit/Stück Pack/Quantity
24 343 01 <sup>1</sup>	0,1	0,01	0,001	3 x grün   green	360	12
24 343 03 <sup>1</sup>	0,2	0,01	0,001	3 x blau   blue	360	12
24 343 06	0,5	0,008	0,01	3 x gelb   yellow	360	12
24 343 11	1	0,008	0,01	2 x gelb   yellow	360	12
24 343 16	2	0,015	0,02	2 x schwarz   black	360	12
24 343 23	5	0,04	0,05	2 x rot   red	360	12
24 343 29	10	0,08	0,1	2 x orange   orange	360	12
24 343 34	25	0,15	0,1	2 x weiß   white	450	12

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<sup>1</sup> Non-ISO size, calibrated to contain („Ex“).

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*Numbering from the top down. Calibration is based on the poured out volume (“Ex”) at a + 20 °C reference temperature. Due to the scale, variable volumes can be held and then dispensed in the same or differing increments.*

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*Numbering from the top down. Due to the scale, variable volumes can be held and then dispensed in the same or differing increments.*

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

**MEASURING PIPETTE**  
**FOR COMPLETE OUTFLOW**  
**FROM SODA-LIME GLASS**  
**class B, graduated**

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
**MEASURING PIPETTE FOR**  
**COMPLETE AND PARTIAL OUTFLOW**  
**FROM SODA-LIME GLASS**  
**class AS, main graduations as circular divisions,**  
**5 seconds waiting time**

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**BULB PIPETTE**  
**FROM SODA-LIME GLASS**  
**class B**

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
*Calibration is based on the poured out volume (“Ex”) at a +20 °C reference temperature. Calibrated to measure and discharge a single volume.*

*Typical applications: accurate measurement and decanting of liquids.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Farbkennzeichnung Colour code DIN 12 621	Ablaufzeit Run-out time s	Gesamtlänge Overall length mm	Verpackungseinheit/Stück Pack/Quantity
24 337 01	1	0,01	blau   blue	5-20	270	12
24 337 02	2	0,015	orange   orange	5-25	330	12
24 337 07	5	0,02	weiß   white	7-30	400	12
24 337 08	10	0,03	rot   red	8-40	440	12
24 337 12	20	0,05	gelb   yellow	9-50	510	6
24 337 14	25	0,05	blau   blue	10-50	520	6
24 337 17	50	0,08	rot   red	13-60	540	6
24 337 24	100	0,12	gelb   yellow	25-60	585	6

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**BULB PIPETTE**  
**FROM SODA-LIME GLASS**  
**class AS, 5 seconds waiting time**

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*Calibrated to measure and discharge a single volume (“Ex”) at a +20 °C reference temperature. Calibrated to measure and discharge a single volume.*

*Typical applications: accurate measurement and decanting of liquids.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Fehlergrenze Accuracy limits ± ml	Farbkennzeichnung Colour code DIN 12 621	Ablaufzeit Run-out time s	Gesamtlänge Overall length mm	Verpackungseinheit/Stück Pack/Quantity
24 338 01 <sup>1</sup>	1	0,007	blau   blue	5-9	270	12
24 338 02	2	0,01	orange   orange	5-9	330	12
24 338 07	5	0,015	weiß   white	7-11	400	12
24 338 08	10	0,02	rot   red	8-12	440	12
24 338 12	20	0,03	gelb   yellow	9-13	510	6
24 338 14	25	0,03	blau   blue	10-15	520	6
24 338 17	50	0,05	rot   red	13-18	540	6
24 338 24	100	0,08	gelb   yellow	25-30	585	6

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“MY LONG-TERM PROJECT ... TO SEPARATE  
THE IMPORTANT FROM THE UNIMPORTANT ...  
IF ONLY EVERYTHING WAS AS **EASY** AS  
SEPARATING USING DURAN®.”

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## 6 | GLASS FILTRATION APPARATUS AND ACCESSORIES

*Due to their high chemical and thermal shock resistance DURAN® filters and the corresponding filter plates are ideal for separations, e. g. with strong acids or alkalis. Thus they offer advantages in comparison with other materials such as plastic or paper. DURAN® filter products have a maximum operating temperature of +450 °C and are therefore far superior to other materials.*

*The corresponding filtration vessels are specially optimised to the matching filtration apparatus and are vacuum-tight due to their special geometry and high wall thickness. This characteristic has been approved by the German TÜV and marked with the "GS"-indication.*

*The glass filters are classified as being in porosity classes 0 to 5 according to their nominal maximum pore size. The following table shows the corresponding porosity range. The specified pore sizes always relate to the largest pore in the plate. This specification also characterises the minimum nominal size of particles which may be retained by the filtration.*

### Usage tips:

- The maximum permissible operating temperature is +450 °C.

- Uniform heating is recommended to avoid thermal stresses and resultant breakages.
- Heat glass filtration apparatus with disk diameters of more than 20 mm in initially cold ovens or sterilisers only.
- The heating or cooling rate should not exceed 8 °C/min.
- When filtering hot substances observe the thermal shock resistance and, if necessary, preheat the filtration apparatus in a drying cabinet.
- Wet filtration apparatus should be heated slowly up to +80 °C and dried for one hour before increasing the temperature further.

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Particular suited to use with hot or aggressive media thanks to the thermal shock and chemical resistance of DURAN®. Also available in **soda-lime glass<sup>2</sup>**.

Typical applications: decanting and filtration of substances.

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From DURAN® glass with its good thermal shock and chemical resistance. Also available in **soda-lime glass<sup>1</sup>**.

Typical applications: decanting of powdered substances and granulated material.

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


**DURAN® POWDER FUNNEL**  
with short, wide stem

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From DURAN® glass with its good thermal shock and chemical resistance. Also available in **soda-lime glass<sup>3</sup>**.

Typical applications: filtering and decanting of liquids of different densities.

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
<sup>1</sup> Non-DIN ISO size.

<sup>2</sup> Details apply only to DURAN® version.

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**DURAN® ANALYTICAL FUNNEL**  
 for quick filtration

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From DURAN® glass with its good thermal shock and chemical resistance.

Typical application: for rapid liquid filtration.

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*From DURAN® glass with its good thermal shock and chemical resistance. The ribbed form is ideal for filtering with round-paper.*

*Typical application: filtering of liquids.*

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
*Made from DURAN® glass with its good thermal shock and chemical resistance. Filter funnels mate to the filtering flask via a conical rubber seal (GUKO).*

*Typical application: filtering of liquids.*

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

From DURAN® glass with its good thermal shock and chemical resistance. Filter funnels mate to the filtering flask via a conical rubber seal (GUKO).

Typical applications: qualitative inorganic analysis and preparative chemistry.

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
From DURAN® glass with its good thermal shock and chemical resistance. The Buechner funnel features a glass support for membrane and paper filters.

Best.-Nr. Cat. No.	Inhalt Capacity ml	Platten Disc d mm	d <sub>i</sub> mm	h mm	passende Papier- Rundfilter matching filter paper d mm	AD OD mm	Ver- packungs- einheit/Stück Pack/Quantity
21 341 22	70	48	10	132	45	57	1
21 341 28	125	60	10	140	55	72	1
21 341 34	220	73	18	190	70	90	1
21 341 44	500	95	22	240	90	106	1
21 341 54	1 000	120	22	270	110	136	1

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<sup>1</sup> Substitute required porosity (see page 239) for X in the Cat No.


Best.-Nr. Cat. No.	d mm	passender Filtriervorstoß Best.-Nr. suitable filter adapter Cat. No.	Verpackungseinheit/Stück Pack/Quantity
29 201 14	26	24 316 16	10
29 201 21	33	24 316 22	10
29 201 26	41	24 316 26	10
29 201 31	49	24 316 32	10

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
**RUBBER SLEEVE from EPDM<sup>1</sup> for filter crucibles**

Best.-Nr. Cat. No.	d mm	d <sub>i</sub> mm	h mm	passende Gummimanschette Best.-Nr. suitable rubber sleeve Cat. No.	Verpackungseinheit/ Stück Pack/Quantity
24 316 16	27	10	108	29 201 14	10
24 316 22	34	10	110	29 201 21	10
24 316 26	41	10	125	29 201 26	10
24 316 32	50	10	132	29 201 31	10

**DURAN® FILTRIERVORSTOSS**


**DURAN® FILTER CRUCIBLE / FUNNEL ADAPTER**

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**DURAN® FILTER FUNNEL HEAD  
 WITH PP<sup>1</sup> FUNNEL AND TWO  
 FKM<sup>1</sup> SEALS**

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*Interchangeable filter disks. Available in three filter diameters, each of which is available in four different porosities. Important: insert filter disc between two FKM seals. After filtration, disk can be removed to allow simple and safe removal of the filtrand. Long filter disk service life, as disks are not damaged when the filtrand is scraped off. Easy cleaning of both sides is possible. Cost-effective as components and disks can be ordered individually as required.*

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
From DURAN® glass with its good thermal shock and chemical resistance. Fused glass rim.

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
**DURAN® FILTER DISK WITH GLASS RIM**  
**patented**

Best.-Nr. Cat. No.	Porosität Porosity	für Platte for disc d mm	Verpackungseinheit/ Stück Pack/Quantity
25 202 4X <sup>1</sup>	1-4	24	1
25 205 0X <sup>1</sup>	1-4	50	1
25 209 0X <sup>1</sup>	1-4	90	1

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

**FUNNEL FOR FILTER FUNNEL HEAD from PP<sup>1</sup>**

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
*Typical application: in-line filtration of gas lines to remove solid impurities (e.g. dust).*

Best.-Nr. Cat. No.	Porosität Porosity	Platte Disc d mm	d <sub>1</sub> mm	Verpackungseinheit/Stück Pack/Quantity
25 855 0X <sup>1</sup>	I-4	30	10	I
25 855 1X <sup>1</sup>	I-4	60	15	I
25 855 2X <sup>1</sup>	I-4	90	16	I

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
**DURAN® IMMERSION FILTER**  
for reverse filtration

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*Typical application: extraction of clear filtrate (not the filtrand).*

Best.-Nr. Cat. No.	Porosität Porosity	Platte Disc d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 855 6X <sup>1</sup>	I-4	35	10	210	I

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
*Typical application: reactions between gases and liquids.*

Best.-Nr. Cat. No.	Porosität Porosity	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 856 0X <sup>1</sup>	0-2	6	22	250	5
25 856 1X <sup>1</sup>	0-2	9	25	250	5
25 856 2X <sup>1</sup>	I	11	34	250	5

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*With screw-connection system. The insertion height of the head is adjustable. Individual parts can also be ordered separately (see page 140).*

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
*Typical applications: cleaning (“washing”) of gases with solvents.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	DIN-Gewinde DIN Thread GL	Olive Hose connection d mm	Verpackungseinheit/ Stück Pack/Quantity
24 713 00	500	45	9	1

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**DURAN® GAS WASHING BOTTLE**  
**head Drechsel-type, without filter disk**

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*With screw-connection system. The insertion height of the head is adjustable. The distribution of the gas in the liquid is carried out through the filter disk and thus the absorption efficiency significantly increased. Gas washing bottles also work reliably at high flow rates. Individual parts can also be ordered separately (see page 140).*

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
**DURAN® GAS WASHING BOTTLE**  
**head Drechsel-type, with filter disk**

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
**DURAN® GAS WASHING BOTTLE**  
**head with filter disk, with standard ground joint**

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
*Typical application: cleaning (“washing”) of gases with solvents.*

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*DURAN® GAS WASHING BOTTLE*  
 with fused-in filter disk,  
 with standard ground joint and cap

Best.-Nr. Cat. No.	Porosität Porosity	Inhalt Capacity	NS Joint	d	h	h <sub>1</sub>	Olive Hose connection d mm	Verpackungseinheit/Stück Pack/Quantity
25 703 01	I	350 ml	29/32	60 mm	250 mm	180 mm	10	I

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
*DURAN® MICRO IMMERSION FILTER*  
 for reverse filtration

Best.-Nr. Cat. No.	Porosität Porosity	d <sub>1</sub> mm	d mm	h mm	Verpackungseinheit/Stück Pack/Quantity
25 857 1X <sup>1</sup>	I-4	10	6	100	10

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
<sup>1</sup> Substitute required porosity (see page 239) for X in the Cat. No.


**DURAN® MIKRO-FILTERKERZE**  
 ohne Rohr

*DURAN® MICRO FILTER CANDLE*  
 without tube

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**DURAN® MICRO FILTER CANDLE**  
**with narrow tube**

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

**DURAN® MICRO FILTER FUNNEL**  
**Pregl type**

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

*Heavy walled for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". For suitable rubber gaskets, see page 147.*

*Typical applications: separations by vacuum filtration.*

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**DURAN® FILTERING FLASK WITH SIDE-ARM SOCKET**  
Erlenmeyer shape

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
21 183 36	250	85	34	155	10
21 183 44	500	105	34	185	10
21 183 54	1 000	135	45	230	10
21 183 63	2 000	166	60	255	1

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
**DURAN® FILTERING FLASK WITH SIDE-ARM SOCKET**  
*bottle shape*

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
**Note:** These filtering flasks have a ground socket 17.5/26 for vacuum tube of 15 to 18 mm OD (e.g. 6 x 5 mm or 8 x 5 mm, DIN 12 865).

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*Provision of filtration flasks with a socket has not only made work in preparation and analytical laboratories easier and simpler, but has also reduced the risk of accidents.*

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**DURAN® FILTERING FLASK WITH GLASS HOSE CONNECTION**  
*Erlenmeyer shape*

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*Due to the heavy wall thickness the apparatus is vacuum-tight. For suitable rubber gaskets, see page 147.*

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*Heavy walled for vacuum use. For suitable rubber gaskets, see page 147.*

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
<sup>1</sup> Does not conform to the "equipment and product safety regulations".

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Flaschenform

DURAN® FILTERING FLASK WITH GLASS HOSE CONNECTION  
bottle shape

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*Heavy walled for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". The plastic hose connection is replaceable. For suitable rubber seals, see page 147.*

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
*Typical application: Separations by vacuum filtration.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	d mm	d <sub>1</sub> mm	h mm	Olive Hose connection d mm	Verpackungs- einheit/Stück Pack/Quantity
21 204 24 5	100	64	24	105	9	10
21 204 36 5	250	85	34	155	9	10
21 204 44 5	500	105	34	185	9	10
21 204 54 5	1 000	135	45	230	9	10
21 204 63 5	2 000	166	60	255	9	1

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
**Flaschenform**

*DURAN® FILTERING FLASK  
 WITH KECK™ ASSEMBLY SET  
 bottle shape*

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*Due to the heavy wall thickness the apparatus is for vacuum use. These filtering flasks fulfil the regulations of the "equipment and product safety regulations". The plastic hose connections can be replaced. For suitable rubber gaskets, see page 147.*


*Typical application: separations by vacuum filtration.*

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*With removable plastic hose connection (PBT<sup>1</sup>), short and long screw (PP<sup>1</sup>), seals (VMQ<sup>1</sup>, EPDM<sup>1</sup>).  
 Suitable for filtering flasks 100 - 20 000 ml.*

Best.-Nr. Cat. No.	Olive Hose connection d mm	Verpackungseinheit/Stück Pack/Quantity
29 258 54	9	10

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
*EPDM<sup>1</sup> GUKO SET, FROM EPDM<sup>1</sup>*  
*conical rubber gaskets for filtering flasks*

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
*GUKO FROM EPDM<sup>1</sup>*  
*conical rubber gaskets for filtering flasks*

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
*Virtually universal applications, as the medium only comes into contact with glass and PTFE. The scaled funnel simplifies dosing and analysis. With PTFE plate holder. Filter paper, membrane filters (47 mm) or glass filters can be used for filtration. Plates and PTFE adapters are replaceable. Easy and fast cleaning. All components are available as spare parts.*

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*Application examples: Coarse and fine filtration, filtration of HPLC media, residue analysis*

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
<sup>1</sup> Comprises the following order numbers:  
 21 204 445, 29 202 27, 24 316 26, 29 201 26,  
 25 851 23 and 25 851 24.

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
**DURAN® LID FOR FILTER**  
**APPARATUS, WITT TYPE**  
**with tubulature, standard ground joint 29/32**

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*DURAN® FILTER APPARATUS,  
 WITT TYPE*

*complete with interchangeable lid and KECK™  
 assembly set<sup>1</sup>, standard ground joint 29/32*

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
*DURAN® FILTER APPARATUS,  
 WITT TYPE*

*complete with interchangeable lid and KECK™  
 assembly set<sup>1</sup>, standard ground joint 45/40*

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


*DURAN® LID FOR FILTER APPARATUS,  
WITT TYPE*  
*wide rough-ground tubulature, standard  
ground joint 45/40*

Best.-Nr. Cat. No.	Wasserverbrauch mind. Water flow rate min. l/h	Wasserdruck mind. Water pressure min. bar	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 362 99	300	1,2	275	1

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*DURAN® WATER JET PUMP  
with non-return valve*

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*Throughput: 400 l/h at 3.5 bar water  
pressure und 12 °C water temperature.*

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
*WATER JET PUMP  
from plastic (PP<sup>1</sup>), with non-return valve,  
hose connection and adapters for 1/2" and 3/4"*

Best.-Nr. Cat. No.	Olive passend für Schlauch ID Connection suitable for hose ID mm	Wasserverbrauch mind. Water flow rate min. l/h	Wasserverbrauch mind. Water flow rate min. bar	h mm	Verpackungs- einheit/Stück Pack/Quantity
29 250 01	9-12	170	1	235	1

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“FLEXIBLE, ACCURATE, HIGHLY VERSATILE ...

IT LOOKS AS IF I SHOULD FOLLOW THE EXAMPLE  
OF MY LABORATORY GLASSWARE.”

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## 7 | INTERCHANGEABLE GLASSWARE

*DURAN® interchangeable glassware is indispensable for laboratory work. DURAN GROUP offers a wide range of bottles and flasks with standard ground necks, vessels with flat flanges, condensers and stirrer shafts.*

*The DURAN® flat flange reaction vessels are valued for their universal suitability for use in the laboratories of a wide range of specialisations.*

*Whether reaction, distillation, evaporation or desiccation, DURAN® offers a wide range of unfinished and finished parts which always provide the optimum solution for the particular application. Due to the pure glass-glass connections, reactions with highly corrosive or highly chemically reactive substances can be carried out without problems.*

*The vessels are notable due to a robust glass flange design with an optimum external flange angle of 45°. Due to the precisely ground joint, the vessels can be closed tightly when using a sealing ring.*

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*Matching stainless-steel quick release clamps, with three retaining clips, ensure easy and safe handling.*

*All individual parts and a wide range of accessories such as lids, seals, quick-release clamps etc. are compatible and can be interchanged as required. Vessels and lids can be matched using their DN (nominal diameter) number.*

### Usage tips:

- All components are suitable for use under full vacuum and approved for operating over-pressures (see product related pages).
- Before use, it is recommended that the glass surfaces of the vessels be checked for damage such as scratches, cracks or nicks.
- Damaged glass vessels should not be used for safety reasons.
- Due to the high wall thickness and reduced thermal shock resistance under pressure loading, the flat flange vessels should be heated uniformly and gradually.

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*Thanks to the uniform wall thickness, round bottom flasks are ideal as heating vessels. The geometry permits very uniform heating.*

*Typical applications: distillation, extraction.*

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
*The conical geometry makes them ideal for small-scale reactions.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 195 08	10	14/23	30	75	10
24 195 14	25	14/23	40	90	10
24 195 20	50	14/23	50	110	10
24 195 25	100	14/23	62	125	10

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
### DURAN® ROUND BOTTOM FLASK with standard ground joint

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### DURAN® PEAR SHAPE FLASK with standard ground joint

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
**DURAN® ERLNMEYER FLASK**  
*with standard ground joint*

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
*With easy-to-read scale and large labelling field for easy marking. The conical shape makes these flasks ideal for mixing liquids and, due to the even wall thickness, also suitable for use as heating glassware.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 193 13	25	14/23	42	75	10
24 193 20	50	14/23	51	85	10
24 193 17	50	19/26	51	85	10
24 193 18 <sup>1</sup>	50	24/29	51	85	10
24 193 19	50	29/32	51	85	10
24 193 24	100	19/26	64	105	10
24 193 26 <sup>1</sup>	100	24/29	64	105	10
24 193 27	100	29/32	64	105	10
24 193 32 <sup>1</sup>	200	29/32	79	131	10
24 193 36	250	24/29	85	140	10
24 193 37	250	29/32	85	140	10
24 193 38 <sup>1</sup>	250	45/40	85	140	10
24 193 39 <sup>1</sup>	300	29/32	87	155	10
24 193 44	500	24/29	105	175	10
24 193 46	500	29/32	105	175	10
24 193 47 <sup>1</sup>	500	45/40	105	175	10
24 193 54	1 000	24/29	131	220	10
24 193 56	1 000	29/32	131	220	10
24 193 57 <sup>1</sup>	1 000	45/40	131	220	10

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*See pages 176–178 for the Keck™ clips for conical joints.*

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
*Due to the flat bottom the flask can be set upon a bench without a support ring.*

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**DURAN® FLAT BOTTOM FLASK  
with standard ground joint**

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hals Neck NS	d mm	h mm	Verpackungseinheit/ Stück Pack/Quantity
24 171 19	50	29/32	51	85	10
24 171 24	100	19/26	64	103	10
24 171 26	100	24/29	64	103	10
24 171 27	100	29/32	64	103	10
24 171 36 <sup>1</sup>	250	24/29	85	130	10
24 171 37	250	29/32	85	130	10
24 171 44 <sup>1</sup>	500	24/29	105	160	10
24 171 46	500	29/32	105	160	10
24 171 54 <sup>1</sup>	1 000	24/29	131	187	10
24 171 56	1 000	29/32	131	187	10
24 171 63 <sup>1</sup>	2 000	29/32	166	230	10

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

**DURAN® EVAPORATING FLASK  
with standard ground joint, pear shape**

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**DURAN® TWIN-NECK, ROUND  
BOTTOM FLASK**  
with standard ground joint, inclined side neck

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*Thanks to the uniform wall thickness, round bottom flasks are ideal as heating vessels. The geometry permits very uniform heating. Depending upon the application, accessories, columns, thermometers, dropping funnels, boiling capillaries, etc. can be fitted.*

*Typical applications: distillation, extraction.*

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
**DURAN® TRIPLE-NECK, ROUND  
BOTTOM FLASK**  
with standard ground joint, inclined side necks

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
*Typical application: distillation.*

Best.-Nr. Cat. No.	Nutzlänge Effective length mm	Hülse Socket size NS	Kern Cone NS	Gesamtlänge Overall length mm	Verpackungs- einheit/Stück Pack/Quantity
24 240 71	300	24/29	24/29	450	1
24 240 72	300	29/32	29/32	450	1
24 240 87	500	24/29	24/29	650	1
24 240 88	500	29/32	29/32	650	1

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**DURAN® TRIPLE-NECK, ROUND  
BOTTOM FLASK**


*with standard ground joint, parallel side necks*

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**DURAN® VIGREUX COLUMN**

*with 2 standard ground joints, complete with  
slide-on glass jacket*

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**DURAN® WOLFF BOTTLE**  
**3 standard ground necks**

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
*Vacuum resistant due to the wall thickness and geometry.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Hälse Necks NS	d mm	Verpackungseinheit/ Stück Pack./Quantity
24 709 44	500	19/26	87	1
24 709 54	1 000	24/29	113	1
24 709 63	2 000	29/32	135	1
24 709 73	5 000	34/35	185	1

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**DURAN® WOLFF BOTTLE**  
**3 standard ground necks and bottom tubulature**

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*Pressure and vacuum resistant due to the wall thickness and geometry. Groove for O-ring seal (see pages 180–181 for suitable sealing rings). For suitable KPG stirrer bearings and shafts, see pages 171–173.*

*Typical applications: reactions under pressure and/or high temperature.*

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


**Note:** At the maximum usage temperature of 250 °C and the maximum operating pressure, the temperature difference in the glass wall of the flat flange reaction vessels must not exceed 30 °C.

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*Pressure and vacuum resistant due to the wall thickness and geometry. Groove for O-ring seal (see pages 180–181 for suitable sealing rings).*

*Typical applications: reactions under pressure and/or high temperature.*

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<sup>1</sup> Suitable for Witt-type filter apparatus.

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**Note:** only heat flat-flange beakers in water or oil baths. At the maximum usage temperature of 250 °C and the maximum operating pressure, the temperature difference in the glass wall of the flat flange reaction vessels must not exceed 30 °C.

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*Pressure and vacuum resistant due to the wall thickness and geometry. For suitable KPG stirrer bearings and shafts, see pages 171–173.*

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

**DURAN® FLAT FLANGE LID**  
with centre neck

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Pressure and vacuum resistant due to the wall thickness and geometry.

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
Typical applications: reactions under pressure and/or high temperature.

Best.-Nr. Cat. No.	DN	Flansch AD Flange OD mm	h mm	Hals Neck NS	max. Betriebsüber- druck bei 250 °C max. operating pressure at 250 °C	Verpackungs- einheit/Stück Pack/Quantity
NS 29/32						
24 398 46	100	138	76	29/32	1 bar	1
24 398 51	120	158	105	29/32	1 bar	1
24 398 57	150	184	102	29/32	1 bar	1
24 398 61	200	242	126	29/32	1 bar	1
NS 45/40						
24 450 46	100	138	84	45/40	1 bar	1
24 450 57	150	184	112	45/40	1 bar	1

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Accessories for flat flange vessels, comprising an elastic, silicone core with seamless FEP<sup>1</sup> coating that encloses the ring. The combination of these high-quality materials achieves good elasticity in conjunction with outstanding chemical resistance.

Best.-Nr. Cat. No.	passend für Flansch suitable for flange DN	d x s mm	Verpackungseinheit/ Stück Pack/Quantity
29 222 34	60	75 x 4	1
29 222 46	100	110 x 4	1
29 222 51	120	133 x 4	1
29 222 57	150	157 x 5	1
29 222 61	200	215 x 5	1

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Accessories for flat flange vessels. From silicone (VMQ<sup>1</sup>), so highly elastic. The chemical resistance is reduced in comparison with FEP<sup>1</sup> coated O-rings.


**O-RING, TRANSPARENT**  
aus Silikon (VMQ<sup>1</sup>)

**O-RING, TRANSPARENT**  
from silicone (VMQ<sup>1</sup>)

Best.-Nr. Cat. No.	passend für Flansch suitable for flange DN	d x s mm	Verpackungseinheit/ Stück Pack/Quantity
29 225 34	60	75 x 4	5
29 225 46	100	110 x 4	5
29 225 51	120	133 x 4	5
29 225 57	150	157 x 5	5
29 225 61	200	215 x 5	5

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
Accessories for flat flange vessels.

Best.-Nr. Cat. No.	passend für Flansch suitable for flange DN	Verpackungseinheit/Stück Pack/Quantity
29 071 34	60	
29 071 46	100	
29 071 51	120	
29 071 57	150	
29 071 61	200	

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**QUICK RELEASE CLAMP**  
from stainless steel, with 3 retaining clips

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
For secure fitting of the lid or the reaction vessel using two clamping rods.

Best.-Nr. Cat. No.	passend für Flansch suitable for flange DN	Verpackungseinheit/Stück Pack/Quantity
29 073 46	100	
29 073 57	150	

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
**HOLDING DEVICE FOR REACTION VESSELS**  
from chrome-nickel steel

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**DURAN® LIEBIG CONDENSER**  
**(WEST CONDENSER)**  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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*Relatively small heat exchange surface and thus relatively low cooling capacity.*

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*Typical application: product condenser for distillate separations.*


Best.-Nr. Cat. No.	Mantellänge Jacket length mm	Hülse Socket size NS	Kern Cone NS	Gewinde Thread GL	Verpackungs- einheit/Stück Pack/Quantity
24 251 61	160	14/23	14/23	14	1
24 251 70 <sup>1</sup>	250	14/23	14/23	14	1
24 251 71 <sup>1</sup>	250	24/29	24/29	14	1
24 251 72	250	29/32	29/32	14	1
24 251 81 <sup>1</sup>	400	24/29	24/29	14	1
24 251 82	400	29/32	29/32	14	1

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**DURAN® BULB CONDENSER**  
**(ALLIHN CONDENSER)**  
with 2 standard ground joints and  
2 screw-on plastic hose connections

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*The bulb condenser has a greater cooling surface than the Liebig condenser and thus higher cooling capacity.*

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*Typical application: reflux condenser for condensation and feedback of the (solvent) vapour to the reaction mixture.*

Best.-Nr. Cat. No.	Mantellänge Jacket length mm	Hülse Socket size NS	Kern Cone NS	Gewinde Thread GL	Verpackungs- einheit/Stück Pack/Quantity
24 252 71 <sup>1</sup>	250	24/29	24/29	14	1
24 252 72	250	29/32	29/32	14	1
24 252 81 <sup>1</sup>	400	24/29	24/29	14	1
24 252 82	400	29/32	29/32	14	1

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
**DURAN® COIL DISTILLATE  
CONDENSER**

*with 2 standard ground joints and  
2 screw-on plastic hose connections*

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*The Dimroth condenser comprises a coil condenser located within a tube. This condenser type has a larger heat exchange surface and thus a better cooling effect than the Liebig or Allihn condenser.*

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*Typical applications: product and reflux condenser.*

Best.-Nr. Cat. No.	Mantellänge Jacket length mm	Hülse Socket size NS	Kern Cone NS	Gewinde Thread GL	Verpackungs- einheit/Stück Pack/Quantity
24 254 61 <sup>1</sup>	160	14/23	14/23	14	1
24 254 71 <sup>1</sup>	250	24/29	24/29	14	1
24 254 72	250	29/32	29/32	14	1
24 254 82	400	29/32	29/32	14	1

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

**DURAN® JACKETED COIL**  
**CONDENSER**

with 2 standard ground joints and  
2 screw-on plastic hose connections

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*This type of condenser has a very large heat exchange surface due to its coil condenser and double jacket and is thus especially suited to working with low-boiling point media.*

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*Typical application: use as a reflux condenser for condensation and feedback of the (solvent) vapour to the reaction mixture.*


Best.-Nr. Cat. No.	Mantellänge Jacket length mm	Hülse Socket size NS	Kern Cone NS	Gewinde Thread GL	Verpackungs- einheit/Stück Pack/Quantity
24 255 71 <sup>1</sup>	250	24/29	24/29	14	1
24 255 72	250	29/32	29/32	14	1
24 255 81 <sup>1</sup>	400	24/29	24/29	14	1
24 255 82	400	29/32	29/32	14	1

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**DURAN® DROPPING FUNNEL**  
cylindrical, with scale and ground joint

DIN ISO  
4800

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*With standard ground stopcock and retaining device.*

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*Typical applications: uniform and metered liquid supply to a reaction mixture. The rate of supply can be adjusted.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Teilung Scale ml	d mm	h mm	NS-Glasküken Nenngröße Standard solid key nominal size	Hülse Socket size NS	Ver- packungs- einheit/ Stück Pack/ Quantity
24 122 17	50	1	9	279	3 NS	19/26	1
24 122 24	100	2	9	299	3 NS	19/26	1
24 122 36	250	5	10	381	4 NS	29/32	1
24 122 44	500	10	10	431	4 NS	29/32	1
24 122 54	1 000	20	13	506	6 NS	29/32	1

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*With standard ground joint, standard ground stopcock and retaining device.*

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*Typical application: uniform and metered liquid supply to a reaction mixture. The rate of supply can be adjusted.*

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
**DURAN® DROPPING FUNNEL**  
**cylindrical, with scale and ground joint**

Best.-Nr. Cat. No.	Inhalt Capacity ml	Teilung Scale ml	Kern Cone size NS	h mm	NS-Glasküken Nenngröße Standard solid key nominal size	Hülse Socket size NS	Ver- packungs- einheit/ Stück Pack/ Quantity
24 124 20	50	1	14/23	220	3 NS	19/26	1
24 124 25	100	2	14/23	240	3 NS	19/26	1
24 124 24	100	2	19/26	240	3 NS	19/26	1
24 124 36 <sup>1</sup>	250	5	24/29	320	4 NS	29/32	1
24 124 37	250	5	29/32	320	4 NS	29/32	1
24 124 44 <sup>1</sup>	500	10	24/29	400	4 NS	29/32	1
24 124 46	500	10	29/32	400	4 NS	29/32	1
24 124 56	1 000	20	29/32	480	6 NS	29/32	1

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*With pressure equalisation tube, standard ground cone, standard ground stopcock and retaining device.*

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
**DURAN® DROPPING FUNNEL**  
**cylindrical, with scale**

Best.-Nr. Cat. No.	Inhalt Capacity ml	Teilung Scale ml	Kern Cone size NS	h mm	NS-Glasküken Nenngröße (DIN 12 541) Standard solid key nominal size (DIN 12 541)	Hülse Socket size NS	Ver- packungs- einheit/ Stück Pack/ Quantity
24 125 20	50	1	14/23	240	3 NS	19/26	1
24 125 25	100	2	14/23	270	3 NS	19/26	1
24 125 24	100	2	19/26	270	3 NS	19/26	1
24 125 36 <sup>1</sup>	250	5	24/29	350	4 NS	29/32	1
24 125 37	250	5	29/32	380	4 NS	29/32	1
24 125 44 <sup>1</sup>	500	10	24/29	430	4 NS	29/32	1
24 125 46	500	10	29/32	430	4 NS	29/32	1

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
**DURAN® SEPARATING FUNNEL**  
**conical shape**

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
*With standard ground stopcock, retaining device and plastic stopper. The conical shape makes it highly suited to phase separation (See pages 50–51 for replacement stoppers.).*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Stiel Stern  d <sub>1</sub> mm	h mm	NS-Glasküken Nenngröße (DIN 12 541) Standard solid key nominal size (DIN 12 541)	Hülse Socket size  NS	Verpackungs- einheit/Stück Pack/Quantity
Mit Glasküken With solid key						
24 294 17	50	9	190	3 NS	19/26	1
24 294 24	100	9	230	3 NS	19/26	1
24 294 36	250	10	280	4 NS	29/32	1
24 294 44	500	10	320	4 NS	29/32	1
24 294 54	1 000	13	380	6 NS	29/32	1
24 294 63	2 000	13	430	6 NS	29/32	1
Mit PTFE-Küken With PTFE key						
10 648 05	100	9	230	3 NS	19/26	1
10 648 06	250	10	280	4 NS	29/32	1
10 648 07	500	10	320	4 NS	29/32	1
10 648 09	1 000	13	380	6 NS	29/32	1

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**DURAN® SEPARATING FUNNEL**  
**spherical, heavy-duty version**

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*With standard ground stopcock and standard ground stopper (See pages 50–51 for replacement stoppers.)*

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*Typical application: Phase separation.*

Best.-Nr. Cat. No.	Inhalt Capacity ml	Stopfen Stopper NS	h mm	d mm	Verpackungs- einheit/Stück Pack/Quantity
24 291 36	250	24/20	235	90	1
24 291 44	500	24/29	276	115	1
24 291 54	1 000	29/32	295	132	1
24 291 66	2 500	45/40	370	182	1
24 291 73	5 000	45/40	425	222	1
24 291 86	10 000	45/40	490	286	1

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

*Ground and polished bearing surface, with standard ground cone.*

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
*Ground and polished bearing surface.*

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

*Ground and polished bearing surface,  
with GL screw thread tube and standard  
ground cone.*

Best.-Nr. Cat. No.	Kennzeichnung Designation	d mm	h mm	DIN-Gewinde DIN thread GL	Kern Cone size NS	Verpackungs- einheit/Stück Pack/Quantity
24 750 08	HB 10	10	75	32	24/29	1
24 750 09	HB 10	10	75	32	29/32	1

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*Bearing surface ground and polished.*

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
Bearing surface ground and polished.

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
**DURAN® DRYING TUBE, BENT**  
with standard ground cone

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

*DURAN® CONNECTION PIECE  
with standard ground cone, 90° angle*


Best.-Nr. Cat. No.	Kern Cone size NS	Verpackungseinheit/Stück Pack/Quantity
24 300 06	14/23	1
24 300 08	24/29	1
24 300 09	29/32	1

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*DURAN® RECEIVER ADAPTER  
bent, with ground socket*

Best.-Nr. Cat. No.	Hülse Socket size NS	Verpackungseinheit/Stück Pack/Quantity
24 310 06	14/23	1
24 310 08	24/29	1
24 310 09	29/32	1

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See pages 180–182 for screw caps and  
hose connections.

Best.-Nr. Cat. No.	Hülse Socket size NS	Kern Cone size NS	DIN-Gewinde DIN thread GL	Olive Hose connection d mm	Verpackungs- einheit/Stück Pack/Quantity
24 130 21	14/23	14/23	14	8,6	1
24 130 34	24/29	24/29	14	8,6	1
24 130 46	29/32	29/32	14	8,6	1

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**DURAN® VACUUM RECEIVER  
ADAPTER**

bent, with 2 standard ground joints  
and screw-on plastic hose connection

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

**DURAN® VACUUM RECEIVER  
ADAPTER**

straight, with 2 standard ground joints  
and screw-on plastic hose connection

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**DURAN® ADAPTER**  
*with standard ground joint*

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
**KECK™ CLIP**  
*for conical, joints, from POM<sup>1</sup>*

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
**KECK™ CLIP ASSORTMENT**  
for conical, joints, from POM<sup>1</sup>

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**KECK™ CLIP ASSORTMENT**  
for conical, joints from metal<sup>1</sup>

Best.-Nr. Cat. No.	Für Nenngröße For nominal size	Verpackungseinheit/Stück Pack/Quantity
29 030 02	NS 7	10
29 030 03	NS 10	10
29 030 04	NS 12	10
29 030 06	NS 14	10
29 030 07	NS 19	10
29 030 08	NS 24	10
29 030 09	NS 29	10
29 030 11	NS 34	10

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
**KECK™ CLIP ASSORTMENT**  
for conical, joints from stainless steel

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**KECK™ CLIP ASSORTMENT**  
for spherical joints, from POM<sup>1</sup>

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
For flexible connection of two glass screw-thread connections. With integral silicone (VMQ<sup>1</sup>) seal.


Best.-Nr. Cat. No.	DIN-Gewinde DIN thread GL	Verpackungseinheit/Stück Pack/Quantity
29 226 05 5	14	1
29 226 06 5	18	1
29 226 09 5	25	1
29 226 08 5	32	1

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**DURAN® SCREWTHREAD TUBE**  
with DIN thread

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
*DURAN® SCREWTHREAD TUBE  
with DIN thread and standard ground cone*

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*High leak tightness through use of PTFE<sup>1</sup> coated silicone seal<sup>2</sup>. More chemically resistant than PP<sup>1</sup> cap.*

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

*Suitable for PBT<sup>1</sup> screw caps tamper-evident screw caps<sup>2</sup>. Heat resistance: 130 °C (vapour) and 200 °C (dry heat). Good chemical resistance due to PTFE<sup>1</sup> coating. Silicon is peroxidically cured.*

Best.-Nr. Cat. No.	für Schraubverschlüsse, rot for screw-caps, red GL	Verpackungseinheit/Stück Pack/Quantity
29 248 08	14	10
29 248 11	18	10
29 248 13	25	10
29 248 19	32	10
29 248 28	45	10

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
**SCREW CAP WITH APERTURE**  
**from PBT<sup>1</sup>, red**

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*Suitable for silicone seal for piercing (Septa). More chemically resistant than PP<sup>1</sup> cap.*

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*Typical applications: Injection or removal of media.*

Best.-Nr. Cat. No.	DIN-Gewinde DIN thread GL	Bohrung Aperture bore d mm	Verpackungseinheit/ Stück Pack/Quantity
29 227 05	14	9,5	10
29 227 06	18	11	10
29 227 09	25	15	10
29 227 08	32	20	10
29 227 10	45	34	10

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Suitable for PBT<sup>1</sup> screw cap with aperture. Heat resistance: 130 °C (vapour) and 200 °C (dry heat).

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
**SILICONE SEALING RING (VMQ<sup>1</sup>)**  
with bonded PTFE<sup>1</sup> face

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Typical application: fixing of tubes.

Best.-Nr. Cat. No.	für DIN-Gewinde for DIN thread GL	Größe Size d x d <sub>1</sub> mm	für Rohr for tube Ø mm	Verpackungseinheit/ Stück Pack/Quantity
29 234 06	14	12 x 6	5,5 - 6,5	10
29 235 06	18	16 x 6	5,5 - 6,5	10
29 235 08	18	16 x 8	7,5 - 9,0	10
29 235 10	18	16 x 10	9,0 - 11,0	10
29 237 08	25	22 x 8	7,5 - 9,0	10
29 237 10	25	22 x 10	9,0 - 11,0	10
29 237 12	25	22 x 12	11,0 - 13,0	10
29 236 10	32	29 x 10	9,0 - 11,0	10
29 236 12	32	29 x 12	11,0 - 13,0	10
29 236 14	32	29 x 14	13,0 - 15,0	10
29 236 16	32	29 x 16	15,0 - 17,0	10
29 236 18	32	29 x 18	17,0 - 19,0	10
29 238 26	45	42 x 26	25,0 - 27,0	10
29 238 32	45	42 x 32	31,0 - 33,0	10

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Suitable for PBT<sup>1</sup> screw cap with aperture. Heat resistance: 130 °C (vapour) and 200 °C (dry heat). Silicon is peroxidically cured.

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
Typical applications: injection or removal of media.

Best.-Nr. Cat. No.	für DIN-Gewinde for DIN thread GL	d mm	Dicke Thickness mm	Verpackungseinheit/ Stück Pack/Quantity
29 246 05	14	12	2	100
29 246 06	18	16	2	100
29 246 09	25	22	2	100
29 246 08	32	29	2	100
29 246 10	45	42	3	100

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
**KUNSTSTOFF-OLIVE**  
**gerade, aus PP<sup>1</sup>**

*PLASTIC HOSE CONNECTION*  
*straight, from PP<sup>1</sup>*

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*With silicone seal (VMQ<sup>1</sup>). Suitable for GL 14 screw cap (Cat. No. 29 227 05).*


Best.-Nr. Cat. No.	d mm	d <sub>1</sub> mm	Verpackungseinheit/ Stück Pack/Quantity
29 255 06	8,6	5	10

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**KUNSTSTOFF-OLIVE**  
**gebogen, aus PP<sup>1</sup>**

*PLASTIC HOSE CONNECTION*  
*bent, from PP<sup>1</sup>*

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
*For plastic hose connections (Cat. No. 29 255 and 29 247 05). From silicone (VMQ<sup>1</sup>).*

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
*"I USED TO THINK: GLASS IS JUST GLASS...  
TODAY I KNOW THAT YOU CAN'T TELL QUALITY  
AT FIRST GLANCE... BUT IT'S IN THE NAME."*

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## WHAT IS GLASS ?

Glass is an inorganic mixture fused at high temperature which solidifies on cooling but does not crystallize. Its basic components, network formers and modifiers, are present in the common glasses in the form of oxides.

Typical glass formers (network formers) are silicon dioxide ( $\text{SiO}_2$ ), boric acid ( $\text{B}_2\text{O}_3$ ), phosphoric acid ( $\text{P}_2\text{O}_5$ ) and aluminium oxide ( $\text{Al}_2\text{O}_3$ ). These substances are capable of absorbing (dissolving) metal oxides up to a certain proportion without losing their glassy character. This means that the incorporated oxides are not involved in the formation of the glass but modify certain physical properties of the structure of the glass as "network modifiers".

A large number of chemical substances have the property that they solidify from the molten state into a glassy state. The formation of glass depends on its cooling rate and a necessary prerequisite is the existence of mixed types of bond (covalent bonds and ionic bonds) between the atoms or groups of atoms.

As a result, glass-forming products show a strong tendency whilst still in the molten state towards amorphous three-dimensional networking through polymerisation. Crystals are formed when the individual atoms form a regular three-dimensional arrangement in what is known as a "crystal lattice" as soon as the particular substance changes from the liquid to the solid state. Glass, however, forms a largely amorphous "network" when it cools down from the molten state. The components mainly involved in the formation of the glass are therefore described as "network formers". The glass-forming molecules in this network can incorporate ions that open up the network at certain points, changing its structure and thus the properties of the glass. They are therefore called "network modifiers".

## WHAT IS DURAN® ?

### The special features of DURAN®

Very high chemical resistance, nearly inert behaviour, a high usage temperature, minimal thermal expansion and the resultant high resistance to thermal shock are its most significant properties. This optimum physical and chemical performance makes DURAN® the ideal material for use in the laboratory and for the manufacture of chemical apparatus used in large-scale industrial plant. It is also widely used on an industrial scale in all other application areas in which extreme heat resistance, resistance to thermal shock, mechanical strength and exceptional chemical resistance are required.

### Chemical composition of DURAN®

DURAN® has the following approximate composition:

81	% by weight	$\text{SiO}_2$
13	% by weight	$\text{B}_2\text{O}_3$
4	% by weight	$\text{Na}_2\text{O}/\text{K}_2\text{O}$
2	% by weight	$\text{Al}_2\text{O}_3$

DURAN® properties are specified in DIN ISO 3585. In contrast to other borosilicate 3.3 glasses, DURAN® is notable for its highly consistent, technically reproducible quality.


### Chemical properties

The chemical resistance especially of DURAN® glass is more comprehensive than that of all other known materials. DURAN® borosilicate glass is highly resistant to water, acids, saline solutions, organic substances and also halogens such as chlorine and bromine. Its resistance to alkali is also relatively good. Only hydrofluoric acid, boiling phosphoric acid and strong alkalis cause appreciable surface removal of the glass (glass corrosion) at elevated temperatures (>100 °C). Due to the nearly inert behaviour, there are no interactions (e.g. ion exchange) between medium and glass and any spurious influence on experiments is thereby effectively excluded.

#### Hydrolytic resistance


DURAN® corresponds to Class 1 of the glasses that are divided into a total of 5 hydrolytic resistance classes in accordance with ISO 719 (98 °C). The amount of Na<sub>2</sub>O/g glass grain leached out after 1 hour in water at 98 °C is measured. For DURAN® the quantity of Na<sub>2</sub>O leached out is less than 31 µg/g of glass grain. DURAN® also corresponds to Class 1 of the glasses divided into a total of 3 hydrolytic resistance classes in accordance with ISO 720: (121 °C). The quantity of Na<sub>2</sub>O leached out after 1 hour in water at 121 °C is less than 62 µg/g of glass grain. Due to its good hydrolytic resistance, DURAN® meets the requirements of the USP, JP and EP for a neutral glass that corresponds to glass type 1. It can therefore be used in an almost unrestricted way in pharmaceutical applications and in contact with foodstuffs.

Hydrolytic attack on DURAN® as a function of time (100 °C)

#### Acid resistance


DURAN® corresponds to Class 1 of the glasses divided into 4 acid classes in accordance with DIN 12116. As the surface removal after boiling for 6 hours in normal HCl is less than 0.7 mg/100 cm<sup>2</sup>, DURAN® is classed as acid-resistant borosilicate glass. The quantity of alkaline metal oxides leached out in accordance with ISO 1776 is less than 100 µg Na<sub>2</sub>O/100 cm<sup>2</sup>.

Acid attack on DURAN® as a function of acid concentration

Alkali resistance

DURAN® corresponds to Class 2 of the glasses divided into 3 alkali classes in accordance with DIN ISO 695. The surface erosion after 3 hours boiling in a mixture of equal volume fractions of sodium hydroxide solution (concentration 1 mol/l) and sodium carbonate solution (concentration 0.5 mol/l) is only 134 mg/100 cm<sup>2</sup>.

Alkali attack on DURAN® as a function of temperature

Overview of the chemical properties of technical glasses

Description	Chemical resistance class		
	Hydrolytic resistance DIN ISO 719	Acid resistance DIN 12 116	Alkali resistance ISO 695
DURAN®	1	1	2
FIOLAX®	1	1	2
Soda-lime glass	3	1	2
SBW	1	1	1

### Physical properties

#### Temperature resistance when heated and thermal shock resistance

The maximum temperature for short-term use for DURAN® is 500 °C. Above a temperature of 525 °C the glass begins to soften and above a temperature of 860 °C it changes to the liquid state. As it has a very low coefficient of linear expansion ( $\alpha = 3.3 \times 10^{-6} \text{ K}^{-1}$ ), a feature of DURAN® is its high thermal shock resistance up to  $\Delta T = 100 \text{ K}$ . For a temperature change of 1 K, the glass changes by only  $3.3 \times 10^{-6}$  relative length units, resulting in low levels of mechanical strain were a thermal gradient exists. The thermal shock resistance is influenced wall thickness and product geometry.

#### Temperature resistance at low temperatures

DURAN® can be cooled down to the maximum possible negative temperature and is therefore suitable for use with liquid nitrogen (approx.  $-196^\circ\text{C}$ ). During use / freezing special attention should be given to the expansion of the content. In general DURAN® products are recommended for use down to  $-70^\circ\text{C}$ .

When working at low temperatures, the effect of any expansion of a DURAN® vessel's content must be borne in mind. During cooling and thawing ensure that the temperature difference does not exceed 100 K. In practice, therefore, stepwise cooling and heating are recommended. When freezing substances in such items as DURAN® bottles or DURAN® test tubes, the container should only be filled to a maximum of 3/4 of its capacity. Moreover, it should be frozen slanted at an angle of  $45^\circ$  (to enlarge the surface area). The minimum service temperature is dependant upon the properties of any screw caps or other components used. For the blue PP screw cap the minimum temperature is  $-40^\circ\text{C}$ .

#### Use in the microwave

DURAN® laboratory glassware is suitable for use in microwaves. This also applies to plastic coated DURAN® products .

#### Overview of the physical properties of technical glasses

Description	Linear expansion coefficient $\alpha$ (20 °C/300 °C) [ $10^{-6}\text{K}^{-1}$ ]	Transformation temperature [°C]	Density [g/ cm <sup>3</sup> ]
DURAN®	3.3	525	2.23
FIOLAX®	4.9	565	2.34
Soda-lime glass	9.1	525	2.5
SBW	6.5	555	2.45



### Optical properties

In the spectral range from about 310 to 2200 nm the absorption of DURAN® is negligibly low. It is clear and colourless. Fairly large layer thicknesses (axial view through pipes) appear slightly yellow/greenish. Amber-coloured DURAN® products are suited to use with light-sensitive substances (see amber colouring of DURAN®). This results in strong absorption in the short-wave region up to approx. 500 nm. In photochemical processes the light transmission of DURAN® in the ultraviolet range is of particular importance. The degree of light transmission in the UV range indicates the ease with which photochemical reactions can be carried out, for example chlorinations and sulfochlorination. The chlorine molecule absorbs light in the range from 280 to 400 nm and thus serves as a transmitter of the radiation energy.

**Amber colouring of DURAN® laboratory glassware**

Amber colouring enables storage of light sensitive substances in DURAN® products. Light transmission in the wavelength range between 300 and 500 nm is, in comparison with DURAN® clear glass, < 10%. Accordingly amber DURAN® glass corresponds to USP/EP/JJP specifications.

To colour the products, a special diffusion colour is sprayed solely on the outer surface of the clear glass articles with an innovative spraying method. This technology results in high uniform amber coloring. Afterwards, the coating is burned-in and is therefore resistant to chemicals and cleaning in a dishwasher. The proven DURAN® properties on the inner surface remain unaffected; there is no contact or interaction between contents and amber coating. The uniformity of the amber colouring process ensures the quality of the amber colour which is assured by continuous monitoring.

Transmission curves for DURAN®

**CONFORMITY WITH STANDARDS AND GUIDELINES**

Alongside the international standard DIN ISO 3585, in which the properties of borosilicate glass 3.3 are defined, DURAN® laboratory glassware corresponds to the current standards for glass laboratory apparatus. The relevant DIN/ISO standards are given on the product pages of this catalogue. If the standard is changed, e.g. in case of harmonisation to ISO, our dimensions are adjusted accordingly within an appropriate time interval.

DURAN® is a neutral glass of high hydrolytic resistance and thus belongs to glass type 1 in accordance with the European pharmacopeia (EP, chapter 3.2.1), the Japanese pharmacopeia (JP, chapter 7.01) and the United States pharmacopeia (USP section: 660) and National Formulary.

## LABORATORY GLASSWARE AND PLASTICS

### Plastics used with laboratory glass

To complement DURAN® laboratory glassware products, various plastic products such as screw caps are available. Their properties are listed in the following table.

		Temperature resistance range °C
EPDM	Ethylene/propylene-diene-rubber	-45 to +150
ETFE	Partially crystalline ethylene/tetrafluoroethylene copolymer	-100 to +180
EVA	Ethylene-vinyl acetate	-80 to +70
FEP	Tetra-Fluor-Ethylene/Hexafluor-Propylene	-200 to +200
FKM	Fluorinated rubber	-20 to +200
PBT	Polybutylenterephthalat	-45 to +180
PE	Polyethylene	-40 to +80
POM	Polyoxymethylene	-40 to +90
PP	Polypropylene	-40 to +140
PTFE	Polytetrafluoroethylene	-200 to +260
PU	Polyurethane	-30 to +135
TpCh260	Thermoplastic/duroplastic	-196 to +260
VMQ	Silicone rubber	-50 to +200
PSU Compound	Compound polyarylsulfone based	-45 to +180

### Chemical resistance of plastics

Classes of substances +20 °C	PE	PP	PBT	PTFE/FEP	TpCh260	ETFE	VMQ	EPDM	PU	FKM	POM	PSU Compound
Alcohols, aliphatic	+	+	++	++	++	++	+	+	++	-	+	++
Aldehydes	+	+	+	++	++	++	+		++		+	+
Alkaline solutions	++	++	+	++	++	++	-	++	++	-	+	++
Esters	+	+	+	++	++	++	-	++	+	-	-	+
Ethers	-	-	+	++	++	++	-	-	+	-	+	+
Hydrocarbons, aliphatic	-	++	+	++	++	++	-	++	++	++	+	+
Hydrocarbons, aromatic	-	+	+	++	++	++	-	+	++	++	+	-
Hydrocarbons, halogenated	-	+		++	++	++	-	+	-	++	+	-
Ketones	+	+	+	++	++	+	-	++	+	-	+	-
Acids, dilute or weak	+	++	++	++	++	++	-	++	++	++	-	++
Acids, conc. or strong	+	+	+	++	++	++	-	++	+	++	-	++
Acids, oxidising	-	+	-	++	++	+	-	-	+	+	-	+

++ = very good resistance  
 + = good to limited resistance  
 - = low resistance

## CLEANING LABORATORY GLASSWARE

Special glass laboratory apparatus can be washed by hand in a soaking bath or by machine in a lab washer. Laboratory dealers can supply a wide range of detergents and cleaner-disinfectants for both methods. As contamination during the delivery of our laboratory glassware cannot be totally ruled out, we recommend washing laboratory glassware before it is used for the first time. To care properly for laboratory glassware, it should be washed immediately after use at low temperature, on a short cycle and with low alkalinity. Laboratory apparatus that has come into contact with infectious substances or microorganisms should be treated in accordance with the current guidelines. Dependent on the substance, autoclaving (e.g. to kill microorganisms) may be necessary prior to cleaning, but it is generally recommended that cleaning or washing of glass products be carried out prior to autoclaving or hot-air sterilisation. This prevents dirt or impurities from adhering to the glassware surfaces and prevent damage caused by any possibly adhering chemicals.

### Manual cleaning

The generally recognized method is to wipe and rub the glass with a cloth or sponge soaked in cleaning solution. Abrasive cleaners and abrasive sponges should not be used on laboratory glassware as these can damage the surface of the glass. Surface damage can affect the glass properties and limit further use of the product. When soaking glassware it should generally be left in the cleaning solution for 20 to 30 minutes at room temperature, then rinsed with tap water followed by distilled water. To clean the glass as gently as possible, and thus extend its service life, a prolonged soaking time and higher temperature should only be used for stubborn soiling. Laboratory glassware should not be soaked for long periods in strongly alkaline media at more than 70 °C since this can have an adverse effect on the ceramic printing and may cause glass corrosion. Also to be avoided is severe mechanical action, e.g. scraping using a metal spoon.

### Washer-disinfectors for automatic laboratory glassware reprocessing

Washer-disinfectors for automatic laboratory glassware preparation are available in various sizes and performance classes. The product range extends from compact machines of 60 to 90 cm width up to powerful, large capacity machines. The large capacity machines are specially intended for central reprocessing of large quantities of laboratory glassware and are available as both 1-door and 2-door barrier machines for installation in a diaphragm wall.

60 cm wide compact machine  
Performance/load: e.g. 39 narrow neck glasses, 116 pipettes


115 cm wide large capacity machine  
Performance/load: e.g. 232 narrow neck glasses, 232 pipettes

Before purchasing a washer-disinfector, you must first be clear which laboratory glassware, and how much of it, requires reprocessing on a day-to-day basis within the laboratory. Once the machine size is specified, the appropriate accessories can be individually selected. The accessories include trolleys and inserts for secure support of the laboratory glassware. Inserts are primarily for holding wide-necked laboratory glassware. Special injector trolleys are offered to thoroughly clean laboratory glassware with a narrow internal diameter. These couple directly to the water supply of the machine and thus ensure that even internal cleaning of the laboratory glassware is correctly carried out. This system ensures that even difficult-to-access points are cleaned, which would be very difficult, or even impossible, to clean manually.

### Phases of machine-based reprocessing

Machined-based reprocessing comprises cleaning, rinsing, disinfection (if necessary) and drying of the laboratory apparatus. The following figure shows a typical programme cycle for laboratory glassware reprocessing.

#### Example

- Pre-rinse: cold water without process chemicals
- Cleaning: cold or hot water with alkaline cleaning agent
- Neutralisation: cold or hot water with acidic neutralisation agent

#### Cleaning

Cleaning removes dirt from the surfaces. At this stage, process chemicals (e.g. cleaning agents, surfactants, emulsifiers, neutralizers) are used. Cleaning may comprise several programme blocks, such as pre-rinse, cleaning, neutralisation.

#### Example

- Rinse I: cold water
- Rinse II: deionised or ultrapure water
- Flushing: deionised or ultrapure water at 75 °C

#### Rinsing

During rinsing the dissolved dirt and the process chemicals used are rinsed off. Rinsing can comprise a number of individual programme blocks. The choice of water quality (e.g. tap water, deionised water, ultrapure water) depends on the application (e.g. organic/inorganic analysis, microbiology).

#### Disinfection

During disinfection, infectious contamination is killed/inactivated to such a degree that the laboratory glassware no longer represents an infection risk. On the one hand, disinfection serves to protect personnel within laboratories who work with infectious contamination. On the other hand disinfection prevents transfer of germs from samples and preparations in medical laboratories, hygiene institutes, pharmaceutical laboratories and the food and cosmetic industries. Thus hygienic, problem-free working is guaranteed.

#### Drying

The washer-disinfectors have, dependent on model and construction, an active hot-air dryer which permits not only drying of the external surfaces, but also drying of narrow diameter laboratory glassware. Also laboratory glassware of complex shape is reliably dried using hot-air drying. To effectively protect the laboratory glassware against dust particles and microorganisms, the drying air is passed through a HEPA filter.

### Typical programme using a Miele washer-disinfector for reprocessing of laboratory glassware:

Miele washer-disinfectors for laboratory glassware reprocessing have up to 10 standard programmes. Numerous programme parameters can be adjusted to adapt the standard programmes for particular customer applications. Moreover, customer-specific programmes can be created for special applications.

Inorganic	To remove acid-soluble inorganic residues
Organic	To remove heavy organic residues such as oil, grease, wax, agar
Standard	Simple standard programme for slightly soiled glassware with a low final-rinse requirement
Universal	To remove organic residues (e.g. proteins, oils), for medium-level dirt and a medium final-rinse requirement
Intensive	To remove organic residues (e.g. proteins, cell and tissue cultures, oil), for heavy levels of dirt and a high final-rinse requirement
Plastic	For temperature-sensitive laboratory equipment (e.g. plastic bottles) with a low to medium level of dirt and a medium final-rinse requirement
Vario TD	For cleaning and heat disinfection at 93 °C with 5 minutes temperature-holding time, in accordance with EN ISO 15883-1, disinfection in the last rinse block
Special 93°C-10	For cleaning and heat disinfection at 93 °C with 10 minutes temperature-holding time, disinfection in the first rinse block, used in the case of an outbreak of a notifiable disease.

### Analysis purity through conductivity measurement in the final rinse

The requirements for analysis purity depend largely on the application of the laboratory glassware. To ensure analysis purity, washer-disinfectors for laboratory glassware reprocessing can optionally be provided with a conductivity measurement module. An integrated conductivity measurement offers the following advantages:

- Detection of undesirable contents in the rinse water (dissolved salts of alkaline or acidic process chemicals)
- Definition of a customer-specific permissible conductivity level

### Process reliability for reproducible results

Automatic preparation is a validatable preparation process that delivers reproducible results. This is one reason why automatic preparation should be favoured over manual processes.

To guarantee the reproducibility of the results, the machines have the following safety installations:

- Temperature monitoring using two redundant temperature sensors
- Automatic liquid dosing including dosing volume control
- Spray arm rotation speed monitoring

### Process documentation

In applications which require high standardisation and reproducibility, process documentation contributes significantly to quality control. Process documentation can take place via documentation software or a printer.

### Economy

Nowadays, laboratory glassware preparation must constantly meet ever higher requirements in respect of performance and economy. Machine-based reprocessing is by comparison with manual cleaning, much more efficient: for example, the economy arises from lower time / personnel expenses, shorter process cycles as well as lower power and water consumption. In particular, the short process cycles mean the laboratory glassware is quickly ready for its next use. Minimal handling of contaminated laboratory glassware simultaneously reduces the potential risk to personnel (injury, chemical burns and risk of infection).

### Value retention through gentle preparation

Automatic laboratory glassware preparation is gentler than manual cleaning. The glass surfaces only comes into contact with the alkalinity of the detergent for a short, defined time interval, so that glass corrosion is minimized. The accessories include special holders and locks so that the laboratory glassware is securely fastened and protected against breakage.

#### DURAN GROUP recommends Miele

To guarantee thorough, gentle and safe laboratory glassware preparation, DURAN GROUP recommends Miele washer disinfectors.

Miele "Made in Germany" quality is notable for its high reliability and efficiency in day-to-day use in the laboratory. Short operating times and reliable results ensure that high-value laboratory glassware is once again ready for use after only a short period.

In addition, the gentle preparation also ensures a long service life for DURAN® laboratory glassware.

DURAN GROUP

empfiehlt  
recommends

**Miele**  
PROFESSIONAL

## STERILISATION

When preparing laboratory glass for sterile applications or as part of the cleaning process, sterilisation is a well-established process. DURAN® laboratory glassware is suitable both for autoclaving as well as for hot air and plasma sterilisation ( $H_2O_2$ ). Laboratory apparatus that has come into contact with infectious substances or microorganisms must be cleaned in accordance with the appropriate guidelines for handling these materials. As the case may be, this may include sterilisation.

When carrying out sterilisation, especially of laboratory glassware, the following instructions should be observed: To avoid overpressure, all vessels should always be kept open. When sterilising media, the use of a membrane cap is recommended. Such a cap permits pressure equalisation through a PTFE membrane and hence the cap can be tightly closed. Consequently, the risk of contamination is greatly reduced.

Alongside the standardised procedures described above, individually modified methods are also applicable to all DURAN® products, for example using higher temperatures. However, you must ensure, especially with bottles (due to the screw caps) that the permissible highest temperatures for the plastic used in the accessories is not exceeded (see page 222).

## WORKING UNDER PRESSURE

Only products whose design includes the appropriate geometry and wall thickness, and which are explicitly designated as such, are suitable for working under pressure and / or vacuum (e.g. filtering flasks, desiccators or flat flange vessels).

When used under positive or negative pressure, and especially when also working with differential temperatures, additional care measures must be taken. Glass apparatus that is under pressure or vacuum should only be subject to further stress (e.g. significant temperature change) with extreme caution, as the individual resulting strains are additive and could readily result in failure.

To guarantee optimum user safety, the following points should be borne in mind:

- To avoid stresses in the glass, evacuated vessels or vessels under pressure should not be heated on one side or heated with an open flame.
- When working under pressure the maximum figures indicated in the catalogue should not be exceeded.
- Before using glass equipment under vacuum or pressure it must always be visually inspected to check that it is in perfect condition (no serious scratches, micro-cracks, abrasions, etc.). Damaged glassware should not be used for work under pressure or vacuum for safety reasons.
- Never subject glassware to sudden pressure changes, e.g. always re-pressurise evacuated glass apparatus slowly.
- Laboratory glassware with a flat bottom (e.g. Erlenmeyer and flat bottom flasks) should not be used under pressure or vacuum.
- The plastic coating of laboratory bottles (DURAN® protect) has no influence on pressure resistance. These products are not designed for use under pressure. For pressure applications using laboratory bottles, the DURAN® pressure plus bottle should be used. The DURAN® pressure plus bottle is pressure resistant from -1 to 1.5 bar due to a modified geometry and increased wall thickness.

## SAFETY INSTRUCTIONS

When used according to our specifications, DURAN® glassware is very safe to use. The appropriate guidelines applicable for the use of special glass in laboratories in the country in question should always be complied with.

The following points should, however, be observed in every case:

- For safety reasons, before DURAN® laboratory glassware is used it should be checked to ensure that it is suitable for the intended purpose and that it can be used without problem.
- Defective laboratory glassware represents a risk (e. g. risk of cuts, burns, infection) that should not be underestimated. If appropriate repairs to any item cannot be carried out or cannot be justified for economic reasons, it must be disposed of in the proper manner.
- Only subject DURAN® glassware to sudden temperature changes within the recommended limit for thermal shock resistance ( $\Delta T = 100K$ ). This means that hot laboratory glassware should not be taken out of a drying cabinet and placed on a cold or even wet laboratory bench. This applies in particular to thick-walled glassware such as filtration flasks and desiccators.
- When assembling apparatus make sure that it stands firmly and is not subjected to stress by using appropriate stands.

## DISPOSAL

DURAN® laboratory glass should under no circumstances be disposed of in the domestic glass recycling stream (e. g. bottle banks), since its high melting point and different chemistry make it incompatible with other glass cullet (soda-lime glass) for recycling. The correct way to dispose of it is, in principle, to include it with general household waste (residual waste) in accordance with the relevant guidelines, provided that the glass is quite free of any harmful contamination (Waste code no: 17 02 04).

## DURAN® LABORATORY GLASS IS ECO-FRIENDLY

DURAN® laboratory glass is made from natural, mineral raw materials. Unlike other materials, laboratory glass, when used properly, will give years of service and this means that it is vastly superior to other materials from an ecological viewpoint too. Depending on its use, DURAN® can be disposed of as household waste and does not need to be dealt with as special waste which may be environmentally harmful. Toxic substances cannot leach out because of the raw materials used.

Production processes in our factories have been consistently optimized over recent years to ensure that they are environmentally friendly during the actual manufacturing stage through the minimum usage of valuable resources. Electrical heating and advanced technology in our melting units ensure that no pollutants are released during manufacture in our ultramodern factories which could harm our workers or people living nearby. In addition energy demand is kept as low as possible. The latest waste gas purification equipment is used to avoid emissions which could pollute the environment. A significant investment has been made in an enclosed cooling water recirculation system to cut the amount of fresh water required to a minimum thus helping to conserve vital water resources. We use packaging made from environmentally harmless, recycled paper which can be returned after use to the resource cycle.


## FURTHER PROCESSING

DURAN® items made of borosilicate glass 3.3 are suitable for further processing such as the addition of screw thread tubes, olives, tubulatures, necks and ground glass joints. Preferred for further processing are round, flat bottom and Erlenmeyer flasks. Certain sections of the temperature / viscosity range

are of particular importance for glassworking. In the transformation range the elastic-brittle behaviour of the glass changes with increasing temperature into a markedly viscous one, so that consequently all its physical and chemical properties change significantly with temperature. The transformation temperature range thus plays an important part in stress relief during heating up and the introduction of stress when the glass is cooled. The position of the transformation range is identified by the transformation temperature "T<sub>g</sub>" DIN 52 324.

**Note:**

DURAN GROUP cannot accept any product liability where items are subjected to further processing. In this case the entire responsibility for quality lies with the glassworker. The latter is therefore responsible for ensuring that the further processed item conforms to current directives and safety requirements.

Normal temperature dependence/viscosity curve of, for example, DURAN®; viscosity ranges of important processing techniques, position of fixed points of viscosity and various limiting temperatures.

**DURAN® WITH INDIVIDUAL LABELLING**

Individual and permanent labelling of glass articles is now possible due to **innovative laser marking**. This system enables flexible labelling depending on the customer's requirements in the form of **texts, consecutive serial numbers, barcodes, logos, names or trade name of the laboratory**, etc. This information is processed with the aid of the common file format .tif. The contents are clearly identified by the labelling. Mix-ups in the laboratory can be ruled out, which is very important for sensitive areas such as the pharmaceutical industry or biotechnology. Laser marking is an ideal solution for labelling products. It enables the labelling of glass containers in different variants depending on requirements and complies with DURAN® quality requirements, as there **is no restriction of the product properties**. New, innovative technology also enables the labelling of **small batches**.

**Laser marking**

The laser marking is **burnt into the label field** and does **not** interact with the glass due to the wavelength used. Only the screen-printing ink is removed so that **the glass surface remains undamaged**. The tried-and-tested DURAN® glass properties such as high continuous usage temperature, resistance to temperature change and chemical resistance remain unchanged. The use of the latest laser technology produces **good print quality** and therefore **good legibility**. The lasered DURAN® glass articles are still **autoclavable/sterilisable and also microwave and dishwasher-safe**.

## BOTTLES

### Laboratory bottles

DURAN® laboratory bottles are chemically resistant and stable. The extensive range of original accessories includes screw caps for the widest possible range of applications. Alongside the standard PP screw cap for everyday laboratory use, further caps made from various plastics and having special properties are available. DURAN® laboratory bottles are completed by suitable pouring rings from different plastics, which enable drip-free working. As almost all GL 45 bottles of 100 ml capacity and above use the same thread size, screw caps and pouring rings are fully interchangeable. The bottles, pouring rings and caps are autoclavable/sterilisable.

#### Properties

##### Light protection

- amber bottles are opaque up to 500 nm
- plastic coated bottles are opaque up to 380 nm
- Application: storage of light sensitive substances

##### High thermal shock resistance

Due to their temperature properties, the bottles are suitable for autoclaving and sterilising (see page 227, general section). Because of the bottom geometry and the wall thickness, direct heating with an unshielded flame is not recommended. When using an electronic heating plate or water bath laboratory bottles should be heated gradually.

#### Recommendations

##### Pressure resistance

DURAN® laboratory bottles are, with the exception of the pressure-resistant DURAN® pressure plus bottles, in general not suitable for use under pressure or in a vacuum. DURAN® pressure plus bottles are pressure resistant from -1 to +1.5 bar (overpressure) due to a modified geometry and increased wall thickness.

##### Sterilisation

When sterilising or autoclaving contents, the screw cap must only be loosely fitted (max. one turn). The contents may expand or boil causing a large pressure difference in a closed vessel, which may well result in explosive failure. Alternatively, a DURAN® membrane cap may be used. Pressure equalisation takes place through the PTFE membrane, while at the same time the membrane cap can remain tightly closed, greatly reducing the risk of contamination. See also page 227, general part.

##### Cleaning

Cleaning should be carried out manually in a soaking bath or automatically in a dishwasher (see page 223, general part). When cleaning in a dishwasher, load so that there is no glass-to-glass contact (especially the threads) to avoid chips or abrasions.

##### Freezing substances

Recommendation: The bottle should be frozen slanted at an angle of 45°, filled to a maximum 3/4 (to enlarge the surface area) and dependent on the properties of any screw caps or other components used. For the blue PP screw cap the minimum temperature is -40 °C. Alternatively the Premium screw cap can be used (min. working temperature: -196 °C).

See also page 220, general part.

**Thawing frozen substances**

Frozen contents can be thawed by immersing the bottle in a liquid bath while taking care that the temperature difference between the contents and the bath does not exceed  $\Delta T = 100\text{K}$ . This will ensure that the frozen material is warmed uniformly from every side without damaging the bottle. The contents can, however, also be thawed slowly from above, so that the surface melts first, allowing the material to expand.

**Laboratory bottles with plastic coating**

The coating of DURAN® Protect bottles is a resistant and transparent plastic coating based on a cross-linked copolymer.

The coating adheres securely to the glass surface and fulfils the following functions:

- Protects the glass surface against mechanical damage (scratch protection).
- Holds the fragments together in the event of the glass breaking (splinter protection).
- Minimises liquid loss if the glass breaks (protects against contents escaping and splash).
- Absorbs UV rays up to a light wavelength of 380 nm (light protection).

**Recommendations**

- The plastic coating does not increase the pressure resistance. These bottles are not designed for pressure or vacuum applications.
- If the plastic coated bottle breaks during use, the contents and the plastic coating are likely to come into contact. A test for any interaction between plastic and contents should be carried out to ensure that the contents remain unchanged and can be further used.

**Temperature resistance**

Do not expose DURAN® protect bottles to open flames or direct heat, e.g. on a laboratory hotplate. The maximum operating temperature is  $+135\text{ }^{\circ}\text{C}$  and thus the bottle is suitable for use in an autoclave. Long-term exposure to temperature ( $> 30$  minutes) should be avoided. DURAN® protect bottles can be used for freezing to  $-30\text{ }^{\circ}\text{C}$  and used in microwaves. Thermal and chemical stresses can result in coating discolouration.

**Autoclaving:**

The following procedure, bearing in mind the maximum temperature resistance, is recommended:

- Steam sterilisation at  $+121\text{ }^{\circ}\text{C}$  or  $+134\text{ }^{\circ}\text{C}$ .
- The cycle duration should not exceed 20 minutes.

(See also page 227, general section)


When sterilising, the screw cap should only be loosely applied (max. one turn - do not tighten), or use a membrane cap that allows pressure equalisation.

### Gas washing bottles

By distributing the gas through the liquid by use of a filter disk, the gas surface is significantly increased and the interchange between gas and medium is improved. DURAN® gas washing bottles also work reliably at high flow velocities. The graph illustrates the effectiveness of gas washing bottles with and without a gas filter disk.

Absorption efficiency of two gas wash bottles: A without gas filter and B with gas filter plate

### Filtering flasks with side-arm socket or plastic hose connection

DURAN® filtering flasks are vacuum tight in accordance with DIN 12 476, ISO 6556. Alongside the filtering flasks with glass hose connections, versions are also available with a side-arm socket or plastic hose connection. The ground side-arm socket with dimensions 17.5/26 is suitable for vacuum hoses from 15 to 18 mm OD (e.g. 6 x 5 mm or 8 x 5 mm, DIN 12 865). The plastic hose connections are suitable for hoses of approx. 9 mm internal diameter. The versions with side-arm socket or plastic hose connection offer improved safety for the user.

### DURAN® SUPER DUTY

The new DURAN® Super Duty articles have greater mechanical stability compared to standard DURAN® articles due to increased glass content. The reinforced rim also increases shock resistance and considerably reduces the risk of breakage. They provide maximum possible safety for users when working under mechanical load e.g. frequent cleaning.

Uniform wall-thickness distribution, tried-and-tested DURAN® properties and increased shock resistance extend their service life and make DURAN® Super Duty glass containers more economical.

#### Recommendations

Uniform and slow heating is recommended for the Super Duty products to avoid thermal stresses in the glass. The standard DURAN® beakers and Erlenmeyer flasks should be used when working at very high temperatures or if rapid temperature changes are expected, as they are characterised by excellent resistance to temperature changes. However, the mechanical stability of these DURAN® products is limited compared to the Super Duty product range.

## DESICCATORS

DURAN® desiccators are used for drying moist substances or as storage vessels for moisture-sensitive products. To accelerate the drying process, the desiccators can be used under vacuum. Due to the high wall-thickness of the vessels and the exact machining of the vacuum-tight ground joints on the lid and base, storage under vacuum is possible - even over extremely long periods.

All individual parts and a wide range of accessories such as lids, stopcocks, bases, etc. are compatible and can be interchanged as required. Always ensure the individual parts have the same DN (nominal diameter in millimetres).

For desiccators, the DN is based on the diameter of the sieve plate; this, or the lip it rests upon in the desiccator base, can be measured directly. For lids, measure the outside diameter of the flange and cross-reference with the tables on the product pages. An overview on page 80 indicates which individual parts are required to assemble the desired desiccator.


### Recommendations

- Designed for use under full vacuum (-1 bar)
- Due to the high wall thickness and the reduced thermal shock resistance under pressure loading, the desiccators must not be heated on one side only or heated using a naked flame.
- Before evacuation, it is recommended that the glass surfaces of the desiccator be checked for damage such as scratches, cracks or nicks.  
Damaged desiccators must not be used for safety reasons.
- Never expose desiccators to abrupt pressure changes (do not suddenly ventilate evacuated vessels).

## CENTRIFUGE TUBES AND CULTURE TUBES

### Centrifuge tubes

DURAN® centrifuge tubes are approved in accordance with DIN 58 970 (Part 2) up to a maximum relative centrifugal acceleration (RZB = 4000) and for filling up to their capacity with contents having a maximum density of 1.2 g/ml.

Calculation:

$$RZA = 1,118 \times 10^{-5} \times r \times n^2$$

$$n = \sqrt{\frac{4000}{1,118 \times 10^{-5} \times r}}$$

Example:  $r = 15 \text{ cm}$

Example in the diagram:

number of revolutions (n)

= 4900  $\text{min}^{-1}$

### Culture tubes


In addition to DURAN® culture tubes, our product range also includes soda-lime culture tubes. This is a glass belonging to the third water resistance class and is one of the soda-lime glasses with a high fraction of alkaline and alkaline earth oxides.

Properties of soda-lime glass:

Physical data		Chemical data							
Linear expansion coefficient		Hydrolytic class		(ISO 719)					3
$\alpha_{20/300}$ to DIN 52328:	$9,1 \times 10^{-6} \text{ K}^{-1}$	Acid class		(DIN 12 116)					1
Transformation temperature Tg:	525°C	Alkali class		(ISO 695)					2
Temperature fixed points at viscosity $\eta$ n dPa x s:		Chemical composition							
$10^{13}$ upper annealing temperature	530 °C	(main components in approx. weight %)							
$10^{7,6}$ softening temperature	720 °C	SiO <sub>2</sub>	B <sub>2</sub> O <sub>3</sub>	K <sub>2</sub> O	Al <sub>2</sub> O <sub>3</sub>	Na <sub>2</sub> O	BaO	CaO	MgO
$10^4$ working temperature	1040 °C	69	1	3	4	13	2	5	3
Density $\rho$ :	2,50 g/cm <sup>3</sup>								

### DURAN® baffled flask with GL 45 thread


Oxygen intake is often the limiting factor for cell growth in the cultivation of microorganisms in Erlenmeyer flasks on a vibrating board. The movement causes a liquid sickle to form when using DURAN® Erlenmeyer flasks on a vibrator. The size of the sickle depends on the speed of the board and the vibration diameter. The greater the surface area of the contents, the greater the gas-exchange area and therefore the potential oxygen intake. The speed and the associated oxygen intake can, however, only be increased to a limited extent. The new DURAN® baffled flask with four baffles on the bottom disrupts the laminar flow and produces a turbulent flow. The surface area of the liquid and the gas-exchange area are increased, thereby increasing the oxygen intake. Laboratory trials have demonstrated that the oxygen intake is doubled by the baffles compared to a standard DURAN® Erlenmeyer flask .


The Erlenmeyer flasks with baffles from the DURAN Group can be geometrically reproduced due to completely automated and mechanical production. The wall thickness of the flasks was increased to achieve an excellent mechanical stability and to guarantee a long service life of the products. The special production process enables the manufacture of the product complete with thread in a two-stage process. The flasks can therefore be sealed with the tried-and-tested membrane screw cap from the DURAN Group. This enables a reproducible gas exchange compared to other sealing mechanisms e.g. sealing with cotton wool.

Liquid movement on a vibrating board:

DURAN® Erlenmeyer flask

DURAN® baffled flask

The movement causes a liquid sickle to form when using DURAN® Erlenmeyer flasks on a vibrator. The DURAN® baffled flask with four baffles on the bottom disrupts the laminar flow and produces a turbulent flow. The surface area of the liquid and the gas-exchange area are increased, thereby increasing the oxygen intake.


## FLAT FLANGE RANGE

The DURAN® flat flange reaction vessels are valued for their universal suitability for use in the laboratories of a wide range of specialisations. Whether for reaction, distillation, evaporation or desiccation, DURAN Group offers a wide range of unfinished and finished parts which always provide the optimum solution for the particular application. Due to the pure glass-glass connections, reactions with highly corrosive or highly chemically reactive substances can be carried out without problem.

The vessels are notable due to a robust glass flange design with an optimum flange angle of 45°. The proven flange design (flat ground) is available with groove, consequently O-Rings can be used. The corresponding stainless-steel quick release clamps with three flexible retaining clips ensure easy and safe handling. All individual parts and a wide range of accessories such as lids, O-Rings and quick-release clamps etc. are compatible and can be interchanged as required. In so doing however, you must always ensure the same DN (nominal diameter) of the individual parts applies.

### Recommendations

- All components are suitable for use under full vacuum (-1 bar). Many are rated for positive pressure operation (see product descriptions for details)
- Before use, it is recommended that the glass surfaces be checked for damage such as scratches, cracks or nicks.
- Damaged glassware should not be used for safety reasons.
- Due to the high wall thickness and reduced thermal shock resistance under pressure loading, the flat flange vessels should be heated uniformly and gradually.

Beaded lid for safer handling of the reaction vessel

## Accessories

Flat flange reaction vessels can be sealed by:

- a) O-rings (see below) for use at positive and negative pressures up to max. 230 °C (O-ring dependent)

Advantages:

- Easy to open
- The lid does not stick, even after operation for long periods under vacuum and at high temperatures
- Reduced need to grease contact surfaces

The stainless steel quick release clamps with three holding segments are optimally designed to provide even distribution of contact pressure. The chromium nickel steel support comprising two clamping rods is designed for secure fitting of the reaction vessels or the lids in support bar. For example, if there is a need to change the lid or the vessel, this can be done without dismantling the entire apparatus.

## Shape retentive O-rings

### FEP seamlessly coated elastomer O-rings with silicone core

Comprising an elastic, silicone core with a seamless FEP coating that encloses the ring. The combination of these high-quality materials ensures good elasticity in conjunction with outstanding chemical resistance. The chemical resistance of FEP (tetrafluoroethylene hexafluor-propylene copolymer) is equal to that of PTFE. Hence the material is resistant to almost all chemicals and is suitable for temperature from – 200 °C to + 200 °C.

### Silicone (VMQ) O-rings

These O-rings are made solely from silicone (VMQ) and therefore are highly elastic. Their chemical resistance, however, is reduced in comparison with FEP coated O-rings. Temperature resistance extends from –50 °C to +230 °C.

	Elasticity/recovery	Temperature resistance	Chemical resistance	Solvent resistance	Physiologically harmless
O-Rings, red FEP coated	+	++	++	++	++
O-Rings, transparent made of silicone (VMQ)	++	++	+	+	++

++ = very good resistance

+ = good resistance

## FILTERS AND FILTRATION APPARATUS

DURAN® filters and the corresponding filter plates are precision manufactured from DURAN® borosilicate glass 3.3 with its high chemical and thermal-shock resistance. They are entirely inorganic and inert in most circumstances. There are therefore no leachable organic or ionic species present that could otherwise contaminate the filtrates. They are ideal for separations, e.g. with strong acids or alkalis and can likewise be readily cleaned and reused. DURAN® filter products have a maximum operating temperature of +450 °C.

DURAN® filtration vessels are specially optimised to the matching filtration apparatus (eg funnels with gunko adapters) and are vacuum-tight due to their special geometry and high wall thickness. Their designs have been approved by the TÜV accreditation body and marked with the "GS" indication where appropriate; see specific products for details.

### DURAN® filtering apparatus

The filter apparatus has virtually universal applications with regard to the chemicals to be filtered because the medium only comes into contact with glass and PTFE. The graduated funnel simplifies dosing and analysis. The tried-and-tested DURAN® filtering flask and PTFE hose connection enable safe working in the laboratory. Thanks to the PTFE plate holder, porous glass plates with different porosities can also be used in addition to the split sieve. Filter paper, membrane filters (47mm) or just glass filters can be used for filtration. The replaceable plates and the PTFE adapter in conjunction with the clamp enable rapid changing of porosities or replacement of filters. Cleaning has been significantly simplified compared to a traditional filter funnel as the filter plate can be cleaned quickly and easily from both sides.

#### Recommendations

Coarse and fine and also analytical filtration can be carried out thanks to the available porosities of 10µm - 160µm. Furthermore, the filtration appliance is also suitable for the filtration of HPLC media, testing for bacterial contamination, residue analysis and the filtration of other media.

### Porosity

Porosity measurement is by the Bechhold bubble pressure method, which is widely described in the literature<sup>1</sup>. In the interests of rapid filtration every effort is made to produce filter disks with as many open pores as possible without blockages or closed cavities. This is one of the areas where DURAN® glass filters stand out.

Prerequisite for the successful use of glass filters is selection of the correct porosity. In this respect, the following table lists details of six porosity ranges with indications of their main areas of application. A point to be borne in mind is that the filtration equipment should ideally be selected to ensure that the nominal size of the largest pore is somewhat smaller than the smallest particles to be filtered out. This will prevent infiltration of particles into of the pores.

For quantitative analysis applications, porosity 3 or porosity 4 glass filtration apparatus is used almost exclusively. Different working methods often contain different porosity indications here for the same materials. This is because different processes used in the production of precipitations for gravimetric analysis often result in different grain sizes.

<sup>1</sup> Frank, W.: GIT (1967) Iss.7 pp. 683-688

Porosity	New identification ISO 4793	Nominal max. pore size $\mu\text{m}$	Areas of application, examples
0	P 250	160–250	<b>Gas distribution:</b> Gas distribution in liquids with low gas pressure. Filtration of coarsest precipitates.
1	P 160	100–160	<b>Coarse filtration, Filtration of coarsest precipitates.</b> <b>Gas distribution in liquids</b> Liquid distribution, coarse glass filters, extraction apparatus for coarse-grained material. Substrates for loose filter layers against gelatinous precipitates.
2	P 100	40–100	<b>Preparative fine filtration:</b> Preparative work with crystalline precipitates. Mercury filtration.
3	P 40	16–40	<b>Analytical filtration:</b> Analytical work with medium fine precipitates. Preparative work with fine precipitates. Filtration in cellulose chemistry, fine glass filters. Extraction apparatus for fine grained material.
4	P 16	10–16	<b>Analytical fine filtration:</b> Analytical work with very fine precipitates. (e.g. $\text{BaSO}_4$ , $\text{Cu}_2\text{O}$ ). Preparative work with correspondingly fine precipitates. Non-return valves and check valves for mercury.
5	P 1,6	1,0–1,6	<b>Ultrafine filtration</b>

### Flow rate

To determine the possible applications of glass filter disks and filtration apparatus, it is necessary to know not only the porosity, but also the flow rates of liquids and gases. These are given in Figures 9 and 10 for water and air. The data applies to 30 mm diameter filter disks.

The flow rates for other disk diameters can be calculated by multiplying the value read off by the conversion factor given in Table the following table:



Filter disk diam. mm	10	20	30	40	60	90	120	150	175
Conversion factor	0,13	0,55	1	1,5	2,5	4,3	6,8	9,7	15

**Example:**

Suction filtration of an aqueous solution under vacuum using a suction filter with a 60 mm disk diameter and porosity 4. Figure 9 gives a flow rate of 200 ml/min for a pressure differential of about 900 mbar. Table 8 gives a flow volume of  $200 \times 2.5 = 500$  ml/min for a 60 mm disk diameter. As the flow rate is heavily dependent on the pore diameter (pore radius to the power of 4), deviations from the values indicated may occur. Flow can also be obstructed by the formation of a filter cake over the surface of the filter disk. Further changes to the flow rate occur if liquids are used whose viscosity differs from that of water. The resultant flow rate is then inversely proportional to the viscosity. Differences for gases result when using filter disks that are coated with water or other liquids (gas flow in washing processes). More detailed information can be found in the literature<sup>1</sup>.

Water flow rate

Dry-air flow rate

Water flow rate through filter discs of various porosities as a function of pressure differential.  
For filter discs with Ø 30 mm

Air flow rate through filter discs of various porosities as a function of pressure differential.  
For filter discs with Ø 30 mm

<sup>1</sup> Frank, W.: GIT (1967) Iss.7 pp. 683-688

## Care and cleaning of filtration apparatus

In addition to the information in the general section on page 223, please also note the following guidelines relating to thermal stresses, which apply specifically to filtration apparatus, in order to avoid glass breakage.

### Temperature changes (thermal shock), drying and sterilisation

- The maximum permissible operating temperature is +450 °C.
- Uniform heating is recommended to avoid thermal stresses and resultant breakages.
- Heat glass filtration apparatus with disk diameters of more than 20 mm in initially cold ovens or sterilisers only.
- The heating or cooling rate should not exceed 8°C/min.
- When filtering hot substances avoid temperature differences of more than 100 K; if necessary, preheat the filtration apparatus in a drying cabinet.
- Wet filtration apparatus should be heated slowly up to 80 °C and dried for one hour before increasing the temperature further.

Whenever possible, filtration apparatus should be stood on its rim (stem upwards) to allow air convection between the inside of the vessel and the oven chamber. If placing the filtration apparatus in the oven at an angle cannot be avoided (as in the case of pipeline filters), any support point close to the position of the filter weld must be protected against heating up prematurely by placing heat-insulating material under it.

### Cleaning new glass filtration apparatus

Before using glass filtration apparatus for the first time, it should be rinsed with water (if applicable, acid), to remove any minor contamination that may be present.

### Mechanical cleaning

In many cases, if no precipitate has infiltrated the pores, simple spraying of the surface (e.g. with a spray bottle) will suffice. Brushes or rubber wipers can also be used to clean the surface of the filter disk. If some precipitate has infiltrated into the pores, then back-flushing of the disk is required.

### Recommendations

- Glass filters should always be cleaned immediately after use.
- Do not use sharp objects to remove the filtrate to prevent damage to the filter surface.

### Chemical cleaning

If some of the pores on the filter disk still remain clogged after mechanical cleaning or if it is desirable to make sure that no residue from previous work remains before filtering a new substance, then thorough chemical cleaning is necessary. The choice of solvent used depends on the nature of the contamination (see example in the following overview).

Barium sulfate	hot conc. sulfuric acid
Silver chloride	hot ammonia liquor
Red copper oxide	hot hydrochloric acid and potassium chlorate
Mercury residue	hot conc. nitric acid
Mercury sulfide	hot aqua regia
Albumen	hot ammonia liquor or hydrochloric acid
Grease, oil	acetone, isopropanol
Other organic substances	hot conc. sulfuric acid with addition of nitric acid, sodium nitrate or potassium dichromate

When chemical cleaning is completed, it should be followed by thorough rinsing with copious amounts of water. Use of hot concentrated phosphoric acid and hot alkali solutions is not recommended, as these may attack the glass surface.

### Screwfilters with interchangeable filter disks

With 3 filter sizes, each having 4 filter disks of varying porosity, 12 different filter rates are available. DURAN® screwfilters have a range of benefits compared with conventional filter apparatus:

- Interchangeable filter disks
- Safe and simple removal of the filtered material
- Disks have longer service life, as no damage is caused by scraping off the filtered material
- Filter disks are easy to clean from both sides
- Slit sieve (Cat. No. 21 340 31) can be used in the medium sized screwfilter to support membrane and paper filters
- Space saving
- Cost-effective; filter disks and apparatus can be ordered individually, as required.

### Recommendations

The filter disk should be located between 2 FKM gaskets.

## VOLUMETRIC PRODUCTS

DURAN® volumetric products have closely calibrated scales that permit very accurate determination and measurement of volumes. They are available in two accuracy classes: class A/AS and class B. The two classes differ in the accuracy of measurement with class A being the highest accuracy, and class B is approximately half that of class A. Class AS has the same tolerances as class A, but is designed to permit more rapid outflow; it is applicable to burettes and pipettes.

### Volumetric flasks

DURAN® volumetric flasks are manufactured from the chemically highly resistant borosilicate glass 3.3. Used for the accurate measurement of specific quantities of liquid they are, like virtually all volumetric glassware, volumetric analysis aids. They are mainly used for preparation and storage of standard solutions. Calibration is based on the amount of fluid contained ("In") at a +20 °C reference temperature, which means that when the circular graduation mark is reached, exactly the specified liquid amount is contained in the vessel. Thus the desired concentration can be precisely set. The volume content tolerances for volumetric flasks conform to accuracy class A, the accuracy limits of the German weights and measures regulation and to DIN and ISO guidelines.

### Measuring and mixing cylinders

DURAN® measuring and mixing cylinders are manufactured from borosilicate glass 3.3 and therefore are very resistant to mechanical and thermal stresses. Measuring cylinders are for holding and simultaneously measuring different liquid amounts. Mixing cylinders are for diluting solutions and mixing several components in a given quantity ratio. Their large hexagonal base prevents the cylinder from rolling. The base is equipped with three knobs that increase its stability. The cylinders have uniform wall thickness over the entire measurement range, so wedge errors are avoided. Calibration is based on contained fluid ("In") at a +20 °C reference temperature, which means that when the circular graduation mark is reached, exactly the specified liquid amount is contained in the vessel. Thus the desired concentration can be precisely set. Volume content tolerances for measuring and mixing cylinders conform to DIN and ISO accuracy limits.

### Burettes

DURAN® burettes are manufactured from chemically highly resistant borosilicate glass 3.3. They are primarily used for titration. The precise scale permits exact reading of the liquid quantity required for the titration. Calibration is based on the released volume ("Ex") at a +20 °C reference temperature. The fluid quantity released can be taken exactly from the scale, as the liquid adhesion to the glass is taken into account in the calibration. This only applies, however, if the specified waiting times for reading the scale are adhered to. Volume content tolerances for burettes conform to DIN and ISO accuracy limits. The DURAN® Class B burettes' accuracy limits are roughly one and a half times the Class AS accuracy limit. The tolerances are thus stricter than specified by DIN.

The tried-and-tested DURAN® burettes are also available with PTFE keys. Work in the laboratory is simplified by the fact that unlike glass keys, these do not have to be lubricated.

By the specification of a class "AS", the German weights and measures regulations have, within the scope of the 15th Amendment Regulations, acknowledged that the great majority of volumetric measurements, especially in clinical laboratories, are carried out with water or dilute aqueous solutions; thus apparatus with considerably shorter draining times than previously required but with the same accuracy limits is now admitted by the calibration regulations.

Capacity ml	Accuracy limits class AS suitable for official calibration DIN 12 700 ± ml	Accuracy limits class B	
		DIN 12 700 ± ml	DURAN ± ml
1	0.01	-	-
2	0.01	-	-
5	0.01	-	-
10	0.02	0.05	0.03
25	0.03	0.05	0.04
50	0.05	0.1	0.08
100 <sup>1</sup>	0.08	0.2	0.15

<sup>1</sup> Non-DIN size.

### Pipettes

Measurement and bulb pipettes are made from soda-lime glass (see page 234, for more information on soda-lime glass). Pipettes are for precise measurement and filling of liquids. Measurement pipettes are graduated to permit the taking up of varying liquid quantities and then dispensing of the same or different amounts. Bulb pipettes are designed to repeatedly take up and discharge a fixed volume for each pipette size. Calibration is based on the released volume (“Ex”) at a +20 °C reference temperature. The fluid quantity released can be taken exactly from the scale, as the liquid adhesion to the glass is taken into account in the calibration. This only applies, however, if the specified waiting times for reading the scale are adhered to. Volume content tolerances for calibrated pipettes conform to DIN and ISO accuracy limits. DURAN® Class B pipettes’ accuracy limits are roughly one and a half times the Class AS accuracy limit. The tolerances are thus stricter than specified by DIN.

By the specification of a class “AS”, the German weights and measures regulations have, within the scope of the 15th Amendment Regulations, acknowledged that the great majority of volumetric measurements, especially in clinical laboratories, are carried out with water or dilute aqueous solutions; thus apparatus with considerably shorter draining times than previously required but with the same accuracy limits is now admitted by the calibration regulations.

Capacity ml	Accuracy limits class AS suitable for official calibration ISO 385 ± ml	Accuracy limits class B	
		ISO 385 ± ml	DURAN ± ml
0.1 <sup>1</sup>	-	-	0.01
0.2 <sup>1</sup>	-	-	0.01
0.5	-	0.01	0.008
1	0.007	0.01	0.008
2	0.010	0.02	0.015
5	0.030	0.05	0.040
10	0.050	0.10	0.080
25	0.100	0.20	0.150

<sup>1</sup> Non-ISO size.

### Recommendations

- To ensure a long service life for your volumetric glassware and to exclude possible volume changes, these products should not be heated above +180 °C in drying cabinets or sterilisers.
- Never heat volumetric glassware on a hot plate.
- Always heat up and cool down volumetric glassware gradually, to avoid thermal stresses and thus any possible breakage of the glass.

## GLASS-CERAMIC LABORATORY PROTECTION PLATES

Due to low thermal expansion stresses, these glass ceramic plates are well suited to heating glassware with a Bunsen burner.

### Energy and time savings

The high transparency to infrared radiation means heat energy is transferred to the material being heated with low losses that shortens heating time and results in energy savings of 20% or more. In addition, several vessels can be placed on the plate's square, stable surface.

### Chemically resistant

When working in the laboratory it is impossible in practice to avoid aggressive media boiling over or spilling. The glass-ceramic laboratory protection plate is resistant even against highly corrosive media.

### Trouble-free cleaning

The pore-free smooth surface of the glass-ceramic laboratory protection plate can be cleaned easily either manually or mechanically.

### High temperature resistance

Service temperature from  $-200\text{ }^{\circ}\text{C}$  to  $+700\text{ }^{\circ}\text{C}$ . The glass-ceramic laboratory protection plate is may be used continuously at high temperatures. Durability at  $700\text{ }^{\circ}\text{C}$ : 6000 h; at  $750\text{ }^{\circ}\text{C}$ : 750 h. Even when a hot plate is quenched with cold water, there is no risk of breakage, since it is resistant to thermal shock even with a  $\Delta T > 650\text{K}$ . To avoid overheating, care must be taken not to exceed the above-mentioned limits when working with a Bunsen burner. The glass-ceramic laboratory protection plate retains its shape, remains flat and does not age.

**Note:** Further information about DURAN® laboratory glassware is available upon request.

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<b>A</b>		<b>C</b>	
Acid bottle	52	Calcium chloride cylinder	174
Adapter		Candles, filter micro	141–142
filter crucible	135	Caps, see screw caps	
receiver	174–175	glass	108
stopper type*		Kapsenberg	107
vacuum receiver	175	metal	110
Allihn		Centrifuge tube	99
condenser	166	Cone*	
filter tube	143	Coplin type staining jar	113
Analytical funnel	132	Clamp	
Aspirator bottle	52–54	closure for rolled bottle	111
Assembly set, for filtering flask, KECK™	146	quick release	165
Automatic burette	123	Cleaning scraper	71
<b>B</b>		Clip, KECK™	176–177
Beaker		Column, Vigreux	159
Berzelius	61	Condenser	
filtering, heavy-wall	59	Allihn	166
flat flange	162	bulb	166
high form	60–61	coil distillate	167
low form	59	Dimroth	167
Super Duty	59	jacketed coil	168
Philips	61	Liebig	166
Bearing, KPG stirrer	170–172	West	166
Bell jar	76–77	Cone	77
Berzelius beaker	61	conical flask	66
Bloom test vessel	62	GL 45	174
Blowball	123	Connection system	
Bottle		GL 45	36
aspirator	53–54	GLS 80®	21
culture	106	DG Safety Caps	40–42
culture media	108–110	DG Safety Waste Caps	41
gas washing	139–141	Coupling, screwthread	178
GL 45 laboratory bottle, amber	25	Crucible, filter	135
GL 45 laboratory bottle, clear	24	Crystallizing dish	72
GL 45 laboratory bottle, Premium	30	Culture bottle	106
GL 45 laboratory bottle, pressure plus	27	Culture flask	
GL 45 laboratory bottle, protect,		Erlenmeyer shape	109
PU-coated	26	Fernbach type	103–104
GLS 80® laboratory bottle, amber	16	Kolle type	105
GLS 80® laboratory bottle, clear	15	penicilin type	106
GLS 80® laboratory bottle, protect,		Roux type	105–106
PU-coated	17	Culture media bottle	108–110
HPLC	28	Culture tube	100
reagent, narrow neck	47, 49	Cylinder	
reagent, wide neck	46, 48	calcium chloride	174
rolled flange	111	measuring	120
Roller	112	mixing	119
screw cap, round, amber	44	multi purpose	76
soda lime	43–44, 48–49, 53	standing	76
stainless steel	31–32	<b>D</b>	
square	29, 43, 108	Desiccator	
weighing	72	base	85–86
Woulff	160	lid	87–88
Breed-Demeter, square bottle	108	O-ring	90
Buechner funnel	134	plate	89
Bulb pipette	126	screw cap	92
Burette	121–122	stopcock	90–91, 93
Burette bottle	124	DG Safety Cap	40, 42
		DG Safety Waste Cap	41

\*This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

Dimroth coil condenser	167	Flat bottom flask	
Dish		narrow neck	68
crystallizing	72	wide neck	68
evaporating	71	with ground joint	155
Petri	97–98	Flat flange*	
watch glass	73	beaker	162
Dispenser	39	lid	163–164
Disposable culture tube	101–103	quick release clamp	165
Drechsel head	139–140	reaction vessel	161
Dropping bottle	53	round bottomed flask	161
Dropping funnel	168–169	Funnel	
Drying tube, bent	173	analytical	132
DUROPLAN® Petri dish	97	Buechner	134
		Bunsen	132
<b>E</b>		dropping	168–169
Engler flask	69	filter	133–134, 142
Erlenmeyer flask		filter adapter	135
for Kapsenberg caps	107	filter, micro	142
narrow neck	63	short stem	131
screwthread (DIN thread)	65	long stem	132
wide neck	64	powder	131
with ground joint	156	ripped	133
Super Duty narrow neck	63	separating	170
Super Duty wide neck	64		
Evaporating dish	71	<b>G</b>	
Evaporating flask	157	Gas distribution tube	138
		Gas washing bottle	139–141
<b>F</b>		Glass box	114
Feigl type spot plate	98	Glass cap	108
Fernbach type culture flask	103–104	Glass ceramic laboratory protection plate	70
Filter apparatus universal	148–149	Glass stopper, amber	50
Filter apparatus, Witt type	150	GL 45, GL 32, GL 25, see bottles	
Filter candle, micro	141–142	GLS 80®	
Filter crucible / funnel adapter	135	bottle	15–17
Filter disk*		connection system	21
with glass rim	137	membrane cap	20
Filter funnel	133	pouring ring	19–20
Filter head	136	quick release screw cap	19
Filter, micro	141–142	stirred reactor	22–23
Filter, immersion	138, 141	GUKO	147
Filter tube, Allihn type	143		
Filtering apparatus	148–149	<b>H</b>	
Filtering flask	143–146	Head, Drechsel	139–140
FKM seals	137	Hellendahl type staining dish	113
Flask		Hexagon base measuring cylinder	119
baffled	103	Holding device for reaction vessel	165
conical	66	Hose connection, plastic	182
culture	106	HPLC bottle	28
Engler	69		
Erlenmeyer	63–65, 107, 156	<b>I</b>	
evaporating	157	Individual labelling	229
flat bottom	157	Immersion filter	138
flat flange	161	Immersion filter, micro	141
filtering	143–146	Iodine flask	66
iodine determination	66		
pear shape	155	<b>J</b>	
penicillin	106	Jacketed coil condenser	168
round bottom	67, 161, 158–159	Jar	74, 113
triple-neck round bottom	158–159		
twin-neck round bottom	158	<b>K</b>	
volumetric	119	Kapsenberg cap	107

\*This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

KECK™			
assembly set	146		
for conical joint*	176–177		
for spherical joint*	177		
Kolle type culture flask	105		
KPG stirrer bearing	171–172		
KPG stirrer shaft	172–173		
<b>L</b>			
Laboratory bottle, see bottles			
Laboratory protection plate, glass ceramic	70		
Levelling bottle	52		
Lid			
for desiccator	87–88		
for filter apparatus, Witt type	149, 151		
for flat flange vessel	163–164		
Liebig condenser (West condenser)	166		
<b>M</b>			
Measuring cylinder*	120		
Measuring pipette	124–125		
Membrane filter	28		
Membrane screw cap			
GL 25, GL 32, GL45	34		
GLS 80®	20		
Metal cap	110		
Micro burette	123		
Micro filter candle	141–142		
Micro filter funnel	142		
Micro immersion filter	141		
Mixing cylinder	119		
Museum jar	75		
<b>N</b>			
Narrow neck reagent bottle	47, 49		
NMR tube	79		
<b>O</b>			
Organ storage jar	73		
O-ring			
for desiccator	90		
for flat flange vessel	164–165		
<b>P</b>			
Pear shape flask	155		
Penicillin flask	106		
Petri dish	98		
Philips beaker	61		
Pipeline filter	138		
Pipette			
bulb	126		
for dropping bottle	53		
measuring	124–125		
Plastic hose connection	182		
Plastic stopper	50		
Plate holder	70		
Porcelain desiccator plate	89		
Pouring ring			
GL 32,	34, 44		
GL 45	33–34, 44		
GLS 80®	19–20		
Powder funnel	131		
Pregl, micro filter funnel	142		
Premium bottle	30		
Premium screw cap	33		
Pressure equalization set	21		
Pressure plus bottle	27		
Production bottle GL 45	18		
Production bottle GLS 80®	18		
Protect bottle	26		
Pump, water jet	151		
<b>Q</b>			
Quick release clamp	165		
Quick release closure, GLS 80®	19–20		
<b>R</b>			
Reaction vessel, flat flange	161		
Reagent bottle			
narrow neck	47, 49		
wide neck	46, 48		
Receiver adapter, bent	174		
Replacement seal	182		
Reservoir bottle	124		
RODAVISS®*			
Rolled flange bottle	111		
Roller bottle	112		
Round bottom flask			
narrow neck	67		
triple-neck	158–159		
twin-neck	158		
wide neck	67		
with ground joint	161, 168–169		
Roux type culture flask	105–106		
Rubber sleeve	135		
<b>S</b>			
Safety stopcock	93		
Schiefferdecker type staining dish	114		
Screw cap			
GLS 80® membrane screw cap	34		
GLS 80® quick release screw cap	19–20		
membrane screw cap	34		
Premium screw cap	33		
screw cap, PBT, PP	34–35		
tamper evident screw cap	33, 44–45		
with aperture	92, 180		
with two hose connectors	37		
Screwthread coupling	178		
Screwthread to ground joint adapters	179		
Screwthread tube	178–179		
Seal, silicone	55, 180–181		
Security screw cap	45		
Sedimentation cone	77		
Separating funnel	170		
Shaft, KPG stirrer	172–173		
Silicone rubber seal	180		
Silicone seal (Septa)	181		
Silicone sealing ring (VMQ)	181		
Socket* for spirit lamps see spirit lamp			
Soda-lime			
disposable culture tubes	101–102		

\*This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

funnel	131–132	Desiccator	83–84
pipette	124–125	receiver adapter, bent, baffled	174–175
reagent bottle, narrow neck	47, 49	Valves	
reagent bottle, wide neck	46, 48	GU*	
screw cap bottle, round	44	PRODURAN®*	
screw cap bottle, square	43	Vigreux column	159
staining dish	113–114	Volumetric flask	119
test tube	78		
Specimen jar	75	<b>W</b>	
Spherical joints*		Watch glass dish	73
Spirit lamp	71	Washing bottle, gas	139
Spot plate, Feigl	98	Water jet pump	151
Square bottle		Weighing bottle	72
Breed-Demeter	108	Wick for spirit lamp	71
laboratory bottle, DURAN®	29	Wide neck	
screw cap bottle, soda lime	43	GLS 80® bottle	15–16
Square quadrupod	70	reagent bottle	46, 48
Staining dish	113–114	Witt, filtration apparatus	150
Staining jar, Coplin type	113	Woulff bottle	160
Staining tray	114		
Stainless steel bottle	31–32		
Stainless-steel desiccator plate	89		
Stainless steel handle	115		
Standing cylinder	76		
STERIPLAN® Petri dish	97		
Stirred Reactor GL 45	38		
Stirred Reactor, GLS 80®	22–23		
Stirrer			
bearings, KPG	170–172		
shafts, KPG	172–173		
Stopcock			
single way*			
three way*			
two way*			
with PTFE-spindle for desiccator	90–91		
with standard ground joint for aspirator			
bottle	54–55		
Stopper			
for desiccator type WERTEX	84, 93		
glass	50–51		
plastic	50		
Super Duty			
beaker	59		
erlenmeyer flask wide neck	64		
erlenmeyer flask narrow neck	63		
<b>T</b>			
Tamper-evident screw cap	33, 44–45		
Test tube	78		
Threaded filter head	136		
Triple-neck round bottomed flask	158–159		
Tube			
centrifuge	99		
culture	100		
NMR	79		
test	78		
Tubing clamps, KECK™*			
Twin-neck round bottomed flask	158		
<b>V</b>			
Vacuum			

\*This article you will find in the DURAN® laboratory glassware for manipulators catalogue.

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Januar 2012

## GENERAL NOTE

The DURAN® laboratory glassware catalogue provides a basic information source for ordering our products. It does not represent a proposal for concluding a concrete agreement and will only serve as the basis for a contract upon explicit inclusion in a contractual relationship. We reserve the right to make changes to technical specifications, article numbers, packaging and design (e.g. due to changes of directives and standards). The contents of the catalogue have been created with the greatest possible care. However, we can accept no liability for the correctness, completeness and actuality of the contents. The presented replicated images provide an illustration of the article, details may however differ from the actual article.

## REGISTERED TRADEMARKS

DURAN®, DUROPLAN®, FIOLAX®, GLS 80®, KPG®, PRODURAN®, RODAVISS® and KECK™ are registered trademarks.

## TERMS AND CONDITIONS FOR DELIVERY AND PAYMENT

The following general terms and conditions for delivery and payment shall be applicable in respect of any and all deliveries and services by the Duran Group GmbH ("us") to customers provided that the customer receives these deliveries or services as part of its commercial or entrepreneurial activity (the "Customer"). Any conflicting general purchasing terms and conditions of the Customer are hereby expressly rejected. Any such general purchasing terms and conditions shall apply only if we expressly confirm them in writing.

### 1. Prices/Terms and Conditions of Payment

**1.1** Unless otherwise agreed, the applicable prices are calculated in EURO (EUR), plus an additional amount for VAT as applicable from time to time. Unless special terms are agreed, the prices should be understood to be prices ex works, with no deduction or discount being granted for immediate payment.

**1.2** If payment deadlines or dates specified in our order confirmation or otherwise agreed are not met, this will automatically give rise to all of the statutory consequences of default, without any special reminder being required. In particular, we reserve the right to charge interest at the applicable rate charged by our bank for utilised credit if such interest exceeds the interest rate prescribed by statute (8 percentage points above the base lending rate). Furthermore, the entire balance shall become due and payable immediately, irrespective of any payment targets.

### 2. Delivery Dates and delivery amounts

**2.1** We will endeavour to adhere to stipulated delivery deadlines. However, due to the hazards and peculiar features of glass processing, delivery deadlines will not be binding unless expressly agreed otherwise. Our delivery times are subject to our suppliers delivering the correct products to us on time. We undertake to inform the Customer of any unavailability of any of our products without undue delay and will reimburse the Customer any amounts paid in respect of the unavailable products.

**2.2** In the case of any custom-made products, we reserve the right to deviate to a reasonable extent from the agreed quantity. The Customer must take delivery of surplus quantities. A variation of +/- 10% in relation to the ordered quantity shall be deemed as the agreed tolerance.

### 3. Place of Performance and Passing of Risk

**3.1** The place of performance for the delivery is the principal place of business of our respective supplying factory. The place of performance for payment is our principal place of business.

**3.2** When goods are transported, the risk (of accidental loss, destruction or deterioration) (the "Risk") shall pass to the Customer as soon as we have delivered the goods to the carrier chosen by us.

### 4. Packaging

Unless otherwise agreed, we will accept the return of packaging only to the extent that we are obliged to do so under the German Packaging Regulation (Verpackungsverordnung) or other mandatorily applicable legal regulations.

### 5. Payment

**5.1** Unless agreed otherwise, our claim for payment of the purchase price becomes due immediately following receipt of the relevant invoice.

**5.2** We reserve the right to assign any claim we may have against the Customer in whole or in part to a third party.

**5.3** Any right for set-off or any right of retention may only be asserted by the Customer in respect of undisputed or finally determined and legally binding claims.

**5.4** The Customer shall, irrespective of any other claims for reimbursement of costs we might have, be obliged to assume any fees, costs and expenses that accrue due to a legally successful enforcement of rights against the Customer outside of the Federal Republic of Germany provided these fees, costs and expenses were required for the enforcement action.

### 6. Warranties in Respect of Defects and Notification of Defects

**6.1** If, despite the greatest of care being taken, the goods give rise to complaints, then, in accordance with § 377 of the German Commercial Code (Handelsgesetzbuch, or "HGB"), obvious defects must be notified without delay, and in any case no later than 14 calendar days after receipt of the goods, and hidden defects must be notified without delay after their discovery, otherwise the goods shall be deemed accepted.

We shall not be liable for damage to deliveries through the breakage of glass during transit ("Breakages") where the cause of the damage arises after the transfer of risk. In such cases any claim for Breakages must be made against the carrier or under the policy of transit insurance. We shall not pay compensation for Breakages where the value of the relevant item is EUR 20.00 or less unless the Breakage is due to intentional conduct or gross negligence on the part of ourselves or our servants or agents.

**6.2** Claims on the basis of defects as to quality shall become time-barred 12 months after delivery of our goods to our Customer. The foregoing provisions shall not apply to the extent that longer limitation periods are mandatorily prescribed by statute pursuant to § 438(1) No. 2 of the German Civil Code (Bürgerliches Gesetzbuch, or "BGB" – Physical Structures and Physical Objects used for Physical Structures), § 479 (1) BGB (Recourse Claim), and § 634a (1) BGB (Construction Defects).

**6.3** Delivered goods which are returned to us because the wrong goods were delivered or due to a defect ("Returned Goods") shall only be accepted if we are notified of the Returned Goods before their dispatch and the following conditions are satisfied:

- a) Upon notification of any Returned Goods, the Customer shall be issued a processing number relating to the Returned Goods; such processing number must be entered on the documentation for the returned items;
- b) Any Returned Goods must be reported to our freight centre by delivering appropriate carriage documents with a reference to the processing number

relating to the Returned Goods attached.

**6.4** If, despite all care being taken, the delivered goods contain a defect that already existed at the time that the Risk passed, then we will, in our sole discretion and subject to receiving notification of the defect within the required time period, repair the goods or deliver substitute goods. We must always be given the opportunity to render supplementary performance (Nacherfüllung) within a reasonable time period.

**6.5** If the supplementary performance fails to rectify the defect, the Customer may – notwithstanding any claims for compensatory damages – rescind the agreement or reduce the amount of the purchase price.

**6.6** The following shall not give rise to any claims based on defects: merely immaterial deviations from the agreed condition of the goods, merely immaterial impairments to their utility, natural wear and tear, or loss or damage that arises after the Risk has passed as a result of incorrect or careless treatment, overuse, unsuitable operating resources, defective building work, unsuitable building foundations or special external influences that are not included or catered for in the contract. In addition, if the Customer or a third party improperly (in a non-workmanlike manner) carries out maintenance work on or makes modifications to the goods, then no further claims based on defects may be made in respect of such goods or the resulting consequences of such maintenance work or modifications.

**6.7** Claims on the part of the Customer for expenses necessary to enable supplementary performance, particularly transport, tolls and other road and transport charges, labour costs and the cost of materials, are excluded to the extent that such expenses are increased because the goods delivered by us were subsequently taken to a location other than the Customer's business premises, unless such transportation is consistent with the authorised use of the goods.

**6.8** Any recourse claims on the part of the Customer against us shall exist only to the extent that the Customer has not entered into any agreements with its customers going beyond the mandatory statutory claims regarding defects. Clause 6.6 shall apply accordingly in respect of any such recourse claim by the Customer against us.

### 7. Industrial Property Rights and Copyright; Title Defects

**7.1** Unless otherwise agreed, we have an obligation (although such obligation exists only in the country in which the place of delivery is located) to deliver the goods free from the industrial property rights and intellectual property rights of third parties (hereinafter referred to as "Proprietary Rights"). In the event that a third party makes legitimate claims against the Customer for infringement of Proprietary Rights based on the goods delivered by the supplier and used in accordance with the contract, we shall be liable to the Customer within the period specified in clause 6.1 above as follows:

a) In our sole discretion and at our own expense, we will either secure a licence for the goods concerned, modify them so that the Proprietary Right is not infringed, or exchange them. If we are unable to do any of the above on reasonable terms, then the Customer shall be entitled to the statutory rights of rescission and reduction of the purchase price.

b) The provisions of clause 8 shall apply to any claims for compensatory damages or claims for the reimbursement of expenses.

c) Our obligations as described above shall exist only on the condition that the Customer notifies us in writing without delay of the claims asserted by the third party, the Customer does not admit to the infringement and leaves in our hands any defence of the claims and settlement negotiations. If the Customer discontinues using the delivered goods in order to mitigate loss or for any other good reason, then the Customer shall notify the third party of the fact that discontinuing use of the goods in no way constitutes an admission of an infringement of Proprietary Rights.

**7.2** Claims on the part of the Customer are excluded if the Customer is responsible for the infringement of the Proprietary Rights.

**7.3** Claims on the part of the Customer shall be further excluded if the infringement of the Proprietary Rights is a result of special instructions issued by the Customer, an application or use of the goods that was not foreseeable by us, or as a result of the Customer modifying the goods or using them together with goods not delivered by us.

**7.4** In the event of an infringement of a Proprietary Right and regarding claims by the Customer arising according to clause 7.1 a), the provisions set forth under clauses 6.3 and 6.7 shall otherwise apply accordingly to the Customer's claims.

**7.5** If other title defects exist, then the provisions of clause 6 shall apply mutatis mutandis.

### 8. Claims for Compensatory Damages; Limitation of Liability

**8.1** In the event of a breach of a pre-contractual, contractual and/or other obligation, including unsatisfactory delivery, tortious conduct and manufacturer's liability, we shall be liable for compensatory damages and the reimbursement of costs – subject to further contractual or statutory liability requirements – only in the case of wilful conduct or gross negligence and in the event of a breach of a material contractual duty only (i.e. being a contractual duty, the infringement of which jeopardises the ultimate purpose of the contract and whose fulfillment the Customer can under regular circumstances expect) also due to ordinary negligence. However, our liability for simple and gross negligence as well as in the event of liability that arises regardless of negligence or fault, shall be limited to typical contractual loss or damage that was foreseeable at the time the contract was entered into.

**8.2** The exclusions and limitations of liability set forth under clause 8.1 shall not apply in the event that a guarantee is given within the meaning of § 443 BGB

with respect to the condition of the goods at the time the Risk passes to the Customer or the durability of the goods (i.e. a declaration by the seller that the object of the purchase as of the time the Risk passes possesses a certain quality or will maintain a certain quality and that the seller is willing to assume responsibility for any consequences arising from the fact that such quality does not exist regardless of negligence or fault), or a defect is fraudulently concealed, in the event of injury to life, physical injury or injury to health, or mandatory liability under the German Product Liability Act (Produkthaftungsgesetz). In the event of fraudulently concealing a defect or in respect of any guarantee pursuant to § 443 BGB, the Customer's rights shall solely be determined according to the statutory law or the content of the guarantee.

**8.3** Irrespective of the Customer's claims regarding compensatory damages and the reimbursement of costs set out in clause 8.1, any further claims or other claims than the rights set out in clauses 6 and 7 regarding any defect or title defects by us or against any of our agents shall be excluded.

### 9. Non-binding Nature of Drawings, Diagrams, Measurements and Weights

Drawings, diagrams, measurements and weights are approximate only, unless they are expressly stipulated to be binding. The Customer must guarantee that working drawings (construction diagrams) supplied by it do not infringe the Proprietary Rights of third parties. The Customer must hold us harmless in the event that rights of recourse are asserted by third parties.

### 10. Documents

Documents supplied by us may not be copied or made available to third parties, or used for any purpose other than the agreed purpose.

### 11. Reservation of Title

**11.1** We shall retain title to the goods until all of our claims, including claims arising in the future, are fully paid. The Customer may process and sell the goods in accordance with the following conditions: If the goods are further processed or remodelled by the Customer, then we shall be deemed the manufacturer within the meaning of § 950 BGB and shall acquire direct title to the intermediate or final products. As a precaution, the Customer hereby assigns and transfers the ownership of any new goods created by further processing or remodeling any goods delivered by us to us subject in each case only to the execution of the relevant purchase contract. In respect of such goods assigned and transferred to us, the Customer shall be merely the custodian or bailee of such goods. If the goods subject to the reservation of title ("Reserved Goods") are mixed or processed with other property not belonging to us, then we shall acquire a co-ownership interest in the new item proportionate to the value of the Reserved Goods to the other property.

**11.2** The goods may be sold only in the normal and ordinary course of business and only if claims deriving

from their resale are not assigned to third parties beforehand. The Customer's claims deriving from a resale of the Reserved Goods are hereby assigned to us subject only to the execution of the purchase agreement between us and the Customer, this assignment shall also include any right arising from the fact and to the extent that these goods are mixed or combined with other property. In such a case, the assigned claims shall serve as our security only up to the value of the Reserved Goods sold in each case. We will not collect on the assigned claims for as long as the Customer complies with its payment obligations. However, the Customer has an obligation to disclose to us the identity of the third party debtor at our request and to notify such debtor of the assignment. The Customer may collect on the claims resulting from sale of the Reserved Goods unless and until it receives instructions from us to the contrary. The Customer must immediately transfer any amounts collected by it to us if, to the extent that and as soon as our claims are due.

**11.3** Pledges or the granting of security interests or any assignment of the Reserved Goods or the assigned claims are not permitted. The Customer must inform us immediately of any action by third parties affecting the Reserved Goods or the assigned claims. We agree to release the assigned claims in our sole discretion if they exceed the value of our claims to be secured by more than 20% and are derived from fully paid deliveries.

**11.4** In the event of a breach of duty by the Customer, particularly in the case of default on payment, we are entitled to rescind the agreement in whole or in part and recover the Reserved Goods. The Customer has an obligation to deliver up the Reserved Goods. The declaration of recovery or the enforcement of the reservation of title or any seizure of the goods by us constitute a declaration of rescission from the agreement with respect to the Reserved Goods.

**11.5** If, in the case of non-domestic sales, the reservation of title agreed under clause 11 is not permitted with the same effect as under German law, then we shall retain title to the goods until payment of all of our claims arising out of the contractual relationship formed through the sale of the goods. If the foregoing reservation of title is not permitted with the same effect as under German law either, but it is permissible to reserve other rights in respect of the goods, then we are authorised to exercise all of these rights. The Customer shall cooperate in all actions we may wish to take in order to protect our ownership interest or alternative right in the goods.

## 12. Return of Goods

Any acceptance of a return of goods and any repayment of the purchase price relating to such goods shall be in our sole discretion and under the proviso that we are not legally obliged to do so. The following rules shall apply to any goods that are returned to us unless the goods are Returned Goods within the meaning of clause 6.3:

a) Any goods that are returned must have been

purchased within 4 weeks in the case of deliveries within the Federal Republic of Germany or within 8 weeks in the case of deliveries to customers situated in Europe or within 12 weeks in the case of deliveries to customers situated outside of Europe. The time limits commence running on the date that the goods have been delivered at the Customer and expire on the date of receipt of the returned goods.

b) The provisions of clause 6.3 shall apply accordingly to the acceptance, notification and labeling of goods that are returned to us.

c) Only unopened and undamaged goods without additional stickers or labeling attached to them shall be accepted. We must be able to resell the goods.

d) Any return of goods shall be at the Customer's sole cost and risk.

e) We shall also charge a handling fee equivalent to 20% of the value of the item returned subject to a minimum charge of EUR 20.00 per return. Such sums shall be deducted from an amount that is being reimbursed to the Customer.

f) Custom-made products may not be returned.

## 13. Applicable Law and Judicial Forum

**13.1** With the exception of conflict of law rules under German private international law and the provisions of the UN Convention on Contracts for the International Sale of Goods ("UN-CISG"), the substantive law of the Federal Republic of Germany shall apply to all legal relationships with the Customer.

**13.2** Sole place of jurisdiction for both parties regarding all legal disputes arising out of the relevant purchase contracts or in connection with the supply relationship, including matters of currency exchange, is our head offices. If we appear as the plaintiff, we are also entitled to bring an action before the court responsible for the Customer's head office.

## 14. Moulds and tools

Moulds and tools produced on behalf of the Customer, whether by us or sourced from third parties shall remain in our ownership and possession. At the start of the contract, the Customer shall pay the agreed mould and tool contribution which grants the right to exclusively be supplied from these moulds. At the end of the contract, or any other discontinuation of the project, no assignment or transfer of the moulds and tools will take place; they will remain our property of, and in our possession. In these cases, however, the Customer shall be entitled to demand that we scrap the moulds and tools at our own expense and provide evidence of the scrapping to the Customer. An obligation by us to store project-related moulds and tools shall end automatically at the end of the contract or project. If there is no written agreement to the contrary, a project shall be deemed to have ended after the expiry of a two-year period after our confirmation of the Customer's last order.

We shall ensure proper storage, handling and maintenance of the moulds and tools within the usual scope, during the term of the project. If the moulds

or tools are destroyed or damaged due to improper storage, handling or maintenance by us then they shall be repaired or newly acquired at our expense. The same applies to loss, destruction or damage as a result of force majeure. In the case that moulds and tools are used beyond their limit of wear and tear, the Customer shall bear the costs of the new moulds and tools to be acquired by us, up to the amount of the originally agreed cost contribution for the worn part. Should the limit of wear and tear be reached prior to reaching an output quantity individually guaranteed, or if the Customer proves that the wear and tear is due to a fault of the mould or tool, or an operating error by us, then we will bear the full cost of replacement. The above provisions shall apply accordingly to the moulds and tools acquired as replacement.

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
“RESEARCH? ... WHAT IT MEANS TO ME:  
**REMAINING CURIOUS**, ALWAYS ASKING NEW  
QUESTIONS ... AND WHEN IT'S A QUESTION OF  
GLASSWARE, THEN I ALWAYS ASK FOR DURAN®.”

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*OUR INPUT: LABORATORY GLASSWARE.*  
***YOUR FEEDBACK?*** *WHATEVER YOU WANT*  
*TO TELL US. WHETHER IDEAS FOR*  
*APPLICATIONS OR QUESTIONS. WE'RE*  
*LOOKING FORWARD TO YOUR CONTRIBUTION:*  
[\*INFO.DURAN@DURAN-GROUP.COM\*](mailto:INFO.DURAN@DURAN-GROUP.COM)

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## WHY DURAN® IS SPECIAL

*Optimum physical and chemical characteristics make DURAN® glass the material of choice for laboratory applications:*

- *Uniform wall thickness distribution throughout the product*
- *Therefore better mechanical stability and greater thermal shock resistance*
- *Prevents material stresses, which may cause the glass to crack when subject to heating or cooling*
- *Benefits: improves safety for laboratory technicians, increases the lifespan of laboratory glassware and protects valuable substances*

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